



## Green Fluorescent PIP<sub>2</sub> Assay

Overview	2
Relevant Products	2
Materials in the Kit	2
Storage	2
Additional Materials Not Supplied	2
BioSafety Considerations	2
Warranty	3
About these Assays	3
Related Products	3
Protocol for Use	4
Suggestions for Adherent Cells	6
Fluorescence Detection	6
Assay Performance Considerations	7
Trouble Shooting	8
Contact Us	8
References	9

## Overview

The green PIP<sub>2</sub> sensor is based on a dimerization-dependent fluorescent protein. This protocol is optimized for imaging rapidly dividing, immortalized cell lines on a 96-well plate and has been validated in live HEK293, CHO and NIH 3T3 cells. For imaging live iPSC-derived or cultured cells, see Additional Information section.

The BacMam vector carrying these sensors is a modified baculovirus, which can be used for delivery to, and expression in, a wide variety of mammalian cell types including primary cultures.

## Relevant Products

Product	Description	Promoter	Recommended Use
D0400G	Green PIP <sub>2</sub>	CMV	Fluorescence imaging

## Materials in the Kit

- PIP<sub>2</sub> sensor BacMam ~ 6x10<sup>9</sup> VG/mL in TNM-FH Insect Culture Medium (Allele Biotech product #ABP-MED-10001).

Green fluorescent sensors that decrease in fluorescence intensity when PIP<sub>2</sub> levels decrease. VG/mL is the titer determined by qPCR, and is the average number of viral genes per mL of the BacMam stock.

- Sodium Butyrate (Sigma Aldrich product # B5887) 500 mM in H<sub>2</sub>O.

Sodium Butyrate is added to the culture to maintain BacMam expression. Other histone deacetylase (HDAC) inhibitors, such as Trichostatin A (TSA) may work too.

- HM1 Muscarinic Receptor BacMam with nuclear red label in TNM-FH Insect Culture Medium (Allele Biotech product #ABP-MED-10001).

A Gq-coupled receptor in BacMam provided as a positive control for the purpose of assay optimization. Contains a separate red fluorescent protein that is targeted to the nucleus.

- Carbachol 25 mM in H<sub>2</sub>O

Carbachol can be used to stimulate Gq signaling through the positive control, the HM1 Muscarinic acetylcholine Receptor.

## Storage

BacMam stocks should be stored at 4°C and protected from light. Avoid freeze/thaw cycles.

## Additional Materials Not Supplied

- Greiner CellCoat (#655946) is our preferred 96-well plate available from VWR.
- Dulbecco's Phosphate Buffered Saline (DPBS) available from VWR [Dulbecco, R. and Vogt, M.1957].

## BioSafety Considerations

BacMam is the modified baculovirus, *Autographa californica*, AcMNPV. The baculovirus contained in this kit is pseudotyped to infect mammalian cells. In mammalian cells, the baculovirus genome is silent and it cannot replicate to produce new virus in mammalian cells. While it should be handled carefully, in a sterile environment, it is classified as a BSL-1 reagent.

This product is for research use only and is not recommended for use or sale in human or animal diagnostic or therapeutic products.

## Warranty

Materials are provided without warranty, express or implied. End user is responsible for making sure product use complies with applicable regulations. No right to resell products or any components of these products is conveyed.

## About these Assays


PIP<sub>2</sub>, also known as Phosphatidylinositol 4,5-bisphosphate or PtdIns(4,5)P<sub>2</sub>, is a substrate for several signaling molecules. It is a precursor for the DAG/IP<sub>3</sub> pathway, but also plays a key role in phospholipid signaling and is involved in the regulation of ion channels and transporters. The fluorescent sensors used in the assays described here can be combined with different colored sensors, such as the red GECO calcium sensor or a red DAG sensor to measure multiple signals simultaneously.

## Related Products

Product	Sensor Description	Promoter	Recommended Use
U0300R	Red Upward DAG	CMV	Fluorescence imaging and Plate reader assay (Z' > 0.5)
D0300R	Red Downward DAG	CMV	Fluorescence imaging and Plate reader assay (Z' > 0.5)
U0600R	Red GECO Ca <sup>2+</sup>	CMV	Fluorescence imaging and Plate reader assay (Z' > 0.5)

## Protocol for Use

This protocol is optimized for rapidly dividing, immortalized cell lines, however, the protocol can be adjusted for transducing non-dividing cells adherent cells such as neurons, islets, cardiomyocytes, and iPSC-derived lines. We recommend that you take the time to optimize the assay for your particular cell type. See our Suggestions for Adherent Cells section following this protocol.



**Example:**

**96 wells** (1 plate)

100 uL cell suspension (500,000 cells/ml) per well.

100 uL cells x **110** (**96 wells** + 10% scale) = **11000 uL** cell suspension.

- When preparing the master mix, scale up by 10-15% to avoid coming up short. To seed a 96-well plate, multiply amounts in Step 1 and Step 2 by 110-120.

**Step 2) Prepare Viral Transduction Reaction (Tube B)**

- Prepare a 500 mM stock solution of sodium butyrate in sterile water (in your kit).
- For each transduction reaction (i.e. one well in a 96-well plate), prepare the transduction solution by mixing 25 uL of the BacMam stock with 0.6 uL of the 500 mM stock solution of sodium butyrate\*\* and 24.4 uL of the complete culture media for your cells, for a total volume of 50 uL. Mix gently.

\*\* Concentration of sodium butyrate should be 6mM in this step. Following Step 3, final concentration of sodium butyrate will be 2mM.

**Example:**

96 wells needed (1 plate). The number of wells desired, in bold, corresponds to the number in Step 1 above.

<i>Single Well</i>	<i>Master Mix</i>
20 uL Sensor	x <b>110</b> = 2000 uL
5 ul HM1 Receptor Control	x <b>110</b> = 550 uL
0.6 uL 500 mM Sodium Butyrate	x <b>110</b> = 66 uL
<u>24.4uL Complete Media</u>	<u>x <b>110</b> = 2684 uL</u>
50 uL Total Volume	x 110 = <b>5500 uL</b> transduction mix (96 wells)

**Step 3) Mix Cells and Transduction Mix from above.**

- Mix Tube A and Tube B (100 uL tube A + 50 uL tube B). Mix gently and then seed 150 uL of mix per well on the 96-well plate.
- Cover plate with aluminum foil to protect from light and incubate at room temperature for 30 minutes.
- Incubate overnight under normal growth conditions, protected from light.

**Example:**

**96 wells** needed (1 plate)

<i>Single Well</i>	<i>Master Mix</i>
100 uL cell suspension	x 110 = 11000 uL
<u>50 uL transduction reaction</u>	<u>x 110 = 5500 uL</u>
150 uL total volume per well	x110 = 16,500 uL total reaction volume

## DAY 2 FLUORESCENCE IMAGING

- Cells are now ready for assay. Prior to imaging, replace culture media with DPBS. Wash gently so as not to dislodge cells. **Cover the cells and allow them to rest at room temperature in DPBS for 20-30 minutes before measuring fluorescence.** Experiments are performed at 25°C using standard GFP excitation and emission wavelengths.
- Add 50uL of 200uM carbachol (50 uM final concentration per well) to activate the DAG/IP<sub>3</sub> pathway in a set of control wells. A decrease in the fluorescence intensity will be observed after addition of the carbachol when PIP<sub>2</sub> is hydrolyzed to produce IP<sub>3</sub> within the cell.

## Suggestions for Adherent Cells

The protocol above is optimized for rapidly dividing immortalized cells, however, this product is compatible with screening primary cultures and iPSC-derived lines, where the cells are plated before transduction. Specific protocols will vary by cell type, so it is important to take the time to titrate the BacMam stock for best results. For expression in rare cell types, or specific cells in mixed cultures, Cre-dependent and specific promoter systems are available for many of our sensors.

- Prepare a 500 mM stock solution of sodium butyrate in sterile water.
- For each transduction reaction (e.g. one well in a 96-well plate, 100 uL per well), prepare the transduction solution by mixing 25 uL of the BacMam stock with 23.5 uL of DPBS and 1.5 uL of the 500 mM stock solution of sodium butyrate for a total volume of 50 uL. Mix the solution gently.
- Sensor expression and cell health can be controlled by titrating the virus, so it is worth taking the time to optimize the assay for your particular cell type.
- Prepare a dilution series of transduction reactions by varying the amount of BacMam ranging from 15 uL to 80 uL and adjusting the amount of DPBS accordingly.
- Add the transduction reaction directly to the plated cells (no aspiration of cell medium necessary). Gently rock the plate 4-5 times in each direction to mix throughout the well. Incubate the cells under normal growth conditions, protected from light, for 6 hours (5% CO<sub>2</sub> and 37°C).
- Aspirate transduction solution and add 100 uL complete growth medium with sodium butyrate at a concentration of 1-2 mM. Return cells to normal growth conditions for approximately 24 hrs.

## Fluorescence Detection

The optimal excitation wavelength for imaging green fluorescent PIP<sub>2</sub> sensor is 480 nm, but the absorption band of this protein is quite broad, so broad bandpass filters that pass 450 to 480 nm light can be used effectively. On the emission side, the green light spans 510 to 550 nm, so broad band pass emission filters can also collect much of the emission. These filter properties are similar to many of the FITC filter sets commonly available on most microscopes and plate readers.




Figure 1. The absorption and emission properties of a typical green fluorescent protein are plotted above as a function of wavelength. Optimal excitation light ranges from 450 to 490 nm, while optimal emission filters should select the emitted light between 510 and 550 nm.

## Timing

The PIP<sub>2</sub> assay measures PIP<sub>2</sub> in living cells, in real time. Be sure to capture fluorescence during the peak response as shown in Figure 2.

## Assay Performance Considerations

### How we measure the infectivity of the viral stock.

Typically, viruses are quantified in terms of plaque forming units (PFU). In the case of BacMam, this would be a measurement of the viruses that are capable of transducing an insect cell, the natural host. Since mammalian cell expression is the goal for this assay, we quantify infectivity by measuring (VG) per milliliter (mL) of the BacMam stock. We use primers specific to the VSVG gene present in the BacMam genome. Viral samples are prepared to release viral genomic DNA, then multiple dilutions of the preparation are run in qPCR against a standard curve to generate an average titer for each viral stock. Check the label on the tube to find VG/mL for your PIP<sub>2</sub> sensor stock.

### Green PIP<sub>2</sub> + hMIRsanta 50 $\mu$ M Carbachol treatment





Figure 2. Green fluorescent PIP<sub>2</sub> sensor in HEK 293 cells, co-transduced with human muscarinic acetylcholine control receptor M1, and a red nuclear label. Fluorescence was captured from cells before the addition of 50  $\mu$ M carbachol, and sampled at regular intervals. The maximal response is reached at  $\sim$ 60 seconds after the addition of the drug, and the response begins its return to baseline  $\sim$ 100 seconds after drug is added.

## Expression levels of the sensor.

To optimize the assay in your particular cell type, it is important to optimize the amount of BacMam virus used in the transduction. Too little virus will produce variable results particularly if the sensor expression levels are low and difficult to detect on your instrument.

## Receptor expression

The magnitude of the sensor response can be affected by the level of GPCR expression in your cells. We have found that low levels of receptor expression produce the largest signals, while high levels of receptor expression often produce smaller responses. This is consistent with the observation that over expression of some receptors can artificially change the resting levels of second messengers.

## Trouble Shooting

Here are a few simple steps that may help you trouble shoot if needed.

### Are your cells fluorescent?

Different types of promoters drive expression in mammalian cells. The CMV promoter in our BacMam vectors is an effective promoter in many cell lines, but not all. If your cells are not fluorescent, check that the CMV promoter works in your cell line. Twenty four hours after transduction with the downward PIP<sub>2</sub> sensor, you should see green fluorescent cells in a typical epifluorescence microscope, or the transduced wells in a 96 well plate should be significantly more fluorescent than untransduced cells in wells on the same plate.

HDAC inhibitors are important to maintain expression of the sensors. While BacMam transduction without the HDAC inhibitor will initially generate low levels of sensor expression, it is important to follow protocols and include sodium butyrate or another HDAC inhibitor such as VPA, or trichostatin A (TSA) to generate optimal levels of sensor expression.

### Is the positive control working?

If the cells are expressing the sensor, and fluorescence is detectable on your instrument, then check the sensor using the positive control receptor included in this kit. Adding 5 uL of the hM1 receptor stock to set of control wells will ensure that a Gq coupled receptor is present in all of the cells. You can double check to make sure the M1 receptor is expressed by examining the cells in a fluorescent microscope with filters for red fluorescence. You should see red nuclear fluorescence that marks the cells that express the exogenous M1 receptor as shown in Figure 2.

Addition of carbachol will cause a change in fluorescence, when the receptor control is present in the cells, as shown in Figure 2. This positive control can be used to optimize three aspects of your assay. First, a serial dilution series of the sensor with a constant amount of PIP<sub>2</sub> sensor can be used to determine the optimal sensor expression for your instrument and cell type. Second, it is important to titrate the amount of BacMam sufficient to transduce all of the cells in the well. Third, it is important to determine whether your instrument can measure the peak response in the appropriate time frame.

## Contact Us

If you have any questions about the protocols described here, or if you have ideas about how we can improve these tools, then we hope to hear from you. Your feedback is extremely valuable. Please send an email to [info@montanamolecular.com](mailto:info@montanamolecular.com), and we'll respond as quickly as we can.

## References

1. Graham FL, Smiley J, Russell WC, Nairn R: Characteristics of a human cell line transformed by DNA from human adenovirus type 5. *J Gen Virol* 1977, 36(1):59-74.
2. Dulbecco R and Vogt M: Plaque formation and isolation of pure lines with poliomyelitis viruses. *The Journal of experimental medicine* 1954.
3. Chalfie M, Tu Y, Euskirchen G, Ward WW, Prasher DC: Green fluorescent protein as a marker for gene expression. *Science* 1994.
4. Kost T, Condreay J, Ames R, Rees S, Romanos M: Implementation of BacMam virus gene delivery technology in a drug discovery setting. *Drug Discovery Today* 2007, 12(9-10): 396-403.
5. Tewson P, Quinn AM, Hughes TE.: A Multiplexed Fluorescent Assay for Independent Second-Messenger Systems: Decoding GPCR Activation in Living Cells. *Journal of Biomolecular Screening* 2013 vol. 18 no. 7 797-806.