



## ELISA-VIDITEST anti-VCA EBV IgA

ODZ-096

### Instruction manual

**PRODUCER:** VIDIA spol.s r.o., Nad Safinou II 365, Vestec, 252 42 Jesenice, Czech Republic, Tel.: +420 261 090 565, www.vidia.cz, info@vidia.cz

#### 1. TITLE:

ELISA-VIDITEST anti-VCA EBV IgA – ELISA kit for the detection of IgA antibodies to Epstein-Barr virus (EBV) capsid antigen (VCA) in serum.

#### 2. INTENDED USE:

ELISA-VIDITEST anti-VCA EBV IgA is intended for in vitro diagnostic procedures in EBV induced and EBV associated diseases, such as infectious mononucleosis (IM) and the chronic active EBV infection. The test is also useful in the diagnosis of Burkitt's lymphoma, nasopharyngeal carcinoma, carcinoma of the Waldeyer's ring and in characterisation of opportunistic lymphomas (oligo- and polyclonal). The other use of the test is in characterisation of chronic fatigue syndrome, in neuroinfections and immunosuppression which is frequently associated with EBV reactivation. The determination of IgA antibodies to VCA is suitable mainly during EBV reactivation, when the levels of IgM antibodies to VCA are low and there are also present IgG antibodies to VCA and EBNA-1.

#### 3. TEST PRINCIPLE:

ELISA-VIDITEST anti-VCA EBV IgA assay is a solid-phase immunoanalytical test. The polystyrene strips are coated with specific antigen that bear immunodominant epitopes of VCA complex. The anti-VCA EBV antibodies, if present in the tested sera, bind to the immobilized antigens and the antibodies being in complexes with antigen are later on recognised by animal anti-human IgA antibodies labelled with horseradish peroxidase. The labelled antibodies are revealed by an enzymatic reaction with a chromogenic substrate. Negative sera do not react and the mild change in colour, if present, may be attributed to the reaction background.

#### 4. KIT COMPONENTS:

ELISA break-away strips coated with specific antigen	STRIPS Ag	1 microplate
1.3 mL Calibrator (Cal) r.t.u. <sup>1)</sup>	CAL	1 vial
1.3 mL Positive control serum (PCS) r.t.u.	CONTROL +	1 vial
1.3 mL Negative control serum (NCS) r.t.u.	CONTROL -	1 vial
15 mL Anti-Human IgA antibodies labelled with horseradish peroxidase r.t.u (Px-conjugate)	CONJ	1 vial
125 mL Wash buffer 10x conc.	WASH 10x	1 vial
125 mL Dilution buffer (DB) r.t.u.	DIL	1 vial
15 mL TMB substrate r.t.u.	TMB	1 vial
15 mL Stop solution r.t.u.	STOP	1 vial

Sealable pouch for unused strips  
Instruction manual  
Certificate of quality  
1) ready to use

## 5. MATERIAL REQUIRED BUT NOT PROVIDED WITH THE KIT:

1. Distilled or deionised water for dilution of the Wash buffer concentrate.
2. Appropriate equipment for pipetting, liquid dispensing and washing.
3. Spectrophotometer/colorimeter (microplate reader – wavelength 450 nm).

## 6. PREPARATION OF REAGENTS AND SAMPLES:

- a. Allow all kit components to reach room temperature.
- b. Vortex samples and the controls in order to ensure homogeneity and mix all solutions well prior use.
- c. **Dilute serum samples 1:100 in Dilution buffer** and mix (e.g. 5  $\mu$ L of serum sample + 500  $\mu$ L of Dilution buffer). Do not dilute the Control sera and Calibrator, they are ready to use.
- d. Prepare Wash buffer by diluting the Wash buffer concentrate 10 times with an appropriate volume of distilled or deionised water (e.g. 100 mL of the concentrated Wash buffer + 900 mL of distilled water). If there are crystals of salt present in the concentrated Wash buffer, warm up the vial to +32 to +37°C in a water bath. Diluted Wash buffer is stable for one week if stored at +2 to +10°C.
- e. Mix well Px-conjugate (CONJ) and TMB substrate r.t.u.
- f. Do not dilute TMB substrate and Stop solution and Px-conjugate they are ready to use.

## 7. ASSAY PROCEDURE:

- a. Allow the microwell strips sealed inside the aluminium bag to reach room temperature. Withdraw an adequate number of strips and put the unused strips into the provided pouch and seal it carefully with the desiccant kept inside.
- b. Start with filling the first well with 100  $\mu$ L of Dilution buffer (DIL) to estimate the reaction background. Fill next two wells with 100  $\mu$ L/well of Calibrator (CAL) and then fill next wells with Negative control serum (CONTROL -). It is also suitable to apply Positive control serum (CONTROL +) for the test control. Fill the remaining wells with 100  $\mu$ L of diluted serum samples (S1, S2, S3,...). It is sufficient to apply one serum into one well, however, if you wish to minimize laboratory error, apply the Calibrator in triplet and samples and Controls in doublets.
- c. Incubate **30 minutes** ( $\pm$ 2 min) at room temperature.
- d. Aspirate the liquid from wells into a waste bottle containing an appropriate disinfectant (see Safety Precautions). Wash and aspirate the wells four times with 250  $\mu$ L/well of Wash buffer. Avoid cross-contamination between wells! If any liquid stays trapped inside the wells, invert the plate and tap it on an adsorbent paper to remove the remaining drops.
- e. Add 100  $\mu$ L of Px-conjugate r.t.u. (CONJ) into each well. **Incubate 60 minutes ( $\pm$ 5 min) at room temperature.**

- f. Aspirate and wash 4x with 250  $\mu$ L/well of Wash buffer. Tap the plate on an adsorbent paper.
- g. Dispense 100  $\mu$ L of TMB substrate into each well. **Incubate for 10 minutes (+/-5 seconds) in dark at room temperature.** The time measurement must be started at the beginning of TMB dispensing. Cover the strips and keep them in the dark during the incubation with TMB substrate.
- h. **Stop the reaction** by adding 100  $\mu$ L of Stop solution. Use the same pipetting rhythm as with the TMB substrate to ensure the same reaction time in all wells. Tap gently the microplate few times to ensure complete mixing of the reagents.
- i. Measure the absorbance at 450 nm with a microplate reader within 10 minutes. It is recommended to use reference reading at 620-690 nm.

Figure 1: Pipetting scheme

	1	2	3	4	5	6	7	8	9	10	11	12
a	<b>DIL</b>	<b>S5</b>										
b	<b>CAL</b>	<b>S...</b>										
c	<b>CAL</b>											
d	<b>CONTROL -</b>											
e	<b>S1</b>											
f	<b>S2</b>											
g	<b>S3</b>											
h	<b>S4</b>											

## 8. PROCESSING OF RESULTS:

Begin the processing of results with subtraction of the background absorbance (absorbance of the DIL well) from the absorbances of all other wells.

### 8.1. Processing of results for the qualitative interpretation

1. Compute the absorbance mean of the wells with Calibrator (CAL). If the Calibrator was applied in three parallels and in one of them the absorbance is different from the mean in more than 20% then exclude the deviating well from the calculation and compute a new absorbance mean with using the other two wells.
2. **Compute the cut-off value** by multiplying the mean of Calibrator (CAL) with a correction factor. **The correction factor value for Calibrator for this Lot is written in enclosed Quality control certificate.**
3. Serum samples with absorbances lower than the 90% cut-off value are considered negative and samples with absorbances higher than the 110% of the cut-off value are considered positive.

## 8.2 Processing of results for the semiquantitative interpretation

Determine Positivity Index for each serum sample as follows:

1. Compute the cut-off value (see the previous paragraph)
2. Compute the Positivity Index according to the following formula:

$$\text{Sample Positivity Index} = \frac{\text{Sample absorbance}}{\text{Cut-off value}}$$

3. Express the serum reactivity according to Table 1 (Semiquantitative interpretation of results)

Table 1: Semiquantitative interpretation of result:

Positivity Index	Interpretation
< 0.9	Negative
0.90 – 1.10	+/-
1.11 - 4.00	+
4.01 – 6.50	++
6.51 – 12.00	+++
> 12.00	++++

*Note! An equivocal sample reactivity, i.e. interpreted as +/-, requires repetition of the sample testing. If the result is again indifferent then it is recommended to use an alternative testing method or to obtain another sample from the patient, usually withdrawn 1-2 weeks later.*

Example of calculation:

Calibrator absorbances = 0.986; 0.996; 0.998  
Mean Calibrator absorbance = 0.993  
Correction factor = 0.19  
Cut-off value = 0.993\*0.19 = 0.189  
Sample absorbance = 0.800  
Sample Positivity Index = 0.800 / 0.189 = 4.23

## 9. INTERPRETATION OF RESULTS:

IgA anti-VCA is considered as a marker of active EBV infection (either the primary infection or reactivation). IgA anti-VCA positivity, resulting from asymptomatic EBV reactivation, can be observed in some healthy seropositive individuals, too. Result of IgA anti-VCA antibody assay must be interpreted only in the context with patient's symptoms and with the results of the other complementary EBV serological examinations – see following Table:

Interpretation of the anti-EBV antibodies pattern in serum.

VCA IgG	VCA IgM	VCA IgA	EA (D) IgG	EBNA-1 IgG	EBNA-1 IgM	EBV status
-	-	-	-	-	-	EBV negative
-	+	+	+	-	+	EBV primoinfection (acute phase)
+	+	+	-	-	+	
+	+	+	+	+	+	EBV primoinfection (post-acute phase )
+	+	+	+	-	+	
+	-	-	-	-	-	
-	+	+	-	-	-	EBV primoinfection (convalescent phase)
-	+	-	-	+	-	
+	-	+	-	+	-	Seropositivity without signs of the active EBV infection
+	-	-	-	+	-	
+	-	-	+	+	-	EBV reactivation
+	+	+	+	+	-	

For the determination of the antibodies to other EBV antigens, VIDIA offers these kits:  
ELISA-VIDITEST anti-EBNA-1 EBV IgG and IgM; ELISA-VIDITEST anti-VCA EBV IgG and IgM; ELISA-VIDITEST anti-EA (D) EBV IgG.

As a confirmatory tests can be used IF-VIDITEST anti-VCA EBV IgG and IgM; IF-VIDITEST anti-EA (D+R or D) EBV IgG.

## 10. VALIDITY, SPECIFICITY AND SENSITIVITY OF THE TEST:

### 10.1 Validity of the test

The test is valid if:

The background of the reaction (the absorbance of the Dilution buffer) is less than 0.100.

The absorbance of Calibrator should be in range that is written in enclosed Quality control certificate.

The Controls keep the order that: CONTROL - < CAL (< CONTROL +).

### 10.2 Precision of the test

The intraassay variability (within the test) and the interassay variability (between tests) were performed with samples of variable absorbance values.

### 10.2.1 Intraassay variability

The coefficient of intraassay variability is max. 5%. It is measured for each particular Lot at least on 12 parallels of the same microtitration plate.

Example: (N = number of parallels of the same microtitration plate, SD = standard deviation)

N	Mean absorbance	SD	CV%
16	1.335	0.050	3.8 %
16	0.614	0.023	3.7 %

### 10.2.2 Interassay variability

The coefficient of interassay variability is max. 15%. It is measured for each particular Lot as comparison of the OD values of the same serum sample in several consecutive tests.

Example: (N = number of an independent examination of the same serum sample, SD = standard deviation):

N	Mean Absorbance	SD	Range (min-max)	CV%
18	1.369	0.064	1.223-1.476	4.7 %
18	0.463	0.060	0.337-0.569	12.9 %
14	1.128	0.093	0.945-1.319	8.2 %

### 10.2.3. Recovery test

Measured values of recovery test for every Lot are between 80-120% of expected values.

## 10.3 Diagnostic sensitivity and specificity

Diagnostic sensitivity and specificity was performed using the samples from patients with infectious mononucleosis, patients with serological proved EBV reactivation and healthy persons. The result was confirmed with other commercial test. The diagnostic sensitivity is 85%, in patients with IM diagnosis it is 92%. The diagnostic specificity is 100%.

## 10.4. Interference

Haemolytic and lipemic samples have no influence on the test results up to concentration of 50 mg/mL of haemoglobin, 5 mg/mL of bilirubin and 50 mg/mL of triglycerides.

## **11. SAFETY PRECAUTIONS:**

All ingredients of the kit are intended for laboratory use only.

Controls contain human sera that has been tested negative for HBsAg, anti-HIV-1,2 and anti-HCV. However they should be regarded as contagious and handled and disposed of according to the appropriate regulations.

Autoclave all reusable materials that were in contact with human samples for 1 hour at 121°C, burn disposable ignitable materials, decontaminate liquid wastes and nonignitable materials with 3% chloramine. Liquid wastes containing acid (Stop solution) should be neutralized in 4% sodium bicarbonate solution.

The waste from strip washing disinfect in waste container using appropriate disinfection solution (e.g. Incidur, Incidin, chloramin,...) in concentrations recommended by the producer.

Handle Stop solution with care. Avoid contact with skin or mucous membranes. In case of contact with skin, rinse immediately with plenty of water and seek medical advice.

Do not smoke, eat or drink during work. Do not pipette by mouth. Wear disposable gloves while handling reagents or samples and wash your hands thoroughly afterwards. Avoid spilling or producing aerosol.

## **12. HANDLING PRECAUTIONS:**

Manufacturer guarantees performance of the entire ELISA kit.

Follow the assay procedure indicated in the Instruction manual.

Wash buffer, TMB substrate, Stop solution and Dilution buffer are interchangeable between ELISA-VIDITEST kits unless otherwise stated in the instruction manual.

Calibrator, Control sera, TMB substrate, Dilution buffer and Px-conjugate contain preservative ProClin 300<sup>®</sup>.

Avoid microbial contamination of serum samples and kit reagents.

Avoid cross-contamination of reagents.

Avoid contact of the TMB substrate with oxidizing agents or metal surfaces.

Variations in the test results are usually due to:

- Insufficient mixing of reagents and samples
- Inaccurate pipetting and inadequate incubation times
- Poor washing technique or spilling the rim of well with sample or Px-conjugate
- Use of identical pipette tip for different solutions

## **13. STORAGE AND EXPIRATION:**

Store the kit and the kit reagents at +2 to +10°C, in a dry place and protected from the light. Store unused strips in the sealable pouch and keep the desiccant inside.


Store serum samples at +2 to +10°C up to one week. For longer period make aliquots and keep them at -18°C to -28°C. Avoid repeated thawing and freezing.

Do not store diluted serum samples. Always prepare fresh.

Kits are shipped in cooling bags, the transport time of 72 hours have no influence on expiration. If you find damage at any part of the kit, please inform the manufacturer immediately.

Expiration date is indicated at the ELISA kit label and at all reagent labels.

#### 14. FLOW CHART:



#### References:

Wakiguchi H., Fujieda M., Kubota H., Matsumoto K., Wakiguchi A., Kurashige T., Oda M.: Elevated IgA antibodies to Epstein-Barr virus in children with chronic active Epstein-Barr virus infection; *Acta Med Okayama*. 1989 Jun; 43 (3): 193-6

Marklund G., Henle W., Henle G., Ernberg I.: IgA antibodies to Epstein-Barr virus in infectious mononucleosis; *Scand J Infect Dis*. 1986; 18(2): 111-9

Roubalova K., Roubal J., Skopovy P., Fucikova T., Domorazkova E., Vonka V.: Antibody response to Epstein-Barr virus antigens in patients with chronic viral infection; *J Med Virol*. 1988 May; 25(1): 115-22

Chalupa P., Slesinger P., Krchnakova A., Prasek J.: IgA antibodies against Epstein-Barr virus capsid antigen; *Cas Lek Cesk.* 1992 Sep 10; 131(17): 526-9

Nikoskelainen J., Neel EU, Stevens DA: Epstein-Barr virus-specific serum immunoglobulin A as an acute-phase antibody in infectious mononucleosis; *J Clin microbiol.* 1979 Jul; 10(1): 75-9

Edwards JM, Woodroof M.: EB virus specific IgA in serum of patients with infectious mononucleosis and of healthy people of different ages; *J Clin Pathol.* 1979 Oct; 32(10): 1036-40

Hadar T., Rahima M., Kahan E., Sidi J., Rakowsky E., Sarov B., Sarov I.: Significance of specific Epstein-Barr virus IgA and elevated IgG antibodies to viral capsid antigens in nasopharyngeal carcinoma patients; PMID: 3025351 [PubMed]

Date of the last revision of this manual: 11/2012