



# PRODUCT INFORMATION SHEET

## IMMOBILIZED PROTEINASE K G3M

# P3502

Proteinase K from *Tritirachium album*. Chromatographically purified, free of ribo- and deoxyribonucleases.

Proteinase K is an unspecific serine protease with strong proteolytic activity on denatured (in SDS) and high molecular weight native proteins. It cleaves peptide bonds mostly after the carboxyl group of N-substituted hydrophobic, aliphatic and aromatic amino acids.

**G3m:** 25 µg proteinase K immobilized on matrix G3m per CR-column.

0.7 mAnson units immobilized per CR-column.

This CR-column cuts at least 370 µg BSA per application.

**Nr. 5 Storage buffer:** 50 mM Tris/HCl, pH 7.5

**Nr. 16 Reaction buffer:** 50 mM Tris/HCl, 5 mM NaCl, pH 8.0

**Nr. 17 Washing buffer:** 50 mM Tris/HCl, 1.0 M NaCl, pH 8.0

The columns are more active in 0.1% SDS and at 40°C. Also active in PBS buffer (20 mM Na-phosphate, 150 mM NaCl at pH 7.6).

## Protocol

For more details see MoBiTec-CRC-Handbook.

### 1. Dilute delivered buffers (at least 2 ml each) with sterile doubly distilled water.

For 1 application you need:

1 ml 10x reaction buffer and 9 ml doubly distilled water

2 ml 5x washing buffer and 8 ml doubly distilled water

1 ml 10x storage buffer and 9 ml doubly distilled water

The substrate should be in reaction buffer

### 2. Equilibrate the CR-column with 10 ml reaction buffer.

Fill 10 ml reaction buffer into a syringe, let the reaction buffer run through the column by gravity to the upper filter. In case the buffer runs very slowly, apply pressure by a syringe.

### 3. Load substrate solution in reaction buffer.

Small volumes (< 70 µl): spin the CR-column 5 seconds in a benchtop centrifuge (2000 rpm are sufficient). Let the substrate solution enter the matrix material.

Larger volumes: Let the substrate solution run through the column.

Flow-rate: up to 70 µl/minute

Keep the substrate in the column for about 1 minute at room temperature. Higher turn-over is obtained when the substrate is applied to the column again or incubated for longer times.

### 4. Elute the product solution.

Small volumes (< 70 µl): centrifuge the product out of the column.

Larger volumes: Let the substrate run through the column and spin the residual solution out of the matrix

Notice: Molecules < 700 Dalton have to be eluted with 7 ml reaction buffer.

It does not harm the columns if they run dry.

### 5. Wash the column with 10 ml washing buffer.

### 6. Equilibrate the column with 10 ml storage buffer.

Store the column at 4°C. Never freeze a CR-column!