

hp-Vector Expression System  
for *Bacillus megaterium*\*



**Mo Bi Tec**  
MOLECULAR BIOTECHNOLOGY



\*Reclassified as *Priestia megaterium*, Gupta *et al.* 2020



## Content

1. Features .....	4
1.1. Recombinant protein production with <i>Bacillus megaterium</i> .....	4
1.2. High performance (hp) vectors for xylose-inducible gene expression .....	4
2. Introduction.....	4
2.1. The hp-Vector Expression System .....	5
3. Expression Vectors and Strains .....	6
4. Protocols .....	8
4.1. Cloning the DNA fragment of interest.....	8
4.2. General remarks on the handling of <i>B. megaterium</i> .....	8
4.3. Transformation of <i>B. megaterium</i> protoplasts .....	8
4.4. Recombinant protein production and secretion .....	9
4.5. Analysis of proteins from cell extract (intracellular and membrane proteins) .....	10
4.6. Precipitation of proteins in the cell free supernatant.....	10
4.7. Scale up protein production.....	11
5. Materials.....	11
6. References .....	13
7. Vector Maps .....	15
7.1. Vector map of p3STOP1623hp .....	15
7.2. Vector map of pC-HIS1623hp (#BMEG31) .....	16
7.3. Vector map of pN-HIS-TEV1623hp .....	17
7.4. Vector map of pSP <sub>LipA</sub> -hp .....	18
7.5. Vector map of pSP <sub>YocH</sub> -hp .....	19
7.6. Vector map of p3STOP1623-2RBSHp.....	20
7.7. Vector map of pC-STREP1623hp .....	21
7.8. Vector map of pN-STREP-Xa1623hp.....	22
7.9. Vector map of pN-STREP-TEV1623hp .....	23
7.10. Vector map of pMGBM19 .....	24
7.11. Vector map of helper plasmid .....	25
7.12. Map of the control vector pGFP1624hp .....	26
8. Order Information, Shipping, and Storage .....	27
9. Contact and Support .....	28



## 1. Features

### 1.1. Recombinant protein production with *Bacillus megaterium*

- Stable, high-yield production of recombinant proteins using *Bacillus megaterium*
- Growth rate comparable to that of *Escherichia coli* but reaching higher final optical density (OD)
- No alkaline proteases activity even up to 5 hours after induction
- No endotoxins
- MoBiTec host strains are asporogenic on common media
- *B. megaterium* is not pathogenic

### 1.2. High performance (hp) vectors for xylose-inducible gene expression

- Tightly regulated and efficiently xylose-inducible promoter system
- Optimized sequence of promoter and ribosome binding site (RBS)
- Up to tenfold enhanced protein yields in comparison to basic expression system
- hp vectors for intracellular or extracellular protein production
- hp vectors encoding for C- or N-terminal His<sub>6</sub>- or Strep-tag for easy purification of recombinant proteins.
- Removable N-terminal tag versions (TEV or Xa protease cleavage site included)
- 2RBS vector for simultaneous dual expression
- Convenient cloning due to *B. megaterium* / *E. coli* shuttle vectors with large MCS
- System suitable also for (other) *Bacillus* ssp. (Lakowitz *et al.* 2017)

## 2. Introduction

**Preliminary note:** Due to recent findings (Gupta *et al.* 2020) *Bacillus megaterium* was reclassified as *Priestia megaterium* and belongs now to the novel genera *Priestia*. However, this handbook will stick for the present to the former and more familiar name: *Bacillus megaterium*.

*B. megaterium* is a rod-shaped, non-pathogenic Gram-positive bacterium that is able to grow on a wide variety of carbon sources. Like other soil bacteria, it is able to secrete high amounts of proteins into the culture medium (up to 0.5 g/L of recombinant protein). This property, together with its ability to grow in many ecological niches such as waste from meat industry, molasses or petrochemical effluents makes it an interesting organism for the industrial large-scale protein production. Further, it was studied (next to *B. subtilis*) as Gram-positive model organism regarding biochemistry, sporulation and bacteriophages.

The name “megaterium” is due to its large size of vegetative cells (>10 µm). Cells often occur in pairs of chains, joined together by polysaccharides on the cell walls. An overview about the features of this unique organism is given in several review articles listed in chapter 6.


One of the genetic regulatory elements for carbon metabolism of *B. megaterium* is the xylose operon that has been described by Rygus and Hillen (1991). It is used in the vector-based and easy-to-handle expression system described here. Within this system, MoBiTec




is offering diverse xylose inducible *E. coli* / *B. megaterium* shuttle vectors and three different *B. megaterium* strains, provided as protoplasts, ready for transformation. All vectors are multi-copy plasmids for cloning in *E. coli* and xylose-induced gene expression in *B. megaterium*. Different signals for protein secretion and affinity tags for purification (including removable tag versions) arranged in a high variety of combinations are available. By using this system, recombinant proteins can be produced, secreted and purified in the g/L scale.

## 2.1. The hp-Vector Expression System

With our improved hp- (high performance) vectors the yields of recombinant proteins can be more than tenfold enhanced in comparison to the basic plasmids carrying the native promoter system. All plasmids have an established multiple cloning site (MCS) for versatile cloning and an optimized promoter and ribosome binding site (RBS). We offer vectors encoding C- or N-terminal His<sub>6</sub>- or Strep-tags for easy purification, a tag-free plasmid, and a vector carrying two ribosomal binding sites (2RBS) for simultaneous dual expression. The protein secretion with signal peptides of the lipase A (LipA) or of the YochH (putative exported cell wall-binding protein) is increased up to ninefold, while the secreted protein can also be fused to a C-terminal His<sub>6</sub>-tag. Induction of protein production of all vectors is achieved by the tightly regulated and efficiently inducible xylose promoter.



**Fig. 1.** Soluble protein fractions 6 h after induction of heterologous gene expression



**Fig. 2.** Relative fluorescence mediated by the optimized hp-vector (red; square) and the basic plasmid (blue; circle), respectively, over time [h] after induction



### 3. Expression Vectors and Strains

#### Vectors of the 1623hp series

All vectors of the 1623hp series were constructed from the basic vector pSTOP1622 (referring to the manual on *Bacillus megaterium* Protein Production System) which is used as a shuttle vector for xylose-inducible production of recombinant target proteins in *Bacillus megaterium*. The pSTOP1622 contains one stop codon downstream of the multiple cloning site (MCS). The plasmid p3STOP1623 has been created by introducing two additional stop codons downstream of the MCS into the pSTOP1622 allowing cloning in all three reading frames.

The high performance (hp)-expression vectors as indicated below were generated by optimizing the promoter  $P_{xyIA}$  (optimized -35 box) and ribosomal binding site (optimized RBS) of the basic vector p3STOP1623.

#### **p3STOP1623hp**

This high-performance plasmid contains two additional stop codons downstream of the MSC covering all open reading frames, an optimized -35 region of  $P_{xyIA}$  and an optimized RBS. All following plasmids are based on this plasmid.

#### **pC-HIS1623hp**

This plasmid allows a C-terminal 6xHis-tag (3' of MCS) fusion including a stop codon downstream of the tag.

#### **pN-HIS-TEV1623hp**

The plasmid encodes an N-terminal 6xHis-tag fusion including a TEV (tobacco etch virus) recognition site that is cleavable with TEV protease treatment.

#### **pSP<sub>LipA</sub>-hp**

This plasmid contains the coding sequence of signal peptide of *B. megaterium* protein LipA for secretion and it encodes for a His<sub>6</sub>-tag for C-terminal fusion.

#### **pSP<sub>Yoch</sub>-hp**

The plasmid encodes the signal peptide of *B. megaterium* protein Yoch which enables secretion of proteins fused to a C-terminal His<sub>6</sub>-tag.

#### **p3STOP1623-2RBShp**

This plasmid carries an additional optimized ribosomal binding site. It is suitable for simultaneous dual expression.

#### **pC-STREP1623hp**

This plasmid allows a C-terminal StrepII-tag fusion to the protein of interest.

#### **pN-STREP-Xa1623hp**

pN-StrepXa1623hp encodes an N-terminal StrepII-tag fusion followed by a Factor Xa recognition site which is cleavable with Factor Xa protease.

#### **pN-STREP-TEV1623hp**

This plasmid carries the coding region of an N-terminal StrepII-tag fusion with TEV (tobacco etch virus) recognition site which is cleavable with TEV protease.



## Vectors for special requirements

### **pMGBm19**

pMGBm19 is an *E. coli*/*Bacillus megaterium* shuttle vector with xylose-inducible P<sub>xyIA</sub> promoter that is designed for co-expression studies. It can be used in combination with any other vector of the 1520, 1622, and 1623hp series, since it contains an origin of replication of a different compatibility group (*repM100* replicon).

### **pMMEc4**

Since the xylose-inducible P<sub>xyIA</sub> promoter is not tightly controlled in *E. coli*, cloning toxic genes into vectors of the 1520, 1622, and 1623hp series, respectively, using *E. coli* as host may be difficult. In such cases, we recommend using the pMMEc4 helper plasmid. This *E. coli* vector (not replicating in *Bacillus*!) encodes for the xylose repressor XylR and is designed for blocking any expression starting from the P<sub>xyIA</sub> promoter while cloning the gene of interest within *E. coli* (Jordan et al., 2007).

In pMMEc4 the expression of *xyIR* is controlled by the arabinose-dependent promoter P<sub>BAD</sub> and the AraC protein. In the presence of 0.2% arabinose, the AraC protein binds to the operator sequence that activates the expression of the *xyIR* gene, and additionally upregulates its own expression. The vector pMMEc4 carries the p15A origin of replication that is compatible with vectors from other incompatibility groups, like ColE1 or the origin of pBR322.

### **pGFP1624hp**

pGFP1624hp is a positive control vector expressing *gfp* under control of the xylose-inducible P<sub>xyIA</sub> promoter. This vector belongs to the high-performance expression vectors series of *B. megaterium* (1623hp series) with an optimized promoter region that leads to improved protein yield. pGFP1624hp is an *E. coli*/*B. megaterium* shuttle vector. The vector can be propagated within *E. coli* by selection with 100 µg/ml ampicillin. For maintenance of the vector in *B. megaterium* use 10 µg/ml of tetracycline.

## **Bacillus megaterium strains**

**1) The strain WH320 (#BMEG02)** is a chemical mutant of strain DSM319 which is deficient in the production of β-galactosidase ( $\Delta lacZ$ ). It was described by Rygus and Hillen (Rygus and Hillen, 1991).

**2) The strain MS941 (#BMEG50)** was generated from the wild-type strain DSM319 by deletion of the gene *nprM* encoding the major extracellular protease (Wittchen and Meinhardt 1995). Because of reduced extracellular protease activity this strain is well suited for extracellular protein production.

**3) The strain YYBm1 (#BMEG04)** carries the *nprM* deletion (like strain MS941) and an additional deletion of the xylose isomerase gene *xyIA*. It is thus unable to metabolize xylose which is used as inducer for gene activation (Yang et al. 2006).



## 4. Protocols

### 4.1. Cloning the DNA fragment of interest

The *E. coli* / *B. megaterium* shuttle vectors are supplied as lyophilized DNA. Green and Sambrook (2012) describe all standard protocols for propagation of the plasmid in *E. coli*, plasmid DNA preparation, restriction endonuclease cleavages, and ligation of the DNA fragment of interest into the vector.

After ligation of the insert the vectors should be propagated in *E. coli* (use 100 µg/ml ampicillin for selection) before transforming *B. megaterium* protoplasts (use 10 µg/ml tetracycline for selection). For all cloning purposes and propagation of plasmids we recommend to use *E. coli* strain DH10B.

### 4.2. General remarks on the handling of *B. megaterium*

*B. megaterium* strains grow well on rich media such as Luria-Bertani (LB) broth at 30 °C and 37 °C. Make sure to aerate liquid cultures well in baffled flasks by vigorous agitation.

We found MS941, WH320, YYBm1, and derived strains to be asporogenic on common growth media - they will die on the plates stored at 4 °C within two weeks, so prepare glycerol stocks (30% w/v) as a backup and streak the working cultures on fresh plates every 7 to 10 days.

Positive clones carrying the plasmid with a gene of interest can be selected by adding 10 µg/ml of tetracycline to the growth medium.

To prove successful overexpression of the target gene, harvest small samples of the culture just before, and at certain time intervals after induction of recombinant gene expression with xylose. To obtain crude extracts for protein analysis, cells have to be lysed using lysozyme. Simple boiling of cells in sample buffer (Laemmli, 1970), which is quite convenient for *E. coli*, **does not** lyse *B. megaterium* cells.

### 4.3. Transformation of *B. megaterium* protoplasts

For recombinant protein production, *B. megaterium* protoplasts are transformed with the plasmids coding for the protein of interest. After transformation it is advisable to screen at least three different clones for protein production as the yield can vary among different clones.

Since intact *B. megaterium* cells cannot be transformed, MoBiTec conveniently provides protoplasts of *B. megaterium* cells, which are ready for transformation (strains MS941, WH320, and YYBm1). They can be used at least 2 months after date of arrival and have to be stored at -80 °C. The protoplast suspension is supplied in 5 aliquots of 500 µl each to prevent multiple freezing and thawing of protoplasts that are not used immediately. **One aliquot is provided per transformation.** It is advisable to use two of the vials for the control experiments as described below.

#### Control experiments:

**Note:** Each lot of protoplasts undergoes this test during our quality control as well.

#### Negative control: protoplasts without DNA

This is a test reassuring that the protoplasts are not only fully viable but also free of contaminations before using them for transformation. Perform this test according to the transformation protocol as demonstrated below but without plasmid DNA donation. After incubation in SMMP (see below, step 10), apply CR5-top agar to the protoplasts and split the sample in two portions. You may plate one sample on a LB plate with antibiotic such





as tetracycline, and another one on a plate without any drug. In this case, bacterial colonies will grow only on a solid medium without antibiotic.

**Positive control: protoplasts transformed with an empty plasmid (not included)**

This is a control for a successful transformation and should yield lots of colonies on the plates supplemented with an antibiotic (here: tetracycline). If this transformation works well, but you have problems with the plasmid containing your target gene, the problem is most likely associated with your construct.

Essential buffers (listed in chapter 5) should be kept at room temperature (RT) at least 30 minutes before using them!

**Transformation procedure:**

1. For each transformation thaw one aliquot of protoplast (500  $\mu$ l) on ice.
2. Mix 3-5  $\mu$ g of plasmid DNA (dissolved in H<sub>2</sub>O<sub>dest.</sub>) with 10 $\mu$ l SMMP and incubate for 20 min. at 37°C.
3. Combine the plasmid DNA/SMMP mix with one aliquot of thawed protoplasts in a 15 ml tube.
4. Add 1.5 ml of PEG-P, mix gently, and incubate for 2 minutes at RT.
5. Add 5 ml SMMP and mix carefully by rolling the tube.
6. Harvest cells by gentle centrifugation (1,300 x g for 10 minutes at RT), discard the supernatant immediately after centrifugation. Supernatant does not have to be removed completely.  
*(Note: do not check for a cell precipitate - most of the time it will be invisible)*
7. Add 500  $\mu$ l of SMMP to remaining supernatant (containing bacterial cells) and transfer to a 1.5 ml microcentrifuge tube.
8. Incubate at 30 °C for 45 min without shaking followed by another 45 min shaking at 300 rpm at 37°C in a thermomixer.
9. Prepare 2.5 ml aliquots of CR5-top agar in sterile tubes.
10. After the incubation step in SMMP add all cells to 2.5 ml top agar, mix gently by rolling the tube between both hands (do not vortex!), and pour onto a pre-warmed plate of LB containing desired antibiotic.
11. Incubate overnight (up to 24 h) at 30 °C - expect colonies of varying diameter because some will be covered with agar and others have easier access to air.  
*(Note: bacterial colonies on the top of the agar surface will be shiny)*
12. Transfer several single colonies (at least 3) on fresh plates within two days.

**Note:** *protein production may vary among the single colonies due to yet unknown reasons.*

**4.4. Recombinant protein production and secretion**

1. Grow the recombinant *B. megaterium* cells in 50 mL suited medium (e.g., LB medium) including tetracycline (10 mg/ml) in 300 ml baffled flasks to an optical density (OD<sub>578nm</sub>) of 0.3-0.4 at 37 °C under vigorous shaking at 180-250 rpm  
The main culture should be set up with an overnight culture in a dilution of 1:100.
2. Take a sample as control before induction.
3. Induce the xylose-inducible promoter by addition of 0.5% (w/v) of D-xylose.
4. Incubate at 37 °C and shake vigorously (180-250 rpm).
5. Withdraw samples every 2 hours for 6 to 8 hours (now, cells should have entered the stationary phase). An additional sample should be taken after around 24 h. Take samples for the OD<sub>578nm</sub> measurement and protein analysis. For extracellular protein



analysis take 2 ml of cell culture. Intracellular protein analysis requires a smaller volume (0.5 ml).

6. Centrifuge each sample at 4°C and 13,000 rpm for 10 min to separate cells and cell-free supernatant.
7. For extracellular protein analysis retain the supernatant and store at 4 °C, for intracellular protein analysis completely discard the supernatant and store the cells at -20 °C.

#### **4.5. Analysis of proteins from cell extract (intracellular and membrane proteins)**

1. Resuspend cells in 60 µl of lysis buffer.
2. Incubate at 37 °C for 30 min with vigorous shaking at 1,000 rpm in a thermomixer. An effective cell lysis can be obtained by whirling the samples every 10 minutes.
3. Add 60 µl of 2x Laemmli buffer
4. At this step, you may store your sample at -20°C
5. Heat the samples at 95 °C for 10 min and load an appropriate volume on SDS-PAGE. The amount of loaded sample should be calculated according to the OD<sub>578nm</sub> of the main culture. This allows to compare protein amounts of similar cell counts.
6. Perform Western blot using appropriate antibodies (not included in the kit).

#### **4.6. Precipitation of proteins in the cell free supernatant**

##### **Option 1 (non-denaturing): Ammonium sulfate precipitation**

1. Add 600 mg of ground ammonium sulfate to 1.5 ml of cell-free supernatant and incubate for two hours at 4 °C while shaking.
2. Centrifuge at 13,000 rpm and 4 °C for 30 minutes.
3. Completely remove the supernatant, centrifuge again for 1 min, and make sure the protein precipitate is dry.
4. Add 10 µl of 8 M urea (in 50 mM Tris-HCl, pH 7.5) and 10 µl 2x Laemmli buffer to dissolve proteins again.
5. Spin briefly at 13,000 rpm, heat at 95°C for 10 minutes, and load an appropriate volume onto SDS-PAGE for analysis. The amount of loaded sample should be calculated according to the OD<sub>578nm</sub> of the main culture. This allows to compare protein amounts of similar cell counts.
6. Perform Western blot using appropriate antibodies (not included in the kit).

##### **Option 2 (denaturing): TCA precipitation**

1. Add 200 µl of 100% TCA solution to 1,8 ml of cell-free supernatant and incubate on ice overnight
2. Centrifuge at 13,000 rpm and 4°C for 15 min.
3. Wash the pellet with 1 ml acetone, centrifuge at 13,000 rpm and 4°C for 15 min, discard the supernatant.
4. Dry the precipitate at 60°C.
5. Dissolve the proteins in 100 µl of 2x Laemmli buffer
5. At this step, you may store your sample at -20°C
6. Heat the samples at 95 °C for 10 min and load an appropriate volume on SDS-PAGE. The amount of loaded sample should be calculated according to the OD<sub>578nm</sub> of the main culture. This allows to compare protein amounts of similar cell counts.
7. Perform Western blot using appropriate antibodies (not included in the kit).



#### 4.7. Scale up protein production

1. Grow larger culture and induce as indicated above.
2. Harvest cells at the time point of maximal protein overproduction, as determined by the test experiments.

#### 5. Materials

##### 2 × AB3 (Antibiotic Medium No. 3, DIFCO)

- 7 g AB3 (Difco) in 200 ml deionized water
- Autoclave for 20 min

##### 2 × SMM

**Solubilize in the given order!**

- 1.16 g maleic acid (40 mM)
- 800 mg NaOH (80 mM)
- 2.03 g MgCl<sub>2</sub> × 6H<sub>2</sub>O (40 mM)
- 85.58 g sucrose (1 M)
- Solubilize each component in deion. water
- mix and fill with deion. water to 250 ml
- sterilize by filtration

##### SMMP

- 2 × AB3 and 2 × SMM 1:1 (freshly prepared!)

##### PEG-P

- solubilize 20 g PEG-6000 with 1 × SMM and fill to 50 ml
- autoclave for 15 min

##### CR5-top agar

Prepare separately for 500 ml:

###### Solution A

- 51.5 g sucrose
- 3.25 g MOPS
- 300mg NaOH
- add deionized water to 250 ml
- sterilize by filtration

###### Solution B

- 2 g agar
- 100 mg casamino acids
- 5 g yeast extract
- add deionized water to 142.5 ml
- autoclave for 20 min

###### 12% proline

- 3 g L-proline
- add deionized water to 25 ml
- sterilize by filtration

**8 × CR5-salts**

- 1.25 g K<sub>2</sub>SO<sub>4</sub>
- 50 g MgCl<sub>2</sub> × 6 H<sub>2</sub>O
- 250mg KH<sub>2</sub>PO<sub>4</sub>
- 11 g CaCl<sub>2</sub>
- solubilize in 625 ml deion. water
- autoclave for 15 min

**20% glucose**

- 5 g glucose
- add deionized water to 25 ml
- sterilize by filtration or autoclave

**for a 2.5 ml portion of CR5-top agar add the following (in the given order!):**

- 1.25 ml Solution A
- 288 µl CR5-salts
- 125 µl 12% proline
- 125 µl 20% glucose

**90 minutes after transformation:**

- Boil Solution B
- add 713 µl to the provided CR5-top agar
- immediately add the regenerated protoplasts and put onto prewarmed agar plates containing the corresponding antibiotic (here: tetracycline)

**Lysis buffer**

- 100 mM Natrium-phosphate-buffer (pH 6.5-7.0)
- 2 mM MgCl
- 5 mg/ml lysozyme
- 2 µl HS-Nuclease (5 U/µl, cat. #GE-NUC10700-01) per ml lysis buffer shortly before use

**2x Laemmli buffer**

- 4 ml of 10% SDS (final 4%)
- 2 ml 100% glycerol (final 20%)
- 1.25 ml 1M Tris-Cl, pH 6.8 (final 125 mM)
- 1ml β-mercaptoethanol (final 10%)
- 1.75 ml H<sub>2</sub>O<sub>deion.</sub>
- Add bromphenol blue to a final concentration of 0.02% (w/v)



## 6. References

Biedendieck, R., Yang, Y., Deckwer, W., Malten, M. & Jahn, D. (2007). Plasmid system for the intracellular production and purification of affinity-tagged proteins in *Bacillus megaterium*. *BiotechnolBioeng*96, 525-537.

Jordan E, Hust M, Roth A, Biedendieck R, Schirrmann T, Jahn D, Dübel S. (2007). Production of recombinant antibody fragments in *Bacillus megaterium*. *Microb Cell Fact*. Jan 15;6:2.

Lakowitz, A., Krull, R. & Biedendieck, R. (2017) Recombinant production of the antibody fragment D1.3 scFv with different *Bacillus* strains. *Microb Cell Fact*. 16, 14

Malten, M., Biedendieck, R., Gamer, M., Drews, A., Stammen, S., Buchholz, K., Dijkhuizen, L. & Jahn, D. (2006). A *Bacillus megaterium* plasmid system for the production, export, and one-step purification of affinity-tagged heterologous levansucrase from growth medium. *Appl Environ Microbiol*72, 1677-1679.

Rygus, T. & Hillen, W. (1991). Inducible high-level expression of heterologous genes in *Bacillus megaterium* using the regulatory elements of the xylose-utilization operon. *ApplMicrobiolBiotechnol*35, 594-599.

Stammen, S., Müller, B.K., Korneli, C., Biedendieck, R., Gamer, M., Franco-Lara, E. & Jahn, D. (2010) High yield intra- and extracellular protein production using *Bacillus megaterium*. *Appl Environ Microbiol*76, 4037-4046.

Wittchen, K. & Meinhardt, F. (1995). Inactivation of the major extracellular protease from *Bacillus megaterium* DSM319 by gene replacement. *ApplMicrobiolBiotechnol*42, 871-877.

Yang, Y., Malten, M., Biedendieck, R. et al. (2006). *Bacillus megaterium* as a recombinant protein production host. *Microb Cell Fact* 5, P74

## Reviews

Biedendieck, R., Bunk, B., Fürch, T., Franco-Lara, E., Jahn, M. & Jahn, D. (2010) Systems biology of recombinant protein production in *Bacillus megaterium* in "Advances in Biochemical Engineering / Biotechnology" (Scheper, T., Ed.); *Biosystems Engineering* (Wittmann, C., Ed.) Springer, Berlin, Heidelberg 120, 133-161

Bunk, B., Biedendieck, R., Jahn, D. & Vary, P.S. (2010) *Bacillus megaterium* and other Bacilli: Industrial Applications. In Flickinger M. C. (Ed.) *Encyclopedia of Industrial Biotechnology: Bioprocess, Bioseparation, and Cell Technology*. Volume 1 pp 429-443. John Wiley & Sons. Inc., Hoboken, NJ.

Bunk, B., Schulz, A., Stammen, S., Münch, R., Warren, M. J., Rohde, M., Jahn, D. & Biedendieck, R. (2010). A short story about a big magic bug. *Bioeng Bugs* 1:2, 1-7.

Eppinger, M., Bunk, B., Johns, M.A., Edirisinghe, J.N., Kutumbaka, K.K., Koenig, S.S.K., Creasy, H.H., Rosovitz, M.J., Riley, D.R., Daugherty, S., Martin, M., Elbourne, L.D.H., Paulsen, I., Biedendieck, R., Braun, C., Grayburn, S., Dhingra, S., Lukyanchuk, V., Ball, B., ul-Qamar, R., Seibel, J., Bremer, E., Jahn, D., Ravel, J. & Vary, P.S. (2011) Genome



sequences of the biotechnologically important *B. megaterium* strains QM B1551 and DSM319. *JBacteriol* 193, 4199-4213.

Vary, P. (1992). Development of genetic engineering in *Bacillus megaterium*. *Biotechnology* 22, 251-310.

Vary, P. S. (1994). Prime time for *Bacillus megaterium*. *Microbiology (Reading, Engl.)* 140 (Pt 5), 1001-1013.

Vary, P. S., Biedendieck, R., Fuerch, T., Meinhardt, F., Rohde, M., Deckwer, W.-D. & Jahn, D. (2007) *Bacillus megaterium* – from simple soil bacterium to industrial protein production host. *ApplMicrobiolBiotechnol* 76, 957-967.

## Methods and protocols

Biedendieck, R., Borgmeier, C., Bunk, B., Stammen, S., Scherling, C., Meinhardt, F., Wittmann, C. & Jahn, D. (2011) Systems biology of recombinant protein production using *Bacillus megaterium*. *Methods Enzymol* 500, 165-195.

Green, M. R. & Sambrook, J. (2012). *Molecular cloning: a laboratory manual*. Cold Spring Harbor, N.Y., Cold Spring Harbor Laboratory Press. 4. Edition, 3 vol. (1890 pages).

Laemmli, U. K. (1970). Cleavage of structural proteins during the assembly of the head of bacteriophage T4. *Nature* 227, 680-685.

Sambrook, J. & Russell, D.W. (2001). *Molecular cloning: a laboratory manual*. (New York: Cold Spring Harbor Laboratory Press).

## Further literature

Biedendieck, R., Beine, R., Gamer, M., Jordan, E., Buchholz, K., Seibel, J., Dijkhuizen, L., Malten, M. & Jahn, D. (2007). Export, purification, and activities of affinity tagged *Lactobacillus reuteri* levansucrase produced by *Bacillus megaterium*. *ApplMicrobiolBiotechnol* 74, 1062-1073.

Gupta RS, Patel S, Saini N, Chen S. (2020) Robust demarcation of 17 distinct *Bacillus* species clades, proposed as novel *Bacillaceae* genera, by phylogenomics and comparative genomic analyses: description of *Robertmurraya kyonggiensis* sp. nov. and proposal for an emended genus *Bacillus* limiting it only to the members of the *Subtilis* and *Cereus* clades of species. *Int J Syst Evol Microbiol*, 70:5753-5798.


Rygas, T., Scheler, A., Allmansberger, R. & Hillen, W. (1991). Molecular cloning, structure, promoters and regulatory elements for transcription of the *Bacillus megaterium* encoded regulon for xylose utilization. *ArchMicrobiol* 155, 535-542.



## 7. Vector Maps

All vector maps and the DNA sequences are available for download on our web page at <http://www.mobitec.com>.


### 7.1. Vector map of p3STOP1623hp (#BMEG30)



	Type	Start	End	Name	Description
	Region	1	10	opt. RBS	Optimized ribosomal binding site
	MCS	27	80	MCS	Multiple Cloning Site
	Origin of replication	175	837	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group
	Selectable genetic marker	1859	1001	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	3413	2130	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	4707	3706	repU	Gene encoding replication protein RepU
	Origin of replication	4859	4885	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	6285	5164	xylR	Gene encoding xylose repressor protein XylR
	Promoter	6431	6465	opt. PxyIA	Xylose-inducible promoter P <sub>xyIA</sub> including an optimized -35 region
	Region	6467	6495	XylR BS	XylR binding-site



## 7.2. Vector map of pC-HIS1623hp (#BMEG31)




	Type	Start	End	Name	Description
	Selectable genetic marker	541	6232	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	2097	814	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	3391	2390	repU	Gene of replication protein RepU
	Origin of replication	3543	3569	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	4969	3848	xylR	Xylose repressor gene
	Promoter	5115	5149	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	5151	5179	XylR BS	XylR binding-site
	Region	5220	5229	opt. RBS	Optimized ribosomal binding site
	MCS	5246	5299	MCS	Multiple Cloning Site
	Tag	5312	5329	His-Tag	6x histidine tag
	Origin of replication	5406	6069	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group






### 7.3. Vector map of pN-HIS-TEV1623hp (#BMEG32)



	Type	Start	End	Name	Description
	Selectable genetic marker	542	6274	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	2189	814	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	3391	2390	repU	Gene of replication protein RepU
	Origin of replication	3543	3569	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	5011	3845	xylR	Xylose repressor gene
	Promoter	5120	5138	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	5151	5179	XylR BS	XylR binding-site
	Region	5220	5229	opt. RBS	Optimized ribosomal binding site
	Tag	5258	5275	His-Tag	6x histidine tag
	MCS	5303	5354	MCS	Multiple Cloning Site
	Origin of replication	5448	6111	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group




#### 7.4. Vector map of pSP<sub>LipA</sub>-hp (#BMEG33)



	Type	Start	End	Name	Description
	Selectable genetic marker	542	6304	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	2097	814	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	3391	2390	repU	Gene of replication protein RepU
	Origin of replication	3543	3569	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	4969	3848	xylR	Xylose repressor gene
	Promoter	5115	5149	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	5151	5179	XylR BS	XylR binding-site
	Region	5220	5229	opt. RBS	Optimized ribosomal binding site
	Signal peptide	5237	5314	SPLipA	Signal sequence of <i>lipA</i> gene
	MCS	5315	5372	MCS	Multiple Cloning Site
	Tag	5384	5401	His-Tag	6x histidine tag
	Origin of replication	5493	6139	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group




## 7.5. Vector map of pSP<sub>Yoch</sub>H-hp (#BMEG34)



	Type	Start	End	Name	Description
	Selectable genetic marker	542	6304	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	2097	814	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	3391	2390	repU	Gene of replication protein RepU
	Origin of replication	3543	3569	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	4969	3848	xylR	Xylose repressor gene
	Promoter	5115	5149	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	5151	5179	xylR BS	xylR binding-site
	Region	5220	5229	opt. RBS	Optimized ribosomal binding site
	Signal peptide	5237	5308	SPYochH	Signal sequence of <i>yochH</i> gene
	MCS	5312	5365	MCS	Multiple Cloning Site
	Tag	5378	5395	His-Tag	6x histidine tag
	Origin of replication	5487	6133	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group




## 7.6. Vector map of p3STOP1623-2RBSHp (#BMEG35)



	Type	Start	End	Name	Description
	Origin of replication	692	1338	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group
	Selectable genetic marker	1589	732	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	3144	1861	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	4438	3437	repU	Gene of replication protein RepU
	Origin of replication	4590	4616	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	6016	4895	xylR	Xylose repressor gene
	Promoter	6162	6196	opt. PxyIA	Xylose-inducible promoter including an optimized -35 region
	Region	6198	6226	XylR BS	XylR binding-site
	Region	6267	6275	opt. RBS1	Optimized ribosomal binding site 1
	Region	6311	6321	opt. RBS2	Optimized ribosomal binding site 2
	MCS	6323	6366	MCS	Multiple Cloning Site




## 7.7. Vector map of pC-STREP1623hp (#BMEG36)



	Type	Start	End	Name	Description
	Selectable genetic marker	542	6256	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	2097	814	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	3391	2390	repU	Gene of replication protein RepU
	Origin of replication	3543	3569	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	4969	3848	xyIR	Xylose repressor gene
	Promoter	5115	5149	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	5151	5179	XylR BS	xylR binding-site
	Region	5220	5229	opt. RBS	Optimized ribosomal binding site
	MCS	5246	5302	MCS	Multiple Cloning Site
	Tag	5303	5326	Strep-Tag	Streptavidin tag II
	Origin of replication	5433	6096	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group




## 7.8. Vector map of pN-STREP-Xa1623hp (#BMEG37)



	Type	Start	End	Name	Description
	Tag	5258	5281	Strep-Tag	Streptavidin tag II
	MCS	5289	5338	MCS	Multiple Cloning Site
	Origin of replication	5433	6095	CoIE1 ori*	origin of replication ( <i>E. coli</i> ); CoIE1 incompatibility group
	Selectable genetic marker	542	6256	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	2097	811	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	3391	2387	repU	Gene of replication protein RepU
	Origin of replication	3543	3569	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	4969	3845	xylR	Xylose repressor gene
	Promoter	5115	5149	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	5151	5179	XylR BS	xylR binding-site
	Region	5220	5229	opt. RBS	Optimized ribosomal binding site




## 7.9. Vector map of pN-STREP-TEV1623hp (#BMEG38)



	Type	Start	End	Name	Description
	Selectable genetic marker	1408	551	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	2963	1677	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	4257	3253	repU	Gene of replication protein RepU
	Origin of replication	4409	4435	pBC16 ori	Origin of replication ( <i>Bacillus</i> )
	Gene	5835	4711	xylR	Xylose repressor gene
	Promoter	5981	6015	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	6017	6045	XylR BS	xylR binding-site
	Region	6086	6095	opt. RBS	Optimized ribosomal binding site
	Tag	6124	6147	Strep-Tag	Streptavidin tag II
	MCS	6182	6225	MCS	Multiple Cloning Site
	Origin of replication	6320	388	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group



## 7.10. Vector map of pMGBM19 (#BMEG39)




	Type	Start	End	Name	Description
	Origin of replication	922	1576	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group
	Selectable genetic marker	2598	1741	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	3733	3086	Cm	Chloramphenicol resistance ( <i>Bacillus</i> )
	Gene	5299	4136	xylR	Xylose repressor gene
	Promoter	5403	5436	PxyIA	Xylose-inducible promoter
	MCS	5525	5567	MCS	Multiple Cloning Site
	Region	5700	5702	Stop	Stop Codon
	Origin of replication	5808	535	pMB100 Replicon	Origin of replication ( <i>Bacillus</i> )






### 7.11. Vector map of helper plasmid pMMEc4 (#PEC04)



	Type	Start	End	Name	Description
	Selectable genetic marker	2648	1989	Cm	Chloramphenicol resistance ( <i>Bacillus</i> )
	Origin of replication	3010	3892	p15A ori	Origin of replication ( <i>E. coli</i> )
	Gene	4971	4093	araC	Gene of AraC regulatory protein
	Operator	5001	5016	Oara	araO2, operator sequence
	Promoter	5122	5150	ParaC	Arabinose inducible promoter
	Operator	5158	5169	Oara	araO1, operator sequence
	Operator	5210	5248	Oara	araI1/2, operator sequence
	Promoter	5272	5302	pBAD	Arabinose inducible promoter
	Gene	5329	6492	xyIR	Xylose repressor gene



## 7.12. Map of the control vector pGFP1624hp (#BMEG40C)



	Type	Start	End	Name	Description
	Gene	1421	255	xylR	Xylose repressor gene
	Promoter	1525	1559	opt. PxylA	Xylose-inducible promoter including an optimized -35 region
	Region	1561	1589	XylR BS	xylR binding-site
	Region	1630	1639	opt. RBS	Optimized ribosomal binding site
	Gene	1647	2390	gfp	Green fluorescent protein (Gfp)
	Origin of replication	2523	3169	ColE1 ori*	Origin of replication ( <i>E. coli</i> ); ColE1 incompatibility group
	Selectable genetic marker	4191	3334	Amp	Ampicillin resistance ( <i>E. coli</i> )
	Selectable genetic marker	5836	4463	Tet	Tetracycline resistance ( <i>Bacillus</i> )
	Gene	7040	6039	repU	Gene of replication protein RepU
	Origin of replication	7192	7218	pBC16 ori	Origin of replication ( <i>Bacillus</i> )



## 8. Order Information, Shipping, and Storage

Order#	Description	Amount
BMEG30	<i>Bacillus megaterium</i> vector p3STOP1623hp, lyophilized	10 µg
BMEG31	<i>Bacillus megaterium</i> vector pC-HIS1623hp, lyophilized	10 µg
BMEG32	<i>Bacillus megaterium</i> vector pN-HIS-TEV1623hp, lyophilized	10 µg
BMEG33	<i>Bacillus megaterium</i> vector pSP <sub>LipA</sub> -hp, lyophilized	10 µg
BMEG34	<i>Bacillus megaterium</i> vector pSP <sub>YocH</sub> -hp, lyophilized	10 µg
BMEG35	<i>Bacillus megaterium</i> vector p3STOP1623-2RBSHp, lyophilized	10 µg
BMEG36	<i>Bacillus megaterium</i> vector pC-STREP1623hp, lyophilized	10 µg
BMEG37	<i>Bacillus megaterium</i> vector pN-STREP-Xa1623hp, lyophilized	10 µg
BMEG38	<i>Bacillus megaterium</i> vector pN-STREP-TEV1623hp, lyophilized	10 µg
BMEG39	<i>Bacillus megaterium</i> vector, pMGBm19, lyophilized DNA	10 µg
BMEG40C	pGFP1624hp, high performance GFP expression vector, positive control, lyophilized DNA	10 µg
PEC04	<i>Escherichia coli</i> P <sub>xyIA</sub> repressing vector, pMMEc4, lyophilized DNA	10 µg
Shipped at RT, store lyophilized vectors at 4 °C and reconstituted vectors at -20 °C Vectors are <i>E. coli</i> / <i>B. megaterium</i> shuttle vectors.		

## Related Products

Order#	Description	Amount
BMEG02	<i>Bacillus megaterium</i> protoplast, strain WH320	5x500 µl
BMEG04	<i>Bacillus megaterium</i> protoplast, strain YYBm1	5x500 µl
BMEG50	<i>Bacillus megaterium</i> protoplast, strain MS941	5x500 µl
PR-ETA10010-01	MobiTEV Protease1, recombinant, His-Tag	1000 U
PR-ETA10010-05	MobiTEV Protease1, recombinant, His-Tag	10 x 1000 U
Shipped on dry ice, store at -80 °C		
PR-ETA10050-04	MoBiTEV Protease2 (TurboTEV), recombinant, GST- & His-Tag	1 mg (10,000 U)
PR-ETA10050-05	MoBiTEV Protease2 (TurboTEV), recombinant, GST- & His-Tag	10 mg (100,000 U)
GE-NUC10700-01	HS-Nuclease, recombinant Endonuclease (encoded by the same gene as Benzonase®)	50,000 U
Shipped on blue ice, store at -20 °C		



Order#	Description	Amount
<b>For His-Tag Purification</b>		
PR-HTK004	MoBiTec Ni-IDA Columns	4 columns
PR-HTK010	MoBiTec Ni-IDA Columns	10 columns
PR-HTK105	MobiSpin Ni-IDA Columns	5 columns
PR-HTK110	MobiSpin Ni-IDA Columns	10 columns
shipped at RT, store columns at RT		

## 9. Contact and Support

MoBiTec GmbH ◆ Lotzestrasse 22a ◆ D-37083 Goettingen ◆ Germany

**Customer Service** - General inquiries & orders

phone: +49 (0)551 707 22 0

fax: +49 (0)551 707 22 22

e-mail: [order@mobitec.com](mailto:order@mobitec.com)

**Technical Service** - Product information

phone: +49 (0)551 707 22 70

fax: +49 (0)551 707 22 77

e-mail: [info@mobitec.com](mailto:info@mobitec.com)

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)