

INSTRUCTIONS

ProFoldin

Ready-to-load DPPC Liposomes with Ammonium Sulfate


CATALOG NUMBER

DPC100AS

INTRODUCTION

Liposomal drug formulations provide great opportunities of improving the drug efficacy and toxicity profiles. Drug molecules with amine groups such as doxorubicin, daunorubicin, topotecan and irinotecan etc. can be loaded by ammonium-based pH gradient. The **Ready-to-load DPPC Liposomes with Ammonium Sulfate (Catalog number DPC100AS)** are high quality non-PEGylated liposomes that are ready to load drug molecules containing amine groups. The drug loading process is completed in about 3 hours after mixing the drug with the liposomes. The composition of the liposomes is 1,2-dipalmitoyl-sn-glycero-3-phophocholine (DPPC) and cholesterol in the weight ratio of 1.9:1. The total lipid concentration is 10 mg/ml. The liposomes are encapsulated with ammonium sulfate for drug remote loading. The average size of liposomes is about 90 nm with poly dispersity index (dpi) of 0.1.

Ready-to load DPPC liposomes with ammonium sulfate



The Ready-to-load DPPC Liposomes with Ammonium Sulfate (Catalog number DPC100AS) includes 10 ml liposomes with 10 mg/ml lipid concentration. The mass ratio of lipids is DPPC: cholesterol = 1.9:1. The concentration of ammonium sulfate encapsulated within the liposomes for drug loading is 200 mM. The buffer is 10 mM histidine, pH 6.5, 9.2 % sucrose. Liposomes are stored in a 2°C to 8°C refrigerator. DO NOT freeze liposomes.

ProFoldin

10 Technology Drive, Suite 40, Number 188 Hudson, MA 01749-2791 USA

Phone: (508) 735-2539 www.profoldin.com FAX: (508) 845-9258 info@profoldin.com

INSTRUCTIONS

DRUG LOADING PROTOCOL

Drug molecules containing at least one amine group are potentially loaded into the liposomes. The drug loading capacity depends on the property of the drug molecules. The typical drug to lipids ratio is from 1:5 to 1:20.

1. Drug loading

Mix the drug solution and the liposome and incubate the mixture in a 48°C water bather for 3 hours. Save the loaded liposomes at 2 to 8°C.

2. Drug encapsulation measurement

The yield of drug encapsulation is measured by a spin-column method using the **Liposome Drug Encapsulation Assay kit (Catalog number LDE10).** A volume of 50 μ l of sample was loaded on the spin column and eluted with 100 μ l of the elution buffer. The encapsulated dug is eluted and the non-encapsulated drug stays on the column. For detailed protocol, please see the Instruction of the **Liposome Drug Encapsulation Assay kit (Catalog number LDE10).**

RELATED PRODUCTS

Liposome products:

SLP20 Spin-columns for Liposome Purification
LDE10 Liposome Drug Encapsulation Assay Kit
LDD05 Liposome Drug Dissolution Assay Kit

LIP1000 MicroGram Lipid Assay Kit

SPS20 Liposome Plasma Stability Test Kit
DPPC002CP DPPC Liposomal Ciprofloxacin- 2 mg

PHPC002CP PEGylated HSPC Liposomal Ciprofloxacin- 2 mg

PHPC002DX PEGylated Liposomal Doxorubicin- 2 mg

DPC001AO Liposomal Acridine Orange Dye
DPC001RG Liposomal Rhodamine G Dye
DPC001FL Liposomal Fluorescein Dye

Nanodisc products:

SMA31-100MG Styrene - Maleic Acid Copolymer 3:1 Free Acid- 100 mg SMA31S-100MG Styrene - Maleic Acid Copolymer 3:1 Sodium Salt- 100 mg

Detergent and lipid products:

DAK1000 Detergent assay kit

LIP1000 MicroGram Lipid Assay Kit

For more information of liposome and nanodisc products, please visit our website at www.profoldin.com.