INSTRUCTIONS

ProFoldin Guanylate Kinase Assay Kit Plus


E. coli Guanylate Kinase Assay Kit Plus E. coli Guanylate Kinase Assay Kit Plus-500

Catalog number: GMK100KE Catalog number: GMK500KE

Introduction

Guanylate kinase is essential for recycling GMP and cGMP. It catalyzes the ATP-dependent phosphorylation of GMP into GDP. Guanylate kinase is an attractive antibacterial target. The guanylate

kinase assay is based on measurement of ADP generated from the kinase reaction. The assay is fluorescence-based and can be carried out using regular black 96-well or 384-well plates or cuvettes.



The *E. coli* GMP Kinase Assay Kit Plus (Catalog No. GMK100KE) includes 600 μl of 10 x assay buffer, 35 μl of 100 x GMP, 35 μl of 100 x E. coli GMP kinase, 35 μl of 100 x ATP, 35 μl of 100 x MUK Reagent A, 35 μl of 100 x MUK Reagent B and 350 μl of 10 x fluorescence dye for 100 assays in a 384-well assay format.

The *E. coli* GMP Kinase Assay Kit Plus-500 (Catalog No. GMK500KE) includes 2000 μl of 10 x assay buffer, 170 μl of 100 x GMP, 170 μl of 100 x E. *coli* GMP kinase, 170 μl of 100 x ATP, 170 μl of 100 x MUK Reagent A, 170 μl of 100 x MUK Reagent B and 1700 μl of 10 x fluorescence dye for 500 assays in a 384-well assay format.



INSTRUCTIONS

Assay Protocol

1. Reagent preparation:

- (1) 1 x reaction buffer: dilute the 10 x Reaction Buffer 10-fold with water. The 1 x reaction buffer is composed of 50 mM Tris-HCl, pH 8.0, 3 mM MgCl₂, 0.2 mM EDTA, 0.5 mM DTT, 50 mM NaCl, 0.003% Brij-35.
- (2) 10 x GMP: dilute the 100 x GMP solutions (5 mM) 10-fold with water.
- (3) 10 x ATP: dilute the 100 x ATP solutions (2 mM) 10-fold with water.
- (4) 10 x guanylate kinase: dilute the 100 x guanylate kinase 10-fold with 1 x reaction buffer.
- (5) 10 x MUK Reagent A: dilute the 100 x MUK Reagent A 10-fold with water.
- (6) 10 x MUK Reagent B: dilute the 100 x MUK Reagent B 10-fold with 1 x reaction buffer.
- (7) 1 x fluorescence dye: dilute the 10 x Fluorescence dye 10-fold with water.

2. Kinase assay

The kinase reaction volume is 30 µl.

- (1) In a standard black 384-well plate (Matrix 4318), mix 12 μl of water, 3 μl of 10 x reaction buffer, 3 μl of 10 x GMP, 3 μl of 10 x ATP, 3 μl of 10 x kinase. Incubate the reaction for 2 min.
- (2) Add 3 µl of 10 x MUK Reagent A, 3 µl of 10 x MUK Reagent B.
- (3) Incubate the reaction mixture at 37°C for 60 min.
- (4) Add 30 µl of the 1 x fluorescence dye into the 30 µl of the reaction mixture.
- (5) Measure the fluorescence intensity at 535 nm with excitation at 485 nm.

Assay Protocol for enzyme inhibition

In order to accurately measure the IC_{50} , it is important to make sure the assay is in the linear range. The linear range can be defined by enzyme concentration-dependence and time-dependence experiments. The enzyme concentration and time in the linear range is selected for the IC_{50} measurement.

For the enzyme inhibition experiment, typically 50 x stock solutions with a 2-fold serial dilution in water or DMSO are prepared. In the wells of the 384-well plate, $0.6~\mu l$ of the 50 x inhibitor is mixed with a premix composed of the buffer and enzyme for 5 min. Then the substrate is added. After incubation for 2 min, 3 μl of 10 x MUK Reagent A and 3 μl of 10 x MUK Reagent B are added. The total assay volume is 30 μl . At the end of the reaction, the dye is added and the fluorescence intensity is measured.

Related Products

S. pneumoniae Guanylate Kinase Assay Kit Plus-500 Catalog number. GMK500KN

E. coli Thymidylate Kinase Assay Kit Plus-500 Catalog number: TMK500KE

E. coli UMP Kinase Assay Kit Plus-500 Catalog number: UMK500KE

Human Thymidylate Kinase Assay Kit Plus-500 Catalog number: HTMK500KE

For more information of drug targets and enzyme assays, please visit www.profoldin.com.