



ProFoldin

10 Technology Drive, Suite 40, Number 188

Hudson, MA 01749-2791 USA

Phone: (508) 735-2539

FAX: (508) 845-9258

www.profoldin.com

info@profoldin.com

INSTRUCTIONS

ProFoldin


E. coli Methylerythritol Phosphate Cytidyltransferase (IspD) Assay Kits

E. coli Methylerythritol Phosphate Cytidyltransferase (IspD) Assay Kit Plus-100
Cat # ISPD100KE

E. coli Methylerythritol Phosphate Cytidyltransferase (IspD) Assay Kit Plus-500
Cat # ISPD500KE

INTRODUCTION

Methylerythritol phosphate cytidyltransferase (IspD) is one of the enzymes in the non-mevalonate pathway for isoprenoid biosynthesis that present in many pathogenic organisms and plants but absent in mammals. IspD is an attractive target for the development of novel antibiotics and herbicides. This enzyme catalyzes the formation of 4-diphosphocytidyl-2-C-methyl-D-erythritol from CTP and 2-C-methyl-D-erythritol 4-phosphate (MEP).



The *E. coli* Methylerythritol Phosphate Cytidyltransferase (IspD) Assay is based on measurement of the pyrophosphate generated from the IspD reaction. The pyrophosphate is converted to phosphate by pyrophosphatase and the phosphate is detected by light absorbance at 650 nm. The assay is in a 384-well or 96-well plate format. It can be used for high throughput screening in drug discovery.

The *E. coli* Methylerythritol Phosphate Cytidyltransferase (IspD) Assay Kit Plus-100 (Catalog # ISPD100KE) contains the reagents for 100 assays in a 384-well plate assay format including 400 μ l of 10 x Buffer, 35 μ l of 100 x MEP, 35 μ l of 100 x CTP, 35 μ l of 100 x *E. coli* IspD (200 nM), 35 μ l of 100 x pyrophosphatase (PPase, 10 U/ml), and 5 ml of Dye MPA3000 for phosphate detection.

The *E. coli* Methylerythritol Phosphate Cytidyltransferase (IspD) Assay Kit Plus-500 (Catalog # ISPD100KE) contains the reagents for 500 assays in a 384-well plate assay format including 2000 μ l of



ProFoldin

10 Technology Drive, Suite 40, Number 188

Hudson, MA 01749-2791 USA

Phone: (508) 735-2539 FAX: (508) 845-9258

www.profoldin.com info@profoldin.com

INSTRUCTIONS

10 x Buffer, 170 μ l of 100 x MEP, 170 μ l of 100 x CTP, 170 μ l of 100 x *E. coli* IspD (200 nM), 170 μ l of 100 x pyrophosphatase (PPase, 10 U/ml), and 25 ml of Dye MPA3000 for phosphate detection.

ASSAY PROTOCOL

The following assay protocol is based on the 384-well plate assay format. The reaction volume is 30 μ l and the final assay volume is 75 μ l. For 96-well plate assays, the reaction volume is 60 μ l and the final assay volume is 150 μ l. For detection using a cuvette, the reaction volume is 400 μ l and the final assay volume is 1000 μ l.

1. Reagent preparation:

For each 10 assay reactions,

- (1) Prepare 297 μ l of premix composed of 257.4 μ l of H₂O, 33 μ l of 10 x Buffer, 3.3 μ l of 100 x *E. coli* IspD and 3.3 μ l of 100 x PPase .
- (2) Prepare 33 μ l of 10 x Enzyme substrate by mixing 3.3 μ l of 100 x MEP, 3.3 μ l of 100 x CTP with 26.4 μ l of water.

2. Reaction:

Mix 27 μ l of the premix with 3 μ l of the 10 x Enzyme substrate in each well. Incubate the reaction mixture at 37°C for 60 min.

3. Detection:

Add 45 μ l of the Dye MPA3000 into the 30 μ l of the reaction mixture. Incubate for 5 min. Measure the light absorbance at 650 nm.

Assay Protocol for enzyme inhibition

The assay can be optimized in terms of assay window, assay linearity and sensitivity to competitive inhibitors. ProFoldin offers HTS assay development service. For more information, please visit our website at <http://www.profoldin.com/services.html>.

Related Products

ProFoldin offers high throughput assays for various drug targets and enzymes including DNA replication enzymes, RNA polymerases, kinases, bacterial cell wall synthesis enzymes and various metabolism enzymes, etc. For more information, please visit www.profoldin.com or send requests to info@profoldin.com.