



## ProFoldin

10 Technology Drive, Suite 40, Number 188  
Hudson, MA 01749-2791 USA  
FAX: (508) 845-9258 Tel: (508) 735-2539  
www.profoldin.com  
[info@profoldin.com](mailto:info@profoldin.com)

# INSTRUCTIONS

## ProFoldin

### M-MLV Reverse Transcriptase Assay Kit

**M-MLV Reverse Transcriptase Assay Kit**


**M-MLV Reverse Transcriptase Assay Kit Plus**

**Catalog No. MLV100K**

**Catalog No. MLV100KE**

### Introduction

The Moloney murine leukemia virus (M-MLV) causes lymphocytic leukemia in mice. M-MLV reverse transcriptase (M-MLV RT) synthesizes a complementary DNA strand from single-stranded RNA or DNA. The M-MLV Reverse Transcriptase Assay is based on measurement of the DNA molecules synthesized by M-MLV RT. The assay can be performed in a 384-well or 96-well plate format for tests of M-MLV reverse transcriptase activities and high throughput screening of inhibitors.



The **M-MLV Reverse Transcriptase Assay Kit** (Catalog number: MLV100K) includes 400  $\mu$ l of 10 x Buffer, 33  $\mu$ l of 100 x Template, 33  $\mu$ l of 100 x dNTPs, 1550  $\mu$ l of 2 x Dye and 1550  $\mu$ l of 50 mM EDTA. It is for 100 assays of M-MLV Reverse Transcriptase. The assay kit includes all reagents except the enzyme.

The **M-MLV Reverse Transcriptase Assay Plus** (Catalog number: MLV100KE) includes all reagents in the **M-MLV Reverse Transcriptase Assay Kit** (Catalog number: MLV100K) plus the enzyme, 7  $\mu$ l 500 x M-MLV RT.

### ASSAY PROTOCOL

The following assay protocol is based on the 384-well plate assay format (plate type: Matrix 4318 or alike). The reaction volume is 30  $\mu$ l and the final assay volume is 60  $\mu$ l. For 96-well plate assays (plate type: Costar 3915 or alike), the reaction volume is 60  $\mu$ l and the final assay volume is 120  $\mu$ l.

#### 1. Reagent preparation:

(1) 10 x Template: Dilute the 100 x Template 10-fold with water. Each assay uses 3  $\mu$ l of 10 x Template.



## ProFoldin

10 Technology Drive, Suite 40, Number 188

Hudson, MA 01749-2791 USA

FAX: (508) 845-9258 Tel: (508) 735-2539

[www.profoldin.com](http://www.profoldin.com)

[info@profoldin.com](mailto:info@profoldin.com)

# INSTRUCTIONS

---

- (2) 10 x enzyme: Dilute the 500 x M-MLV RT (200 Unit /  $\mu$ l) 50-fold with the 1 x assay buffer. Each assay uses 3  $\mu$ l of 10 x enzyme.
- (3) 10 x dNTP mix: Dilute the 100 x dNTP mix (10 mM dATP, 10 mM dGTP) 10-fold with water. Each assay uses 3  $\mu$ l of 10 x dNTP mix.
- (4) 1 x dye: Dilute the 2 x fluorescence dye 2-fold with 50 mM EDTA. Each assay uses 30  $\mu$ l of 1 x dye.

## 2. Reaction:

The total volume of each reaction mixture is 30  $\mu$ l including 18  $\mu$ l of H<sub>2</sub>O, 3  $\mu$ l of 10 x Buffer, 3  $\mu$ l of 10 x template, 3  $\mu$ l of 10 x enzyme, 3  $\mu$ l of 10 x dNTP mix. Incubate the reaction mixture at 37°C for 60 min.

## 3. Detection:

Mix 30  $\mu$ l of the 1 x fluorescence dye with 30  $\mu$ l of the reaction mixture. Measure the fluorescence intensity at 535 nm using the excitation wavelength at 485 nm.

## Assay Protocol for enzyme inhibition

The assay can be optimized in terms of assay window, assay linearity and sensitivity to competitive inhibitors. ProFoldin offers HTS assay development service. For more information, please visit our website at <http://www.profoldin.com/services.html>.

## Related Products

AMV100K	AMV Reverse Transcriptase Assay Kit
AMV100KE	AMV Reverse Transcriptase Assay Kit Plus
HIV100K	HIV Reverse Transcriptase Assay Kit
HIV100KE	HIV Reverse Transcriptase Assay Kit Plus
RPA100KE	<i>E. coli</i> RNA Polymerase Assay Kit Plus
RPA100KSE	<i>S. aureus</i> RNA Polymerase Assay Kit Plus
DPA100KE	<i>E. coli</i> DNA Polymerase III Alpha Assay Kit Plus
DPA100KH	<i>H influenzae</i> DNA Polymerase Assay Kit Plus
DPA100KN	<i>S. pneumoniae</i> DNA Polymerase Assay Kit Plus
HDP A100K	Human DNA Polymerase Alpha Assay Kit
DPB100K	Human DNA Polymerase Beta Assay Kit
DPG100K	Human DNA Polymerase Gamma Assay Kit
DPG100KE	Human DNA Polymerase Gamma Assay Kit Plus
MRPA100K	Human Mitochondrial RNA Polymerase Assay Kit

## More information of drug targets and enzyme assays

For more information of drug targets and enzyme assays, please visit [www.profoldin.com](http://www.profoldin.com) or send emails to [info@profoldin.com](mailto:info@profoldin.com).