ProFoldin 10 Technology Drive, Suite 40, Number 188 Hudson, MA 01749-2791 USA Tel: (508) 735-2539 FAX: (508) 845-9258 www.profoldin.com info@profoldin.com

INSTRUCTIONS


ProFoldin MicroMolar Universal Kinase Assay Kits


MicroMolar Universal Kinase Assay Kit Catalog Number: MUK100K MicroMolar Universal Kinase Assay Kit-1000 Catalog Number: MUK1000K


Introduction

The MicroMolar Universal Kinase Assay Kit is based on measurement of ADP generated from the kinase reaction. It is used to measure the activities of purified kinases producing ADP at concentrations ranging from $0.1 \mu M$ to $10 \mu M$.

The assay is fluorescence-based and can be carried out using regular black or white 96-well or 384-well plates or micro-cuvettes. It is important to avoid using unnecessarily high ATP concentrations in the kinase assays for the following reasons. One is to allow the assay to detect ATP-competitive inhibitors. The other is to avoid a high ADP background from ADP contamination in the ATP samples. It is recommended to use ATP concentrations at 20 µM (micromolar).







ProFoldin 10 Technology Drive, Suite 40, Number 188 Hudson, MA 01749-2791 USA Tel: (508) 735-2539 FAX: (508) 845-9258 www.profoldin.com info@profoldin.com

INSTRUCTIONS

The **MicroMolar Universal Kinase Assay Kit** (Catalog No. **MUK100K**) includes 600 μl of 10 x assay buffer, 35 μl of 100 x ATP, 35 μl of 100 x MUK Reagent A, 35 μl of 100 x MUK Reagent B, 30 μl of 1 mM ADP and 350 μl of 10 x fluorescence dye for 100 assays in a 384-well assay format.

The **MicroMolar Universal Kinase Assay Kit-1000** (Catalog No. **MUK1000K**) includes 4000 µl of 10 x assay buffer, 350 µl of 100 x ATP, 350 µl of 100 x MUK Reagent A, 350 µl of 100 x MUK Reagent B, 200 µl of 1 mM ADP and 3500 µl of 10 x fluorescence dye for 1000 assays in a 384-well assay format.

Assay Protocol

1. Reagent preparation:

- (1) 1 x MUK Buffer: dilute the 10 x MUK Buffer 10-fold with water. The 1x MUK Buffer is composed of 50 mM Tris-HCl, pH 8.0, 50 mM NaCl, 3 mM MgCl₂, 0.5 mM DTT, 0.2 mM EDTA, 0.003% Brij-35.
- (2) 10 x MUK Reagent A: dilute the 100 x MUK Reagent A 10-fold with water
- (3) 10 x MUK Reagent B: dilute the 100 x MUK Reagent B 10-fold with 1x MUK Buffer
- (4) 1 x fluorescence dye: dilute the 10 x fluorescence dye 10-fold with water.

2. Kinase assay

The kinase reaction volume is 30 µl.

- (1) In a standard black 384-well plate, mix 12 µl of water, 3 µl of 10 x reaction buffer, 3 µl of 10 x kinase substrate, 3 µl of 10 x ATP, 3 µl of freshly prepared 10 x kinase. Incubate the reaction for 2 min.
- (2) Add 3 µl of 10 x MUK Reagent A, 3 µl of 10 x MUK Reagent B.
- (3) Incubate the reaction mixture at 37°C for 60 min.
- (4) Add 30 µl of the 1 x fluorescence dye into the 30 µl of the reaction mixture.
- (5) Measure the fluorescence intensity at 535 nm with excitation at 485 nm.

Enzyme inhibition

For the enzyme inhibition experiment, typically $50 \, x$ stock solutions with a 2-fold serial dilution in water or DMSO are prepared. In the wells of the 384-well plate, $0.6 \, \mu l$ of the $50 \, x$ inhibitor is mixed with a premix composed of the buffer and enzyme for $5 \, \text{min}$. Then the substrate is added. After incubation for $2 \, \text{min}$, $3 \, \mu l$ of $10 \, x$ MUK Reagent A and $3 \, \mu l$ of $10 \, x$ MUK Reagent B are added. The total assay volume is $30 \, \mu l$. At the end of the reaction, the dye is added and the fluorescence intensity is measured.

Related products

Thymidylate Kinase Assay Kit Plus-500

Guanylate Kinase Assay Kit Plus-500

Catalog number: TMK500KE

Catalog number: GMK500KE

Catalog number: NAK500KE

Catalog number: NAK500KE

UMP Kinase Assay Kit Plus-500

Catalog number: UMK500KE

Human Thymidylate Kinase Assay Kit Plus-500

Catalog Number: HTMK500KE