

ProFoldin 10 Technology Drive, Suite 40, Number 188 Hudson, MA 01749-2791 USA Phone: (508) 735-2539 FAX: (508) 845-9258 www.profoldin.com info@profoldin.com

INSTRUCTIONS

ProFoldin MurE Assay Kits

E. coli MurE Assay Kit Plus-100 *E. coli* MurE Assay Kit Plus-500

Catalog No. MURE100KE Catalog No. MURE500KE

INTRODUCTION

MurE or UDP-MurNAc-tripeptide ligase is the third amino acid-adding enzymes in the pathway of peptidoglycan biosynthesis in bacteria . It is an essential enzyme and attractive target for antibacterial drug discovery. MurE catalyses the addition of lysine or meso-diaminopimelic acid (DAP) into the MurD product UDP-MurNAc-dipeptide in bacteria generating the UDP-MurNAc-tripeptide. The ligation reaction is coupled to the hydrolysis of ATP forming ADP and inorganic phosphate.



UDP-MurNAc-dipeptide meso-diaminopimelic acid



The *E. coli* MurE Assay is based on measurement of the inorganic phosphate generated from the MurE reaction. The inorganic phosphate is detected by light absorbance at 650 nm. The assay reactions and detection can be performed by using 384-well or 96-well assay plates. Alternatively, the assay reaction can be carried out in Eppendorf tubes and the signal is measured using a cuvette. The high throughput assay can be used for screening inhibitors of *E.coli* MurE in drug discovery research. It may also be used for characterization of *E.coli* MurE.

The *E. coli* MurE Assay Kit Plus-100 (Catalog No. MURE100KE) contains the reagents for 100 assays in a 384-well plate assay format including 400 μ l of 10 x Buffer, 35 μ l of 100 x UDP-MurNAc-L-Ala-D-Glu (UDP-MurNAc-dipeptide, UMAG), 35 μ l of 100 x DAP, 35 μ l of 100 x ATP , 35 μ l of 100 x *E. coli* MurE (5000 nM) and 5 ml of Dye MPA3000 for phosphate detection.

The *E. coli* MurE Assay Kit Plus-500 (Catalog No. MURE500KE) contains the reagents for 500 assays in a 384-well plate assay format including 2000 µl of 10 x Buffer, 170 µl of 100 x UDP-MurNAc-L-Ala-D-Glu (UDP-MurNAc-dipeptide, UMAG), 170 µl of 100 x DAP, 170 µl of 100 x ATP, 170 µl of 100 x *E. coli* MurE (5000 nM) and 25 ml of Dye MPA3000 for phosphate detection.

ASSAY PROTOCOL

The following assay protocol is based on the 384-well plate assay format. The reaction volume is 30 µl

INSTRUCTIONS



and the final assay volume is 75 μ l. For 96-well plate assays, the reaction volume is 60 μ l and the final assay volume is 150 μ l. For detection using a cuvette, the reaction volume is 400 μ l and the final assay volume is 1000 μ l.

1. Reagent preparation:

For each 10 assay reactions,

- Prepare 297 μl of premix composed of 261 μl of H₂O, 33 μl of 10 x Buffer and 3.3 μl of 100 x *E. coli* MurE.
- (2) Prepare 33 µl of 10 x Enzyme substrate by mixing 3.3 µl of 100 x UMAG, 3.3 µl of 100 x DAP, 3.3 µl of 100 x ATP, and 23.1 µl of water.

2. Reaction:

Mix 27 μ l of the premix with 3 μ l of the 10 x Enzyme substrate in each well. Incubate the reaction mixture at 37°C for 60 min.

3. Detection:

Add 45 μ l of the Dye MPA3000 into the 30 μ l of the reaction mixture. Incubate for 5 min. Measure the light absorbance at 650 nm.

ASSAY LINEARITY TEST

Follow the same protocol described above except mixing 27 μ l of the premix with 3 μ l of the 10 x Enzyme substrate at different time points. Plot the reaction signal versus the reaction time to define the linear range.

IC50 MEASUREMENT OF ENZYME INHIBITORS

The concentration range of the inhibitor to be tested depends on the potency of the inhibitor. In general, the maximum concentration is about 10 to 20 fold higher than the IC50 value. The following protocol is for IC50 measurement of one inhibitor with IC50 values around 10 μ M.

- 1. In 8 assay wells, add 0.6 μ l of 2-fold serial dilution solutions of the inhibitor from 5 mM to 0.039 mM in water or DMSO. In one control well, add 0.6 μ l of water or DMSO. In another control well, add 0.6 μ l of 1 M EDTA.
- 2. Prepare 297 µl of premix and 33 µl of 10 x Enzyme substrate as described above.
- 3. Mix 26.4 μ l of the premix and 0.6 μ l of the 50 x inhibitor for 5 min.
- 4. Add 3 µl of the 10 x Enzyme substrate and incubate the assay reaction for the time in the linear range.
- 5. Add 45 μ l of the Dye MPA3000 into the 30 μ l of the reaction mixture. Incubate for 5 min. Measure the light absorbance at 650 nm.
- 6. Calculate IC50s using a computer IC50 fitting software.

Note: If DMSO is used to make the inhibitor solutions, the final concentration of DMSO in the assay is 2 %. It is important to make sure that 2 % DMSO does not affect the enzyme activity. Otherwise, the assay condition should be adjusted accordingly to keep the sufficient signal to background ratio and the assay linearity.