



## ProFoldin

10 Technology Drive, Suite 40, Number 188  
Hudson, MA 01749-2791 USA  
Tel: (508) 735-2539 FAX: (508) 845-9258  
[www.profoldin.com](http://www.profoldin.com)  
[info@profoldin.com](mailto:info@profoldin.com)

# INSTRUCTIONS

## ProFoldin

### Beta Octyl Glucoside ( $\beta$ -OG) Assay Kit


---

CATALOG NUMBER            **OG100K**

#### INTRODUCTION

The Beta Octyl Glucoside ( $\beta$ -OG) Assay Kit is for measurement of detergent  $\beta$ -OG at concentrations below its CMC value. The assay is based on the principle that  $\beta$ -OG interacts with the dye and enhances the light absorbance at 580 nm. It is compatible with most buffers and salts. It is not compatible with phosphate. Molecules such as ATP or ADP that release phosphate may interfere with the assay.

The assay does not detect beta nonyl glucoside ( $\beta$ -NG) or SDS. Please use Beta Nonyl Glucoside ( $\beta$ -NG) Assay Kit (Catalog # NG100K) for measurement of  $\beta$ -NG concentrations; NanoGram SDS Assay Kit (Catalog # SDS200) for measurement of SDS concentrations; and Detergent Assay Kit (Catalog # DAK1000) for other detergents.



The assay kit includes 10 ml of Reagent 1, 5 ml of dilution buffer and 0.1 ml of 100 x Dye for OG. It is for measurement of 100 samples using 96-well plates. A cuvette may also be used for the measurement of detergent concentrations.

#### ASSAY PROTOCOL

The following protocol is for assays using a standard black 96-well plate (Costar 3915 or Greiner 655076). Adjust the assay reagent volumes proportionally for measurements using cuvettes.

- Sample preparation:** Prepare 50  $\mu$ l of standard  $\beta$ -OG solutions with a series of concentrations from zero to 0.2% in 10 mM Tris-HCl, pH 7.5, 100 mM NaCl. Prepare 1 x dye by dilution of the 100 x dye with Reagent 1 (100-fold dilution). Each assay needs 100  $\mu$ l of 1 x dye.
-



## ProFoldin

10 Technology Drive, Suite 40, Number 188  
Hudson, MA 01749-2791 USA  
Tel: (508) 735-2539 FAX: (508) 845-9258  
[www.profoldin.com](http://www.profoldin.com)  
[info@profoldin.com](mailto:info@profoldin.com)

## INSTRUCTIONS

---

- Detection:** Mix 50 µl of the detergent solution with 100 µl of 1 x dye and incubate the solution at room temperature for 5 min. Read the light absorbance at 580 nm ( $A_{580}$ ).
- Data Analysis:** Plot the  $A_{580}$  values and the  $\beta$ -OG concentration [ $\beta$ -OG] to generate the linear standard curve.

$$A_{580} = a [\beta\text{-OG}] + b$$

Where the  $A_{580}$  values are from experimental data, the **a** and **b** values are from the linear fitting between the  $A_{580}$  values and the  $\beta$ -OG concentrations.

### UNKNOWN SAMPLES

Follow the same procedure to measure the light absorbance at 580 nm ( $A_{580}$ ) values from the unknown samples. Calculate the  $\beta$ -OG concentrations in the unknown samples using the  $A_{580}$  values from the unknown samples and the **a** and **b** values from the standard curve.

$$[\beta\text{-OG}] = (A_{580} - b) / a$$

### RELATED PRODUCTS

|         |  |
|---------|--|
| NG100K  | Beta Nonyl Glucoside Assay Kit                           |
| SDS200K | NanoGram SDS Assay Kit                                   |
| DAK1000 | Detergent Assay Kit                                      |
| CMC1000 | Detergent Critical Micelle Concentration (CMC) Assay Kit |
| LIP1000 | MicroGram Lipid Assay Kit                                |
| NPA1000 | NanoMolar Phosphate Assay Kit                            |
| MPA3000 | MicoMolar Phosphate Assay Kit                            |
| PPD1000 | MicroMolar Polyphosphate Assay Kit                       |
| EPA001  | Easy Protein Assay Reagent                               |
| HIS200  | MicroMolar Histidine Assay Kit                           |
| CYS200  | MicroMolar Cysteine Assay kit                            |
| PEP200  | Peptide Assay Kit  |
| PAA100K | MicroMolar Primary Amine Assay Kit                       |
| CAK1000 | Coenzyme A Assay Kit                                     |
| EDTA200 | MicroMolar EDTA Assay kit                                |
| DTT200  | MicroMolar DTT Assay kit                                 |
| MAD100K | MicroMolar ADP Assay kit                                 |
| MUD100K | MicroMolar UDP assay kit                                 |
| MCA1000 | MicroMolar Copper Assay Kit                              |
| NZA1000 | NanoMolar Zinc Assay Kit                                 |
| NMA1000 | NanoMolar Nickel / Cobalt Assay Kit                      |
| CLA100  | MicroMolar Chloride Assay Kit                            |
| MSA200  | MicroMolar Sulfate Assay Kit                             |
| CPT200  | MicroMolar Cisplatin Assay Kit                           |

For more information of concentration assays, please visit [www.profoldin.com](http://www.profoldin.com).