# **INSTRUCTIONS**



 ProFoldin

 10 Technology Drive, Suite 40, Number 188

 Hudson, MA 01749-2791
 USA

 Tel: (508) 735-2539
 FAX: (508) 845-9258

 www.profoldin.com
 info@profoldin.com

# ProFoldin MicroMolar Polyphosphate Assay Kit

# CATALOG NUMBER PPD1000

# **INTRODUCTION**

Inorganic polyphosphate is a linear molecule composed of tens or hundreds of phosphate residues linked together. In bacteria, polyphosphate kinase (PPK) converts polyphosphate and ADP to ATP. Lacking PPK activity resulted in polyphosphate deficiency and failure in expression of RpoS (a sigma factor for RNA polymerase) that leads to cell death. Thus PPK is a potential antibacterial drug target. The MicroMolar Polyphosphate Assay Kit is for measurement of micromolar concentrations of polyphosphate. The assay is based on increase of the fluorescence intensity (emission 550 nm, excitation 415 nm) of the kit fluorescence dye PPD upon binding to polyphosphate. The assay is compatible with regular buffers and various phosphate compounds including inorganic phosphate, pyrophosphate, ATP, ADP and AMP. The assay kit can be used for measurements of polyphosphate in biological samples or environmental water samples.



The kit (catalog number PPD1000) includes  $300 \ \mu l$  of  $100 \ x$  PPD dye and  $30 \ \mu l$  of  $1 \ mM$  sodium polyphosphate (a 45-mer). It is for 1000 assays using 384-well plates (30  $\ \mu l$  of sample volume) or 500 samples using 96-well plates (60  $\ \mu l$  of sample volume). Cuvettes may also be used for measurements.

#### **Reference:**

- 1. Achbergerová L. et al, Degradation of polyphosphates by polyphosphate kinases from Ruegeria pomeroyi. Biotechnology Letters. Volume 36, Issue 10, pp 2029-2035 (2014).
- 2. Zhang J. et al, A Fast Sensor For in Vivo Quantification of Cytosolic Phosphate in Saccharomyces Cerevisiae, Biotechnol Bioeng.Vol 112, pp 1033-1046 (2015).

# **INSTRUCTIONS**



# PROTOCOL

### Standard curve

1. **Sample preparation:** Prepare the polyphosphate solutions, 45-mer sodium polyphosphate, at a series of concentrations ranging from 50  $\mu$ M to zero in a 10 mM HEPES, pH 7.4 buffer. Dilute the 100 x PPD dye 100-fold with water to make the 1 x PPD dye.

Note: The fluorescence intensity has a linear relationship with the polyphosphate concentrations in the range of 1 to 50  $\mu$ M. A polyphosphate concentration higher than 50  $\mu$ M suppresses the fluorescence signal.

2. **Detection:** Mix 30  $\mu$ l of the sample with 30  $\mu$ l of the 1x PPD dye solution for 5 min. Read the fluorescence intensity at 550 nm (excitation 415 nm).

3. **Data Analysis**: Plot the fluorescence intensity **Fc** and the polyphosphate concentration **[Polyphosphate]** to generate the linear standard curve.

# Fc = a [Polyphosphate] + b

Where the **Fc** values are from experimental data, the **a** and **b** values are from the linear fitting between the **Fc** values and the polyphosphate concentrations.

#### **UNKNOWN SAMPLES**

Follow the same procedure to measure the fluorescence intensity Fc values from the unknown samples. Calculate the polyphosphate concentrations in the unknown samples using the Fc values from the unknown samples and the **a** and **b** values from the standard curve.

# [Polyphosphate] = (Fc - b) / a

# **RELATED PRODUCTS**

MPA3000	MicroMolar Phosphate Assay Kit
NPA1000	NanoMolar Phosphate Assay Kit
EPA001	Easy Protein Assay Reagent
HIS200	MicroMolar Histidine Assay Kit
CYS200	MicroMolar Cysteine Assay kit
PAA100K	MicroMolar Primary Amine Assay Kit
EDTA200	MicroMolar EDTA Assay kit
DTT200	MicroMolar DTT Assay kit
DAK1000	Detergent assay kit
LIP1000	MicroGram Lipid Assay Kit

For more information of concentration assays and enzyme essays, please visit www.profoldin.com.