

INSTRUCTIONS



ProFoldin
10 Technology Drive, Suite 40, Number 188
Hudson, MA 01749-2791 USA
FAX: (508) 845-9258 Tel: (508) 735-2539
www.profoldin.com
info@profoldin.com

ProFoldin SARS-CoV-2 RNA Polymerase (RdRp) Assay Kits

SARS-CoV-2 RNA polymerase assay kit-100 (100 assays)

Catalog No. S2RPA100K


SARS-CoV-2 RNA polymerase assay kit plus -20 (20 assays)

Catalog No. S2RPA020KE

SARS-CoV-2 RNA polymerase assay kit plus -100 (100 assays) **Catalog No. S2RPA100KE**

Introduction

The COVID-19 pandemic is caused by infection of SARS-CoV-2 or severe acute respiratory syndrome coronavirus 2, a positive-sense single-stranded RNA virus. The RNA dependent RNA polymerase (RdRp) is a key drug target for this deadly virus. The SARS-CoV-2 RdRP assay is based on detection of the RNA molecules synthesized by the RNA polymerase in the presence of an RNA template and NTPs. The assay can be performed in a 384-well or 96-well plate format for tests of the enzyme activities of the virus RNA polymerase (SARS-CoV-2 RdRp) and high throughput screening of inhibitors.



The **SARS-CoV-2 RNA polymerase assay kit-100** (Catalog No.S2RPA100K) includes 300 μ l of 10 x Buffer, 50 μ l of 50 x RNA template, 60 μ l of 50 x NTPs (ATP, UTP, GTP, CTP) and 700 μ l of 10 x fluorescence dye. It is for 100 assays of SARS-CoV-2 RNA polymerase in a 384-well plate format or 50 assays in a 96-well plate format. The assay kit includes all reagents except the enzyme.

The **SARS-CoV-2 RNA polymerase assay kit plus-20** (Catalog No.S2RPA020KE) includes 75 μ l of 10 x Buffer, 10 μ l of 50 x RNA template, 12 μ l of 50 x NTPs (ATP, UTP, GTP, CTP), 150 μ l of 10 x fluorescence dye and 10 μ l of 50 x SARS-CoV-2 RdRp. It is for 20 assays of SARS-CoV-2 RNA polymerase in a 384-well plate format or 10 assays in a 96-well plate format.

The **SARS-CoV-2 RNA polymerase assay kit plus-100** (Catalog No.S2RPA100KE) includes 300 μ l of 10 x Buffer, 50 μ l of 50 x RNA template, 60 μ l of 50 x mM NTPs (ATP, UTP, GTP, CTP), 700 μ l of 10 x fluorescence dye and 50 μ l of 50 x SARS-CoV-2 RdRp. It is for 100 assays of SARS-CoV-2 RNA polymerase in a 384-well plate format or 50 assays in a 96-well plate format.

INSTRUCTIONS



ProFoldin

10 Technology Drive, Suite 40, Number 188
Hudson, MA 01749-2791 USA
FAX: (508) 845-9258 Tel: (508) 735-2539
www.profoldin.com
info@profoldin.com

ASSAY PROTOCOL

The following assay protocol is based on the 384-well plate assay format (plate type: Matrix 4318 or alike). The reaction volume is 25 µl and the final assay volume is 90 µl. For 96-well plate assays (plate type: Costar 3915 or alike), the reaction volume is 50 µl and the final assay volume is 180 µl.

1. Reaction:

The total volume of each reaction mixture is 25 µl including 21 µl of H₂O, 2.5 µl of 10 x Buffer, 0.5 µl of 50 x template, 0.5 µl of 50 x RNA polymerase and 0.5 µl of 50 x NTPs. Incubate the reaction mixture at 34°C for 2 hours.

2. Detection:

Dilute the 10x fluorescence dye with water to make the 1 x dye. Mix 65 µl of the 1 x dye with 25 µl of the reaction mixture. Measure the fluorescence intensity at 535 nm using the excitation wavelength at 485 nm in 10 min.

Assay Protocol for enzyme inhibition

The assay can be optimized in terms of assay linearity and sensitivity to competitive inhibitors.

ProFoldin offers HTS assay development service. For more information, please visit our website at <http://www.profoldin.com/services.html>.

Related Products

RNA polymerase assay kits:

VRT100K	Viral (Flavivirus) RNA-dependent RNA Polymerase Assay Kit - 100
RPA100KE	<i>E. coli</i> RNA Polymerase Assay Kit Plus
RPA100KSE	<i>S. aureus</i> RNA Polymerase Assay Kit Plus
MRPA100KE	Human Mitochondrial RNA Polymerase Assay Kit Plus
T7RPA100KE	T7 RNA Polymerase Assay Kit Plus
EGA100KE	<i>E. coli</i> DNA Primase Assay Kit Plus
HGA100KE	<i>H. influenzae</i> DNA Primase Assay Kit Plus
AGA100KE	<i>S. aureus</i> DNA Primase Assay Kit Plus
PGA100KE	<i>S. pneumoniae</i> DNA Primase Assay Kit Plus
PNP100KE	<i>E. coli</i> polynucleotide phosphorylase (PNPase) assay plus

DNA polymerase assay kits:

DPA100KE	<i>E. coli</i> DNA Polymerase III Alpha Assay Kit Plus
HDPA100KE	Human DNA Polymerase Alpha Assay Kit Plus
DPB100KE	Human DNA Polymerase Beta Assay Kit Plus
DPG100KE	Human DNA Polymerase Gamma Assay Kit Plus
HIV100KE	HIV Reverse Transcriptase Assay Kit Plus
AMV100KE	AMV Reverse Transcriptase Assay Kit Plus
MLV100KE	M-MLV Reverse Transcriptase Assay Kit Plus

For more information of drug targets and enzyme assays, please visit www.profoldin.com.