INSTRUCTIONS



ProFoldin *E. coli* Thymidylate Kinase Assay Kit Plus

E. coli Thymidylate Kinase Assay Kit Plus *E. coli* Thymidylate Kinase Assay Kit Plus-500


Catalog number: TMK100KE Catalog number: TMK500KE

Introduction

Thymidylate kinase is an important enzyme in the dTTP synthesis pathway for DNA synthesis. The function of this enzyme is to catalyze the phosphorylation of dTMP using ATP to form dTDP and ADP.

 $dTMP + ATP \longrightarrow dTDP + ADP$

Thymidylate kinase is an attractive antibacterial target. The thymidylate kinase assay is based on measurement of ADP generated from the kinase reaction. The assay is fluorescence-based and can be carried out using regular black 96-well or 384-well plates or micro-cuvettes.



The *E. coli* Thymidylate Kinase Assay Kit Plus (Catalog No. TMK100KE) includes 600 µl of 10 x Buffer, 35 µl of 100 x dTMP, 35 µl of 100 x *E. coli* Thymidylate kinase, 35 µl of 100 x ATP, 35 µl of 100 x MUK Reagent A, 35 µl of 100 x MUK Reagent B and 300 µl of 10 x fluorescence dye for 100 assays in a 384-well assay format.

The *E. coli* Thymidylate Kinase Assay Kit Plus-500 (Catalog No. TMK500KE) includes 2000 µl of 10 x Buffer, 170 µl of 100 x dTMP, 170 µl of 100 x *E. coli* Thymidylate kinase, 170 µl of 100 x ATP, 170 µl of 100 x MUK Reagent A, 170 µl of 100 x MUK Reagent B and 1700 µl of 10 x fluorescence dye for 500 assays in a 384well assay format.

INSTRUCTIONS



Assay Protocol

1. Reagent preparation:

- (1) 1 x reaction buffer: dilute the 10 x Buffer 10-fold with water. The 1 x buffer is composed of 50 mM Tris-HCl, pH 8.0, 3 mM MgCl₂, 0.2 mM EDTA, 0.5 mM DTT, 50 mM NaCl, 0.003% Brij-35.
- (2) 10 x dTMP: dilute the 100 x dTMP solutions (5 mM) 10-fold with water.
- (3) 10 x ATP: dilute the 100 x ATP solutions (5 mM) 10-fold with water.
- (4) 10 x thymidylate kinase: dilute the 100 x thymidylate kinase 10-fold with 1 x buffer.
- (5) 10 x MUK Reagent A: dilute the 100 x MUK Reagent A 10-fold with water.
- (6) 10 x MUK Reagent B: dilute the 100 x MUK Reagent B 10-fold with 1 x buffer.
- (7) 1 x fluorescence dye: dilute the 10 x Fluorescence dye 10-fold with water.

2. Kinase assay

The kinase reaction volume is 30 µl.

- In a standard black 384-well plate (Matrix 4318), mix 12 μl of water, 3 μl of 10 x Buffer, 3 μl of 10 x dTMP, 3 μl of 10 x ATP, 3 μl of 10 x kinase. Incubate the reaction for 2 min.
- (2) Add 3 μ l of 10 x MUK Reagent A, 3 μ l of 10 x MUK Reagent B.
- (3) Incubate the reaction mixture at 37° C for 60 min.
- (4) Add 30 μ l of the 1 x fluorescence dye into the 30 μ l of the reaction mixture.
- (5) Measure the fluorescence intensity at 535 nm with excitation at 485 nm.

Assay Protocol for enzyme inhibition

The assay can be optimized in terms of assay window, assay linearity and sensitivity to competitive inhibitors. ProFoldin offers HTS assay development service. For more information, please visit our website at <u>http://www.profoldin.com/services.html</u>.

Reference:

Choi J.Y. et al, Structure Guided Development of Novel Thymidine Mimetics targeting Pseudomonas aeruginosa Thymidylate Kinase: from Hit to Lead Generation, J Med Chem. January 26; 55(2): 852–870 (2012).

Related products

E. coli Guanylate Kinase Assay Kit Plus-500
E. coli NAD+ Kinase Assay Kit Plus-500
E. coli UMP Kinase Assay Kit Plus-500
Human Thymidylate Kinase Assay Kit Plus-500

Catalog number: GMK500KE Catalog number: NAK500KE Catalog number: UMK500KE Catalog number: HTMK500KE

For more information of drug targets and enzyme assays, please visit www.profoldin.com.