

**ProFoldin**

10 Technology Drive, Suite 40, Number 188
Hudson, MA 01749-2791 USA
Tel: (508) 735-2539 FAX: (508) 845-9258
www.profoldin.com
info@profoldin.com


INSTRUCTIONS

ProFoldin *E. coli* DNA Gyrase

***E. coli* DNA Gyrase – for 100 assays**

Catalog No. TOP2-100EC

Protein construct:	Wild-type <i>E. coli</i> DNA gyrase, the (gyrA) ₂ (gyrB) ₂ complex, composed of the gyrA subunit (MW 97 kDa) and gyrB subunit (90 kDa) purified from a bacterial expression system.
MW:	374 kDa
Enzyme concentration:	2 μM
Enzyme activity assay:	The DNA gyrase DNA supercoiling activity is measured by using the DNA Topoisomerase II (Gyrase) Assay Kit (Catalog No. DSA020K)
Storage temperature:	-20 or -80°C. Do not freeze-and-thaw repeatedly.
Enzyme dilution:	Use the 1 x assay to dilute the enzyme just before the assay. Do not store diluted enzyme solution



The *E. coli* DNA Topoisomerase II (Gyrase) – for 100 assays (Catalog No. TOP2-100EC) 50 μl
100 x *E. coli* gyrase (2 μM). It is for 100 assays.

Assay Protocol using DNA Topoisomerase II (Gyrase) Assay Kit

The following assay protocol is based on a 96-well assay plate format using a standard black 96-well plate (Greiner 655076). For 384-well plate assays using a standard black 384-well plate (Matrix 4318), please reduce the reagent volumes proportionally.

1. Reaction:

The total volume of each reaction mixture is 40 μl including: 24 μl of H₂O, 4 μl of 10 x buffer, 4 μl of 10 x relaxed DNA (250 μg/ml relaxed plasmid DNA), 4 μl of 10 x enzyme, 4 μl of 10 mM ATP. Incubate the reaction mixture at 37°C for 60 min.

**ProFoldin**

10 Technology Drive, Suite 40, Number 188
Hudson, MA 01749-2791 USA
Tel: (508) 735-2539 FAX: (508) 845-9258
www.profoldin.com
info@profoldin.com

INSTRUCTIONS

Note: The final concentrations are 20 mM Tris-HCl, pH 8, 35 mM NH₄OAc, 4.6 % glycerol, 1 mM DTT, 0.005% Brij35, 8 mM MgCl₂, 25 µg/ml relaxed plasmid DNA, 1 mM ATP and 20 nM topoisomerase II. The 10 x enzyme is prepared by dilution of the 100x enzyme in 1x assay buffer. A negative control reaction can be the reaction mixture without addition of ATP. Magnesium is essential for the reaction and assay. EDTA should be avoided.

2. Assay

- (1) Make 1 x H19 Dilution Buffer by dilution of the 10 x H19 Dilution Buffer 10-fold with water.
Freshly prepare the H19 dye by dilution of 1 µl the 1500 x H19 dye stock solution with 1.5 ml of 1 x H19 Dilution Buffer (1500 x dilution). Mix the solution evenly by inverting the tube a few times.
- (2) Mix 250 µl of the freshly prepared H19 dye with each reaction solution (40 µl). Incubate the mixture at room temperature for 5 min.
- (3) Measure the fluorescence intensity at 535 nm using the excitation wavelength at 485 nm.

Related Products:

Human Topoisomerase I, 10,000 Units	HTOPI-010
DNA Topoisomerase II (Gyrase) Assay Kit	DSA020K
<i>E. coli</i> DNA Topoisomerase II (Gyrase) Assay Kit Plus-100	DSA100KE
Gel-based <i>E. coli</i> DNA Topoisomerase II (Gyrase) Assay Kit Plus-100	GDSA100KE
Relaxed plasmid DNA - 0.050 mg	RDNA-050UG
Supercoiled plasmid DNA - 0.050 mg	SDNA-050UG
H19 Dye for DNA Relaxation and Supercoiling Assays	DSA1000D
Human Topoisomerase I DNA Relaxation Assay Kit Plus -100	HRA100KE
Human Topoisomerase II DNA Decatenation Assay Kit Plus-100	HDC100KE
<i>E. coli</i> DNA Topoisomerase I Assay Kit Plus-100	DRA100KE
<i>E. coli</i> gyrase ATPase assay Kit Plus	T2A-100KE
96-well Topoisomerase DNA Cleavage Assay Kit	TDC96K
96-Well Topoisomerase DNA Decatenation Assay Kit	TDD96K
96-Well <i>S. aureus</i> Topo IV DNA Decatenation Assay Kit Plus	SDD96KE

Publications

1. Asha M.K. et al, In vitro anti-Helicobacter pylori activity of a flavonoid rich extract of Glycyrrhiza glabra and its probable mechanisms of action, Journal of Ethnopharmacology, Vol.145(2), pp.581-586 (2013).
2. Mora-Pale M et al, Antimicrobial mechanism of resveratrol-trans-dihydrodimer produced from peroxidase-catalyzed oxidation of resveratrol. Biotechnol Bioeng. Vol 112, pp2417-2428 (2015).

For more information of DNA topoisomerase assays and assays for more drug targets and enzymes, please visit www.profoldin.com or send emails to info@profoldin.com.