



ProFoldin


10 Technology Drive, Suite 40, Number 188
Hudson, MA 01749-2791 USA
Tel: (508) 735-2539 FAX: (508) 845-9258
www.profoldin.com
info@profoldin.com

INSTRUCTIONS

ProFoldin *S. aureus* DNA Gyrase

S. aureus DNA Gyrase – for 100 assays Catalog No. TOP2-100SA

Protein construct:	Wild-type <i>S. aureus</i> DNA gyrase, the (gyrA) ₂ (gyrB) ₂ complex, composed of the gyrA subunit (MW 99 kDa) and gyrB subunit (73 kDa) purified from a bacterial expression system.
MW:	344 kDa
Enzyme concentration:	7.5 μM
Enzyme activity assay:	The DNA gyrase DNA supercoiling activity is measured by using the <i>S. aureus</i> DNA Topoisomerase II (Gyrase) Assay Kit (Catalog No. DSA100KS).
Storage temperature:	-20 or -80°C. Do not freeze-and-thaw repeatedly.
Enzyme dilution:	Use the 1 x assay to dilute the enzyme just before the assay. Do not store diluted enzyme solution



The *S. aureus* DNA Gyrase – for 100 assays (Catalog No. TOP2-100SA) includes 43 μl 7.5 μM *S. aureus* gyrase. It is for 100 assays in a 96-well plate format.

Assay Protocol using the *S. aureus* DNA Topoisomerase II (Gyrase) Assay Kit

The following assay protocol is based on a 96-well assay plate format using a standard black 96-well plate (Greiner 655076). For 384-well plate assays using a standard black 384-well plate (Matrix 4318), please reduce the reagent volumes proportionally.



ProFoldin

10 Technology Drive, Suite 40, Number 188
Hudson, MA 01749-2791 USA
Tel: (508) 735-2539 FAX: (508) 845-9258
www.profoldin.com
info@profoldin.com

INSTRUCTIONS

1. Reaction:

The total volume of each reaction mixture is 40 µl including: 19.6 µl of H₂O, 4 µl of 10 x buffer, 8 µl of 2 M potassium glutamate, 4 µl of 10 x relaxed DNA (250 µg/ml), 0.4 µl of 7.5 µM *S. aureus* gyrase, 4 µl of 10 mM ATP. Incubate the reaction mixture at 37°C for 60 min.

Note: The final concentrations are 20 mM Tris-HCl, pH 8, 35 mM NH₄OAc, 4.6 % glycerol, 1 mM DTT, 0.005% Brij35, 8 mM MgCl₂, 400 mM potassium glutamate, 25 µg/ml relaxed plasmid DNA, 1 mM ATP and 75 nM topoisomerase II. A negative control reaction can be the reaction mixture without addition of ATP. Magnesium is essential for the reaction and assay. EDTA should be avoided.

2. Assay

- (1) Make 1 x H19 Dilution Buffer by dilution of the 10 x H19 Dilution Buffer 10-fold with water. Prepare the H19 dye by dilution of 1 µl the 1500 x H19 dye stock solution with 1.5 ml of 1 x H19 Dilution Buffer (1500 x dilution). Mix the solution evenly by inverting the tube a few times.
- (2) Mix 250 µl of the freshly prepared H19 dye with each reaction solution (40 µl). Incubate the mixture at room temperature for 5 min.
- (3) Measure the fluorescence intensity at 535 nm using the excitation wavelength at 485 nm.

Related Products:

<i>S. aureus</i> DNA Topoisomerase II (Gyrase) Assay Kit Plus-100	DSA100KSE
96-Well <i>S. aureus</i> Topo IV DNA Decatenation Assay Kit Plus	SDD96KE
Human Topoisomerase I, 10,000 Units	HTOPI-010
<i>E. coli</i> DNA Topoisomerase II (Gyrase) Assay Kit Plus-100	DSA100KE
Relaxed plasmid DNA - 0.050 mg	RDNA-050UG
Supercoiled plasmid DNA - 0.050 mg	SDNA-050UG
H19 Dye for DNA Relaxation and Supercoiling Assays	DSA1000D
Human Topoisomerase I DNA Relaxation Assay Kit Plus -100	HRA100KE
Human Topoisomerase II DNA Decatenation Assay Kit Plus-100	HDC100KE
<i>E. coli</i> gyrase ATPase assay Kit Plus	T2A-100KE
96-well Topoisomerase DNA Cleavage Assay Kit	TDC96K

Publications

1. Asha M.K. et al, In vitro anti-Helicobacter pylori activity of a flavonoid rich extract of Glycyrrhiza glabra and its probable mechanisms of action, Journal of Ethnopharmacology, Vol.145(2), pp.581-586 (2013).
2. Mora-Pale M et al, Antimicrobial mechanism of resveratrol-trans-dihydrodimer produced from peroxidase-catalyzed oxidation of resveratrol. Biotechnol Bioeng. Vol 112, pp2417-2428 (2015).

For more information of DNA topoisomerase assays and assays for more drug targets and enzymes, please visit www.profoldin.com or send emails to info@profoldin.com.