


# Purification Columns

MOBICOLS | LAB COLUMNS

## Clean-Up

of Proteins, Antibodies, Nucleic Acids & More



Request a  
free sample!



Enjoy short  
processing  
times!



Works with small  
quantities of  
resin & buffer!

**Mobicols**

➤ Unique Empty or Pre-filled Mini Columns

**Lab Columns**

➤ Empty Chromatography Columns  
with Luer-Lock Connections



[www.mobitec.com](http://www.mobitec.com)



Published by MoBiTec GmbH, Goettingen, Germany ©2022

All rights reserved.

# Purification Columns

## MOBICOLS

- ▶ Mini columns (up to 800 µl) - Ideal for small sample volumes
- ▶ Fast and simple purification of biomolecules
- ▶ Practical tool for chemical research, e.g., for chromatographic intragroup lanthanide separation
- ▶ Streamlines the use of resins
- ▶ Gravity-flow columns or spin columns
- ▶ Useful for, e.g., size exclusion or affinity chromatography
- ▶ Empty or pre-packed columns
- ▶ Choose from various filter sizes and accessories
- ▶ Luer-fittings enable connection to various chromatography systems or elution using a standard syringe
- ▶ Small scale laboratory techniques
  - Allows working with small quantities of resin
  - Markedly reduced buffer consumption and processing time
- ▶ Renowned for quality:
  - MoBiTec invented and designed the columns (patented)
  - Highest quality for more than 30 years, Made in Germany
  - Our mini columns have already been sold several million times worldwide. They are an integral part of many laboratory diagnostics kits for clinical routine & research



## Lab Columns – Empty Chromatography Columns

- ▶ Empty chromatography columns for lab scale
- ▶ Syringe-based columns with Luer-Lock adapters at the in- and outlet, providing several handling advantages



# Quick Access: Mobicols & Lab Columns



**BENEFITS**

- Simple and intuitive purification systems for proteins, antibodies, carbohydrates, nucleic acids & more
- From microgram to multiple gram range
- Providing speed, quality, and versatility in research and lab-scale process development
- Increasing efficiency in method development
- Reduction of processing times
- Minimization of error sources

**QUICK ACCESS:**

## MOBICOLS & LAB COLUMNS

**MOBICOLS – Mini Columns, up to 800 µl**

Incredibly versatile chromatography columns. Allows working with small quantities of resin, markedly reduced buffer consumption & processing time.

**UNLIMITED USE (see references)**

Tandem affinity purification • Lanthanide separation • Clotting of intact haemocytes • Desalting of peptides before mass spectrometry analysis • Immobilization of carbohydrates • Identification of RNA aptamers • Isolation of different RNA types • PCR product purification before amplicon sequencing • Microscale immunoprecipitation • Whole-mount in situ hybridization • Extraction of Sialic-acid-containing Glycans and Glycopeptides using Seroitin-bonded silica • Isolation of proteins from polyacrylamide gels • Purification of nuclear signal proteins • Proteolytic enrichment • In vitro pull-down assays for studying protein interactions • Embryonic immunostaining • Desalting of carbohydrates • Immobilization of small molecule drugs • Drug-affinity chromatography • etc.

**Empty Mobicols**

Mini columns can be packed with a variety of chromatographic media. Simply fill in the matrix of choice and create your own purification system.

**Applications**

- Clean-up of proteins, antibodies, nucleic acids & more
- Useful for size exclusion, affinity or ion exchange chromatography
- Featured by a wide range of accessories to customize the columns to specific needs

**Mobicol Gravi**  
Gravity-Flow Columns

Luer-Lock fittings enable connection to many chromatography systems or Luer-Lock syringes; procedure extremely robust, adjustments can be made during the run.

**Classic Variant**  
Bottom plug: separate and manually insertable plug at the outlet (e.g., #M1002).

**Fixed-Type**  
Fixed snap-off plug (leak-proof at varying pressures), reverse side can be used for reclosing (e.g., #M1005).

**Mobicol Spin**  
Spin Columns

Rapid purification by centrifugation; easy handling: load sample, spin, and collect the purified product. Works with standard microcentrifuge.

**Classic Variant**  
Bottom plug: separate and manually insertable plug at the outlet (e.g., #M1005).

**Fixed-Type**  
Fixed snap-off plug (leak-proof at varying pressures), reverse side can be used for reclosing (e.g., #M1005).

**Pre-Filled Mobicols**

Mini chromatography cartridges pre-packed with different resins for separation and purification of various biomolecules

Ready-to-use • Fast procedure • No sample dilution

**Mobicol Spin S**

Fast and reproducible purification of nucleic acids, proteins, polysaccharides, and other macromolecules by size exclusion chromatography.

**Applications**

- Nucleotide removal
- Plasmid & PCR reaction purification
- Removal of unincorporated dyes
- Purification of virus RNA

**Mobicol Spin G**

Includes proven gel filtration resin. Quickly desalts, removes contaminants, and transfers to a new buffer in a single step.

**Applications**

- Buffer exchange
- Nucleotide removal
- Removal of dye terminators
- Detachment of labeling tags

**Mobicol Spin Desalting**

Well-ried silica matrix separates molecules by size exclusion chromatography. Confident handling: silica resins are allowed to fall dry without any loss in purification efficiency. Improved recovery of small-sized nucleic acid fragments.

**Applications**

- Purification of nucleic acids from small contaminants
- Desalting and buffer exchange, especially of small-sized nucleic acid fragments and dNTPs

**LAB COLUMNS – Empty Chromatography Columns**

Easy-to-handle gravity-flow columns. Compatible with any chromatography resin. Syringe-based columns with Luer-Lock connections at both ends. Come with lower filters with 10 µm, 35 µm, or 90 µm pore size.

**2.5 ml Lab Columns**  
FPLC compatible

**5.0 ml Lab Columns**

**10 ml Lab Columns**

**Column Sets**

Mobicol (up to 800 µl)  
2.5 ml Lab Columns  
5 ml Lab Columns  
10 ml Lab Columns

**Applications**

- Affinity chromatography
- Size exclusion chromatography
- Ion exchange chromatography
- E.g., purification of recombinant proteins or plasmid clean-up

**Mobicols with Immobilized Enzymes**

Enzyme immobilization provides enzyme reutilization and prevents the loss of enzyme activity. Immobilized systems provide better enzyme thermostability and pH tolerance. The method shortens redundant downstream and purification processes.

**Applications**

- Dephosphorylation of vector DNA
- Removal of DNases and RNases
- Protein digestions for peptide mapping & protein sequencing

**Mobicol Spin Ni-IDA**

Protein purification by affinity chromatography. Improved target specificity by optimized silica-based Ni<sup>2+</sup>-IDA matrix. Eliminates almost all non-specific interactions of contaminating host proteins with the adsorbent.

**Applications**

- Purification of recombinant polyhistidine-tagged proteins
- Suitable for small proteins, large protein complexes, proteins with low expression rates

**Customized Services**

Custom Filling: Mini Columns packed with the matrix of your choice. (On-source labor-intensive step!)

**Applications**

- Life Science Research
- Laboratory Diagnostics
- Food Safety

Published by MoBiTec - GmbH, Geretsried, Germany © 2022 All rights reserved. [www.mobitec.com](http://www.mobitec.com)



[Download detailed PDF](#)

## Select Citations



[Download Reference List](#)

# Topics in this Issue

▶ Quick Access: Mobicols & Lab Columns	4
▶ Introduction	6
• Purification Columns – Clean-Up of Proteins, Antibodies, Nucleic Acids & More	6
• Mobicols – Versatile Mini Columns	6
• Lab Columns – Empty Chromatography Columns	7
▶ Empty Mobicols	8
• Mobicol Gravi – Empty Gravity-Flow Columns	8
• Mobicol Spin – Empty Spin Columns	10
▶ Mobicol Accessories	12
• Filters	12
• Adapters, Plugs, and More (for Gravity-Flow)	13
▶ Pre-filled Mobicols	15
• Mobicol Spin S-Columns for Nucleic Acid Purification	15
• Mobicol Spin G-50 Columns for Nucleic Acid Purification	18
• Mobicol Spin Desalting Columns – Separation of Small Contaminants	19
• Ni-IDA Columns for Purification of His-tagged Proteins	20
• Mobicols with Immobilized Enzymes	22
▶ Lab Columns – Empty Chromatography Columns	24
• Empty Chromatography Columns for Lab Scale with Luer-Lock Connections	24
• Column Sets	26
▶ Customized Services	27
• Custom Filling: Mini Columns Packed with the Matrix of your Choice	27
▶ Reference List	4

# Introduction

## Purification Columns – Clean-Up of Proteins, Antibodies, Nucleic Acids & More

The purification of biomolecules is an essential part of the daily laboratory routine. Protein purification is vital for the specification of the function, structural studies, and biochemical assays. For example, purified antibodies are essential for many research applications. High purity DNA and RNA is a prerequisite for most molecular biological applications, including enzymatic nucleic acid digestion, transformation, RT-(q)PCR, and NGS library preparation. Our multipurpose chromatography columns make clean-up protocols easier and faster.

MoBiTec offers a variety of columns for the purification of nucleic acids, peptides, proteins, and antibodies. Depending on the matrix, they can be used for affinity, size exclusion, and ion exchange chromatography. They are suitable for small sample volumes of a few  $\mu\text{l}$  up to larger volumes of several ml. Furthermore, they are featured by a wide range of accessories to customize the columns to your specific needs.

Save valuable resources and enhance your purification power with our Mobicols (mini columns) or our Lab Columns (empty chromatography columns).

## Mobicols – Versatile Mini Columns

Mobicols are incredibly versatile mini columns (up to 800  $\mu\text{l}$ ) that are an indispensable addition to every lab. Mobicols can be used both 1) as gravity-flow column (**Mobicol Gravi**) or 2) as spin column (**Mobicol Spin**). Depending on your preference, they can be ordered as empty columns or pre-filled with matrix, e.g., Sephadex® G-50 resin, different Sephacryl® resins, or silica-based resins. Luer-Lock adapters and filters with different diameters and pore sizes are offered to customize the columns to your specific needs. Mobicols are made of polypropylene, the filters of polyethylene. The mini columns can be autoclaved at 120 °C (with inserted filter only at 110 °C for 10 minutes).



### Shared Features of Mobicol Gravi & Mobicol Spin Columns:

- ▶ Fast and simple purification of biomolecules like nucleic acids, peptides, proteins, and other compounds
- ▶ Useful for size exclusion, affinity, or ion exchange chromatography
- ▶ Empty columns - fill in the matrix of choice
- ▶ Also available pre-filled with different matrices
- ▶ Ideal for small sample volumes
- ▶ Plugs included for relocking
- ▶ Filters available with 10  $\mu\text{m}$ , 35  $\mu\text{m}$ , and 90  $\mu\text{m}$  pore sizes
- ▶ Materials: polypropylene (column corpus), polyethylene (filters)
- ▶ Autoclavable at 120 °C (with inserted filter at 110 °C for 10 minutes)

# Introduction

**Additional Feature Mobicol Gravi:**

- ▶ Luer-Lock cap which is compatible to Luer-Lock syringes or other chromatography systems

**Additional Feature Mobicol Spin:**

- ▶ Can be used with a standard microcentrifuge (can be centrifuged in most 1.5 ml or 2 ml collection tubes)

## Lab Columns – Empty Chromatography Columns

These empty syringe-based columns are available in three sizes (2.5 ml, 5 ml, and 10 ml) and come with an inserted lower filter with the pore size of your choice (10  $\mu\text{m}$ , 35  $\mu\text{m}$ , or 90  $\mu\text{m}$ ). Each column has Luer-Lock connectors at both ends for chromatography tubing or syringes. To close the columns, the 5 ml and 10 ml columns include an additional cap and the 2.5 ml columns have a Luer-Lock closing cap. The 2.5 ml columns can be centrifuged and autoclaved (with inserted filters only at 110 °C for 10 minutes). Each package of columns includes two One-Way stopcocks.

Upper filters (used on top of the matrix) and Luer-Lock adapters are available separately for the 2.5 ml, 5 ml, and 10 ml columns.

**Features:**


- ▶ Easy to handle
- ▶ Smallest and largest volumes can be treated straightforwardly
- ▶ Compatible with laboratory standards
- ▶ Each column has Luer-Lock connectors at both ends for chromatography tubing or syringes
- ▶ Two One-Way Stopcocks included

**Applications:**

- ▶ Affinity chromatography
- ▶ Size exclusion chromatography
- ▶ Ion exchange chromatography



# Empty Mobicols



Empty Mobicols are disposable columns that can be packed with a variety of chromatographic resins. You can choose between gravity-flow or spin columns.

## Mobicol Gravi – Empty Gravity-Flow Columns

### Mobicol Gravi – Clean-up by Gravity-Flow Chromatography

Column chromatography can be performed with automated systems which use a pump to force solvent through a packed column at a set flow rate, or can be run by gravity-flow. The gravity method is less expensive, the user has more control, and can make adjustments during the run.

Mobicol Gravi columns are developed for performing size exclusion, affinity, or ion exchange chromatography by gravity-flow. By using the supplied Luer-Lock cap, there is no limitation in volume, i.e., it

is possible to apply small or large volumes of sample or buffer. Using the Luer-Lock cap, a syringe can be connected to the top of the Mobicol Gravi column as a reservoir for sample or buffer. To enable a connection to larger reservoirs, an additional Luer-Lock adapter is available to attach chromatography tubing to the top cap. This even enables a safe connection from a pump to the inlet of a Mobicol Gravi column. Luer-Lock connections for the bottom, a One-Way Stopcock and various filters are also offered to customize the columns to specific needs (Mobicol Accessories, page 13).



Left: large filter, right: small filter



MoBiTec offers small and large filters of different diameters. For gravity-flow chromatography, we strongly recommend using large filters. Otherwise, the flow through the column is too slow, not uniform, and can also break off completely. You can also choose from various filter pore sizes - 10  $\mu\text{m}$ , 35  $\mu\text{m}$ , and 90  $\mu\text{m}$ . The right filter pore size depends on your matrix: the matrix used should not pass through the filter pores.

# Empty Mobicols

## Mobicol Gravi – Features:

- ▶ Mini columns (up to 800 µl) for quick purification of biomolecules and other compounds
- ▶ Fill with the matrix of your choice
- ▶ Supplied with two different caps: screw cap with rubber seal and Luer-Lock cap compatible to, e.g., syringes with Luer-Lock fitting
- ▶ Additional adapters are available for inlet and outlet chromatography tubing
- ▶ Plugs included for relocking
- ▶ Filters with 10 µm, 35 µm, and 90 µm pore sizes available
- ▶ Autoclavable at 120 °C (with inserted filter at 110 °C for 10 minutes)
- ▶ Materials: polypropylene (column corpus), polyethylene (filters)
- ▶ Available as Classic or Fixed type (see box on the right)



### Different Design – Same Application: Mobicol Gravi Classic vs. Mobicol Gravi Fixed

Two types of Mobicol Gravi columns are available: **Mobicol Gravi Classic** and **Mobicol Gravi Fixed**. Both column types differ slightly in shape, but are totally equal regarding sample processing properties.

The Classic variant has a separate and manually insertable plug for the outlet (bottom plug), whereas the Fixed type comes as one piece with a fixed snap-off plug. The Mobicol Gravi Fixed can be packed, stored, or shipped with no leakage. The outlet is 100% tight until bent off. The reverse side of the snap-off plug can then be used for relocking.

**Mobicol Gravi Fixed:**  
Leak-proof outlet plug; tight even at varying pressures (e.g., transportation in airplanes).

## Overview Empty Gravity-Flow Columns

Order No.	Description	Amount
 <b>M1002</b>	Mobicol Gravi Classic with Luer-Lock and screw cap, bottom plug, without filters	50
 <b>M1050</b>	Mobicol Gravi Fixed with Luer-Lock and screw cap, fixed outlet plug, without filters	50
 <b>M105035F</b>	Mobicol Gravi Fixed with Luer-Lock and screw cap, fixed outlet plug and inserted large 35 µm filter	50

**Bottom plug:** separate and manually insertable plug for the outlet.

**Fixed outlet plug:** fixed snap-off plug, reverse side can be used for relocking.


**Filters and other accessories** can be ordered separately.

# Empty Mobicols



## Mobicol Spin – Empty Spin Columns

### Mobicol Spin – Speed-Up your Purification with a Microcentrifuge



Mobicol Spin Columns can be centrifuged in 1.5 ml or 2 ml collection tubes

The Mobicol Spin Columns are developed as fast and convenient routine tool for clean-up of biomolecules by spinning. They can be used without matrix (e.g., to extract trypsinated proteins out of a gel slice for Mass Spec analysis) or filled with a matrix of your choice, e.g., for performing size exclusion or affinity chromatography.

Independent of the preferred chromatography, the handling is easy: load the sample on the column, spin, and collect the flow-through. For collection, the Mobicol Spin column must be placed in a 1.5 ml or 2 ml microtube before centrifugation. Compared to gravitation chromatography, spinning is faster.



Left: large filter, right: small filter

**👍** For spinning, both small and large filters (different diameters) can be selected. If you choose the smaller filters, you will have a larger working volume (minimal death volume). This is due to the fact that the smaller filters are located further down in the column body.

You can also choose from various filter pore sizes – 10 µm, 35 µm, and 90 µm. The appropriate filter pore size depends on your matrix: the matrix used should not pass through the filter pores.

# Empty Mobicols

## Mobicol Spin – Features

- ▶ Mini spin columns (up to 800 µl) for use in a standard microcentrifuge
- ▶ Comes empty – fill in your own material (matrix, gel slice, etc.)
- ▶ Suitable for the purification of nucleic acids, peptides, proteins, and more
- ▶ Easy handling: load sample, spin, and collect the purified product
- ▶ Less than 4 minutes per application
- ▶ Reproducible results
- ▶ No sample dilution
- ▶ Versatile in use



### Different Design – Same Application: Mobicol Spin Classic vs. Mobicol Spin Fixed

Two types of Mobicol Spin columns are available: **Mobicol Spin Classic** and **Mobicol Spin Fixed**. Both column types differ slightly in shape, but are totally equal regarding sample processing properties.

The Classic variant has a separate and manually insertable plug for the outlet (bottom plug), whereas the Fixed type comes in one piece with a fixed snap-off plug. The Mobicol Spin Fixed can be packed, stored, or shipped with no leakage. The outlet is 100% tight until bent off. The reverse side of the snap-off plug can then be used for relocking.

**Mobicol Spin Fixed:**  
Leak-proof outlet plug; tight even at varying pressures (e.g., transportation in airplanes).

## Overview Empty Spin Columns

Order No.	Description	Amount
 <b>M1003</b>	Mobicol Spin Classic with screw cap, bottom plug, without filters	50
 <b>M1053</b>	Mobicol Spin Fixed with screw cap, fixed outlet plug, without filters	50
 <b>M105010S</b>	Mobicol Spin Fixed with screw cap, fixed outlet plug and inserted small 10 µm filter	50
 <b>M105210S</b>	Mobicol Spin Fixed with screw cap, fixed outlet plug and inserted large 10 µm filter	50
 <b>M1002S</b>	Mobicol Spin Classic with screw cap, bottom plug and large 10 µm filter	100

**Bottom plug:** separate and manually insertable plug for the outlet.

**Fixed outlet plug:** fixed snap-off plug, reverse side can be used for relocking.

**Filters and other accessories** can be ordered separately.



For its mini columns MoBiTec offers a variety of filters with different pore sizes and diameters. With this range, the users can select the ideal size for their applications for optimal performance. Multiple accessories designed for gravity-flow chromatography are also available, including Luer-Lock adapters and stopcocks.

## Filters

Filters with different diameters and pore sizes are offered separately to customize the Mobicol columns to your specific needs.

Order No.	Description	Amount	Recommended for	Can be used with (Mobicol Order No.)
 <b>M2110</b>	Filter (small) 10 µm pore size	50	Spinning	M1003, M1053
 <b>M513515</b>	Filter (small) 35 µm pore size	50	Spinning	M1003, M1053
 <b>M2190</b>	Filter (small) 90 µm pore size	50	Spinning	M1003, M1053
 <b>M2210</b>	Filter (large) 10 µm pore size	50	Gravity-flow and spinning	M1002, M1003, M1050 & M1053
 <b>M523515</b>	Filter (large) 35 µm pore size	50	Gravity-flow and spinning	M 1002, M1003, M1050 & M1053
 <b>M2290</b>	Filter (large) 90 µm pore size	50	Gravity-flow and spinning	M1002, M1003, M1050 & M1053

For gravity-flow chromatography, we strongly recommend using large filters. For spinning, both filter sizes (small and large) can be selected. The appropriate filter pore size depends on your matrix: the matrix used should not pass through the filter pores.

## Adapters, Plugs, and More (for Gravity-Flow)




**M3001** Tubing-to-Luer-Lock Adapter



**M3003** Luer-Lock Adapter M



**M3004** Luer-Lock Adapter F




**M3005** Bottom Plugs



**M3006** One-Way Stopcock



**M3007** Luer-Lock Adapter M/M



**M3009** Luer-Lock Caps

The Luer-Lock adapters at the inlet and outlet of the Mobicol Gravi columns enable the purification of larger sample volumes that exceed the actual capacity of the mini column, or even partial automation of the process. For this purpose, MoBiTec offers a variety of adapters to enable the connection to larger reservoirs or even to attach chromatography tubing to the top cap. Without any additional adapter, a syringe can be connected to the Luer-Lock cap. The polypropylene Luer-Lock Adapter M (M3003) for tubing connection to the inlet of the column

allows a flow-through mode at low pressure, whereas the Luer-Lock Adapter F (M3004) is suitable for a tubing connection to the outlet of the column. Furthermore, the Luer-Lock Adapter M/M (M3007) facilitates the connection of the inlet to a Luer-Lock system. The Tubing-to-Luer-Lock Adapter (M3001) enables a safe connection from a pump to the inlet of a Mobicol. The One-Way Stopcock (Luer-Lock F; M3006) is a convenient tool for controlling flow-through at the in- or outlet.

Order No	Description	Amount
<b>M3001</b>	Tubing-to-Luer-Lock Adapter (material outside: metal, inside: glass, with 1 m tubing, for Mobicols or Lab Columns)	1
<b>M3002</b>	Set: Luer-Lock Adapter M, Luer-Lock Adapter F (10 each)	20
<b>M3003</b>	Luer-Lock Adapter M (for Luer-Lock cap)	20
<b>M3004</b>	Luer-Lock Adapter F (for outlet)	20

Order No	Description	Amount
<b>M3005</b>	Bottom Plugs (for Mobicol Gravi Classic or Mobicol Spin Classic)	50
<b>M3006</b>	One-Way Stopcock (Luer-Lock F)	4
<b>M3007</b>	Luer-Lock Adapter M/M	1
<b>M3009</b>	Luer-Lock Caps	50

# Pre-filled Mobicols



## MoBiTec offers different pre-filled mini columns:


- ▶ Spin columns for purification of nucleic acids by size exclusion chromatography
- ▶ Spin columns for purifications of proteins by affinity chromatography
- ▶ Mini columns (for spinning or gravity-flow) with enzymatically active matrices (immobilized enzymes)

## Setup Pre-Filled Mobicols

The pre-packed mini column is closed with a screw cap and contains a specific matrix. The outlet of the mini column is a fixed snap-off plug. Prior to the first centrifugation step, it has to be removed by bending. The reverse side of the plug can be used for relocking.

### Features

- ▶ Columns are pre-packed, equilibrated, and ready-to-use
- ▶ Compatible with laboratory standard
- ▶ Quick procedure: sample purification in less than 4 minutes
- ▶ Easy handling: spin, load sample, spin, and collect the purified product
- ▶ Reproducible results with simplified protocols
- ▶ Numerous samples can be processed simultaneously
- ▶ No sample dilution



## Mobicol Spin S-Columns for Nucleic Acid Purification

### Applications

- ▶ Buffer exchange
- ▶ Nucleotide removal
- ▶ Plasmid, oligonucleotide, and PCR reaction purification
- ▶ Removal of unincorporated dyes or dye terminators


### Description

Mobicol Spin S-columns are pre-packed with S-200, S-300, or S-400 Sephacryl® HR matrices. Sephacryl® High Resolution chromatography resins are highly versatile resins that offer a wide range of fractionation capabilities. They allow fast and reproducible purification of proteins, polysaccharides, and other macromolecules by size exclusion

chromatography. Excellent resolution and flow characteristics as well as long-term physical and chemical stability make Sephacryl® HR the resin of choice for routine purification.

The Mobicol Spin S-columns separate molecules by differences in size. The columns combine the effectiveness of gel filtration with speed of centrifugation. The pore size of the filled-in matrix determines which molecules are small enough to enter the pores of the matrix beads and which molecules are too large to enter. In this way, compounds larger than the pore size are excluded from the resin and located in the void volume. These larger molecules move quickly through the matrix bed when the column is centrifuged. Molecules smaller than the pore size, like hydrated salt ions, do enter the pores of the matrix beads and are held back.

*Sephacryl® is a registered trademark of Cytiva.*



### Size Exclusion Chromatography with Mobicol Spin Columns

Molecules larger than the matrix pore size cannot enter the matrix beads. They appear only in the matrix void volume between the beads. From the inter-bead space, they can be recovered by spinning the column in a collection tube in a benchtop centrifuge. Strong spinning elutes the column void volume without dilution. Molecules smaller than the matrix pore size and hydrated salt ions locate inside the matrix beads and are not eluted by spinning.

# Pre-filled Mobicols

## Benefits of Size Exclusion Chromatography (SEC)

SEC is a liquid chromatography technique which separates molecules by virtue of their size in solution. The size separation step occurs in a column filled with porous particles (stationary phase) and is an entropy-driven process.

Molecules do not bind to the stationary phase. Instead, they are separated by the speed at which they navigate through the inert chromatography resin. Therefore, buffer composition does not directly affect resolution. Conditions can thus be varied to suit the requirements of your sample or downstream processing.

SEC offers mild conditions, simple operation, isocratic elution, and easy scaling up.

SEC is perfectly suited for the separation of biomolecules sensitive to changes in pH, concentration of metal ions or cofactors, or harsh environmental conditions. Separations can be performed in the presence of essential ions, cofactors, detergents, urea, or guanidine hydrochloride at high or low ionic strength.

SEC can be performed within a wide temperature range, provides short analysis times, and avoids sample loss.

## Sephacryl®-Based Columns

Order No	Description	Amount
 <b>SCO200</b>	Mobicol Spin S-200	20 columns
 <b>SCO210</b>	Mobicol Spin S-200	100 columns
 <b>SCO300</b>	Mobicol Spin S-300	20 columns
 <b>SCO310</b>	Mobicol Spin S-300	100 columns
 <b>SCO400</b>	Mobicol Spin S-400	20 columns
 <b>SCO410</b>	Mobicol Spin S-400	100 columns
 <b>SCO234</b>	Mobicol Spin S-200, S-300, S-400	3x 10 columns

Columns are pre-packed with Sephacryl® S-200 HR, Sephacryl® S-300 HR, or Sephacryl® S-400 HR resin and equilibrated in Mobicol Spin Buffer (10 mM Tris/HCl pH 7.6; 1 mM EDTA).

Sephacryl® HR resins are available in different fractionation ranges. This allows for reproducible separations over a wide molecular weight range:

- **Sephacryl® S-200 High Resolution Resin** is designed for purifying antibodies, serum proteins, and mid-size proteins (fractionation range globular proteins  $5 \times 10^3 - 2.5 \times 10^5$ )

- **Sephacryl® S-300 High Resolution Resin** is designed for purifying antibodies, serum proteins, and mid-size proteins (fractionation range globular proteins  $1 \times 10^4 - 1.5 \times 10^6$ )
- **Sephacryl® S-400 High Resolution Resin** is recommended to separate polysaccharides, macromolecules with extended structures, and even small particles such as plasmids (fractionation range globular proteins  $2 \times 10^4 - 8 \times 10^6$ )

Data provided by Cytiva;  
Sephacryl® is a registered trademark of Cytiva.

# Pre-filled Mobicols

## Application Example – Food Safety

### Mobicol Spin S-400 – Recommended by the German Federal Institute for Risk Assessment (BfR) for the Detection of Norovirus in Food



Mobicol Spin S-400: Proven for Norovirus RNA Purification

## Publication Highlight

- ▶ Frozen strawberries can contain different amounts of RT-PCR inhibitors
- ▶ Efficiency of norovirus extraction protocols from strawberries varies remarkably
- ▶ Additional RNA purification with Sephacryl®-based columns (Mobicol Spin S-400) increases detection rate
- ▶ Use of these columns is suggested for norovirus detection in berry samples

*Bartsch C, Szabo K, Dinh-Thanh M, Schrader C, Trojnar E, John R. Comparison and optimization of detection methods for noroviruses in frozen strawberries containing different amounts of RT-PCR inhibitors. Food Microbiol. 2016;60:124-130. doi:10.1016/j.fm.2016.07.005*

# Pre-filled Mobicols



## Mobicol Spin G-50 Columns for Nucleic Acid Purification

### Applications

- ▶ Buffer exchange
- ▶ Nucleotide removal
- ▶ Removal of unincorporated dyes or dye terminators

### Description

Mobicol Spin G-50 columns pre-packed with Sephadex® G-50 resin remove small molecules from DNA/RNA, such as salt ions, and are therefore most suitable for buffer exchange (even from dsDNA fragments as small as 20 bp). They can also be used to remove fluorescent dye dideoxy terminators (e.g., Cy3/Cy5 nucleotides) from cycle sequencing reactions or unincorporated labeled nucleotides (dye- or radiolabeled dNTPs or ddNTPs) from DNA labeling reactions (e.g., PCR probe labeling, Nick Translation, or DNA end-labeling). The purified DNA qualifies for downstream applications like FISH (fluorescence in situ hybridization) or Southern/Northern blotting. The removal of unincorporated labeled nucleotides is a precondition for determining the DNA labeling rate. For good recovery rates, labeled DNA fragments must be at least 20 bp in length.

### Sephadex® G-50 Resin

Sephadex® G-50 is a well-established gel filtration resin for desalting and buffer exchange of biomolecules:

- Quickly desalts, removes contaminants, and transfers to a new buffer in a single step
- Suitable for purification of DNA from small molecules by gel filtration
- Small bead size for shorter diffusion distances

Sephadex® G-50 is a hydrophilic matrix that swells in buffer and acquires its chromatographic properties under wet conditions. In this state, beads have a pore size of 700 Da that allow, e.g., hydrated salt ions or dNTPs (dye-labeled or radiolabeled) to enter into the pores while DNA/RNA > 20 bases and most other biomolecules stay outside.

### Sephadex® G-50 Columns

Order No	Description	Amount
SC0500	Mobicol Spin G-50	20 columns
SC0510	Mobicol Spin G-50	100 columns

*Equilibrated in Mobicol Spin Buffer (10 mM Tris/HCl pH 7.6; 1 mM EDTA)*

*Sephadex® is a registered trademark of Cytiva.*

## Mobicol Spin Desalting Columns – Separation of Small Contaminants

### Applications

- ▶ Purification of nucleic acids from small contaminants
- ▶ Desalting and buffer exchange especially of small-sized nucleic fragments and of dNTPs



### Description

Mobicol Spin Desalting columns are perfectly suited for a fast and efficient routine purification of nucleic acids. The well-trying silica resin of the Mobicol Spin Desalting columns is comprised of uniform beads of a porous, form-stable material. It allows the selective removal of molecules as small as NaCl. Mobicol Spin Desalting columns offer numerous advantages:

- Improved recovery of small-sized nucleic acid fragments
- Desalting efficiency > 90%
- Silica resins are allowed to fall dry without any loss in purification efficiency
- Easy handling and long-term storage stability
- No sample dilution
- Reproducible results with simplified protocols
- Sample loading and recovery in less than 3 minutes
- Spin column chromatography combines effectiveness of gel filtration with speed of centrifugation

Nucleic acids purified with these columns are well suited for use in many downstream applications,

including sequencing, labeling, PCR, cloning, and in vitro transcription. The separation performance of the Mobicol Spin Desalting columns is vigorously tested with a wide range of nucleic acid fragments. In particular, fragments of 100 bp or smaller showed an improved recovery rate.

The desalting matrix separates molecules by differences in size (size exclusion chromatography). Nucleic acid molecules larger than the pore size are excluded from the resin. These molecules quickly move through the matrix bed when the column is centrifuged. Molecules smaller than the pore size, like hydrated salt ions, do enter the pores of the matrix beads and are held back. Thus, nucleic acid molecules are eluted from the column in order of decreasing molecular size.


Since this silica matrix material works without moisture expansion (contrasting Sepharacryl® and Sepharose® resins, which swell in buffer), resins are allowed to fall dry without any loss in purification efficiency. Therefore, Mobicol Spin Desalting columns give reproducible results of high quality and are much easier to handle. Furthermore, the columns are long-term storable in dry state.

### Silica-Based Columns

Order No	Description	Amount
 <b>SCO100</b>	Mobicol Spin Desalting	20 columns
 <b>SCO110</b>	Mobicol Spin Desalting	100 columns

*Including Mobicol Spin Buffer (10 mM Tris/HCl pH 7.6; 1 mM EDTA) for equilibration*

# Pre-filled Mobicols



Structure of IDA in complex with Ni<sup>2+</sup>

## Ni-IDA Columns for Purification of His-tagged Proteins

### Applications

- ▶ Purification of recombinant polyhistidine-tagged proteins
- ▶ Gravity-flow or spinning columns

### First Choice for Higher Yields of His-Tagged Proteins

- Easy, fast, and cost-effective routine purification of recombinant polyhistidine-tagged proteins:
  - under native and denaturing conditions
  - starting from diverse expression systems, e.g., *E. coli*, yeast, insect, and mammalian cells
- Universal use: suitable for small proteins, large protein complexes, proteins with low expression rates
- High binding capacity and high affinity
- Improved target specificity by optimized silica-based Ni<sup>2+</sup>-IDA matrix
- Imidazol-free loading and washing buffer

### His-Tagged Protein Purification

Histidine-tagged protein purification uses affinity chromatography to capture recombinant proteins with 4 - 10 histidine residues. Histidine-tagged


proteins can be purified from a wide range of expression systems, including bacterial, yeast, plant, and mammalian cells under native or denaturing conditions.

Ni-IDA columns from MoBiTec contain a silica-based resin. The form-stable silica matrix is pre-charged with Ni<sup>2+</sup> ions and allows purification on the principle of Immobilized Metal Ion Affinity Chromatography (IMAC). Binding of proteins is based on the interaction between the polyhistidine tag of a recombinant protein and immobilized Ni<sup>2+</sup> ions. The chelating group of the Ni-IDA resin is based on iminodiacetic acid (IDA), which enables strong and efficient binding of target protein to the IMAC matrix.

In contrast to traditional IDA matrices, MoBiTec Ni-IDA matrix is an optimized matrix with low density of IDA ligands. This non-saturating surface concentration of IDA eliminates almost all non-specific interactions of contaminating host proteins with the adsorbent. As a result, Lab Columns Ni-IDA and Mobicol Spin Ni-IDA columns provide much higher target protein purity.

IDA is a tridentate chelator which occupies three of the six binding sites in the coordination sphere of the Ni<sup>2+</sup> ion. The remaining three coordination sites are usually occupied by water molecules and can be exchanged with histidine residues of the recombinant protein.

# Pre-filled Mobicols



**Procedure**

The purification of His-tagged proteins consists of 4 steps: cell lysis, loading & binding, washing, and elution. Cleared cell lysates are loaded onto the matrices. His-tagged proteins are bound while other proteins pass through the matrix. After washing, His-tagged proteins are eluted with buffer from the matrix.

## Mobicol Spin Ni-IDA Columns

- ▶ For protein purification by spinning
- ▶ Compatible with standard microcentrifuges
- ▶ The purification run takes approximately 10 min
- ▶ Very high binding capacity: up to 12 mg protein per spin column
- ▶ Excellent protein recovery rate of > 90%
- ▶ No commercial license required

Mobicol Spin Ni-IDA columns are pre-packed, single-use spin columns for simple, small-scale purification of histidine-tagged proteins and rapid expression screening. This technology allows one-step purification of almost any His-tagged protein from any expression system under native or denaturing conditions.



## Lab Columns Ni-IDA

- ▶ For protein purification by gravity-flow
- ▶ Maximal binding capacity: up to 90 mg protein per column
- ▶ Protein recovery rate of > 80%
- ▶ Columns are long-term storable and reusable

Lab Columns Ni-IDA are designed for fast and simple purification of polyhistidine-tagged proteins without the need for a pump or purification system. Special column frits protect the chromatography medium from running dry during the process. A single column allows purification of approximately 90 mg of protein. The purified protein can be eluted in a small volume, resulting in a highly concentrated target protein.




Order No	Description	Amount
▶ PR-HTK105	Mobicol Spin Ni-IDA	5 columns
▶ PR-HTK110	Mobicol Spin Ni-IDA	10 columns
▶ PR-HTK004	Lab Columns Ni-IDA	4 columns
▶ PR-HTK010	Lab Columns Ni-IDA	10 columns
▶ PEG01*	E. coli pEG-His1 Vector	5 µg DNA

Columns contain dry matrix; storable at RT for at least 1 year.

\* Convenient cloning and subsequent purification of your protein of interest: pEG-His1 vector for IPTG-inducible expression of gene products with C-terminal RGS motif and His<sub>6</sub>-tag. 5' and 3' primers for sequencing are provided in addition.

# Pre-filled Mobicols



## Mobicols with Immobilized Enzymes

### Applications

- ▶ Dephosphorylation of vector DNA to prevent self-ligation
- ▶ Removal of DNases and RNases
- ▶ Protein digestions for peptide mapping and protein sequencing
- ▶ Diverse proteomic applications requiring maximum activity and selectivity
- ▶ Reproducible protein digestion and mass spectrometry-based protein identification
- ▶ Generation of antibody fragments

Enzyme immobilization provides enzyme reutilization and prevents the loss of enzyme activity. The technique further results in increased activity by providing a more suitable microenvironment for the enzyme/substrate reaction. Moreover, immobilized systems can provide better enzyme thermostability and pH tolerance. The method curtails redundant downstream and purification processes.

MoBiTec offers mini columns containing a G3m matrix with covalently immobilized enzymes for quick substrate-to-product turnover. The immobilized enzymes are stable in aqueous media at a pH range of 5 to 10 and column bleeding is negligible. The “stiff” linkers, which keep the enzymes off the matrix

surface, effectively eliminate steric hindrance. This results in high activity of enzymes in the immobilized state.

The enzyme reaction occurs as soon as the substrate is loaded onto the column. The sample is recovered from the column quantitatively either by gravity-flow or by centrifugation. Since the enzymes are covalently bound, they remain in the column after reaction and product elution, making the columns reusable. The eluted product solution is free of enzyme and does not require enzyme removal steps. This saves time and valuable resources.

### Features


- ▶ Covalently bound enzymes enable fast and easy handling
- ▶ Immobilization prevents loss of enzyme activity
- ▶ High enzymatic densities and very high activities
- ▶ Simple product recovery by centrifugation or elution by gravity-flow
- ▶ Suitable for small and large volumes – no sample dilution
- ▶ Treat larger volumes using syringes or chromatography tubing with standard Luer-Lock connections

# Pre-filled Mobicols

- ▶ The eluted product solution is free of enzyme and does not require enzyme removal steps
- ▶ Enjoy short reaction and treatment times
- ▶ Reusability – use the columns multiple times
- ▶ Sample solutions down to 20 µl can be treated and recovered

Order No	Description	Amount
 <b>P3102</b>	Mobicol with immobilized Endoproteinase Glu-C (Prot.V8)	22 U
 <b>P3402</b>	Mobicol with immobilized Papain	0.6 U
 <b>P3122</b>	Mobicol with immobilized Pepsin	0.4 mAnson-U
 <b>P3502</b>	Mobicol with immobilized Proteinase K	0.7 mAnson-U
 <b>P3302</b>	Mobicol with immobilized TLCK-α-Chymotrypsin	1.4 U
 <b>P3702</b>	Mobicol with immobilized TPCK-Trypsin	260 S&T-U
 <b>A3202</b>	Mobicol with immobilized Alkaline Phosphatase (CIP)	100 U
 <b>A3102</b>	Mobicol with immobilized b-Galactosidase	15 U

All columns contain 200 µl immobilized enzyme matrix. Product includes concentrated buffer, a detailed handbook, a data sheet with the description of the immobilized enzyme, the specific buffer composition, and a protocol for the use of the columns.



# Lab Columns – Empty Chromatography Columns

## Empty Chromatography Columns for Lab Scale with Luer-Lock Connections

### Applications

- ▶ Affinity chromatography
- ▶ Size exclusion chromatography
- ▶ Ion exchange chromatography

The empty syringe-based columns are available in three sizes (2.5 ml, 5 ml, and 10 ml) and come with an inserted lower filter with the pore size of your choice (10 µm, 35 µm, or 90 µm). Each column has Luer-Lock connectors at both ends for tubing or syringes. To close the columns, the 5 ml and 10 ml columns include an additional cap whereas the 2.5 ml columns have a Luer-Lock closing cap. The 2.5 ml columns can be centrifuged and autoclaved (with inserted filters only at 110 °C). Each package of columns includes two One-Way Stopcocks.

Upper filters (used on top of the matrix) and Luer-Lock adapters are available as well.

### Shared Features 2.5 ml, 5 ml, and 10 ml Columns

- ▶ Empty gravity-flow columns
- ▶ Compatible with most chromatography resins
- ▶ Luer-Lock connectors at both ends for connecting, e.g., syringes or stopcocks
- ▶ Can be closed with additional closing cap
- ▶ Smallest and largest volumes can be treated easily
- ▶ Come with lower filters with 10 µm, 35 µm, or 90 µm pore size
- ▶ Upper filter (used on top of the matrix) available separately



- ▶ 2 One-Way Stopcocks in each package of 20 columns
- ▶ Luer-Lock adapters can be ordered separately

### Additional Features 2.5 ml Columns:

- ▶ Can be closed tightly with Luer-Lock closing cap
- ▶ Tight up to 4 atmospheres pressure
- ▶ Can also be centrifuged and autoclaved (with inserted filters only at 110 °C)
- ▶ FPLC compatible

# Lab Columns – Empty Chromatography Columns

Order No	Description	Amount
 <b>S10121</b>	2.5 ml Lab Columns - Empty Chromatography Columns (10 µm filter pore size)	20
 <b>S1012</b>	2.5 ml Lab Columns - Empty Chromatography Columns (35 µm filter pore size)	20
 <b>S10129</b>	2.5 ml Lab Columns - Empty Chromatography Columns (90 µm filter pore size)	20
 <b>S10131</b>	5.0 ml Lab Columns - Empty Chromatography Columns (10 µm filter pore size)	20
 <b>S1013</b>	5.0 ml Lab Columns - Empty Chromatography Columns (35 µm filter pore size)	20
 <b>S10139</b>	5.0 ml Lab Columns - Empty Chromatography Columns (90 µm filter pore size)	20
 <b>S10141</b>	10 ml Lab Columns - Empty Chromatography Columns (10 µm filter pore size)	20
 <b>S1014</b>	10 ml Lab Columns - Empty Chromatography Columns (35 µm filter pore size)	20
 <b>S10149</b>	10 ml Lab Columns - Empty Chromatography Columns (90 µm filter pore size)	20
 <b>S1210</b>	Upper filter 10 µm for 2.5 ml Lab Columns	20
 <b>S1235</b>	Upper filter 35 µm for 2.5 ml Lab Columns – (thickness 1.5 mm)	20
 <b>S123532</b>	Upper filter 35 µm for 2.5 ml Lab Columns – (thickness 3.2 mm)	20
 <b>S1290</b>	Upper filter 90 µm for 2.5 ml Lab Columns	20
 <b>S1310</b>	Upper filter 10 µm for 5.0 ml Lab Columns	20
 <b>S1335</b>	Upper filter 35 µm for 5.0 ml Lab Columns – (thickness 1.5 mm)	20
 <b>S133532</b>	Upper filter 35 µm for 5.0 ml Lab Columns – (thickness 3.2 mm)	20
 <b>S1390</b>	Upper filter 90 µm for 5.0 ml Lab Columns	20
 <b>S1410</b>	Upper filter 10 µm for 10 ml Lab Columns	20
 <b>S143532</b>	Upper filter 35 µm for 10 ml Lab Columns – (thickness 3.2 mm)	20
 <b>S1490</b>	Upper filter 90 µm for 10 ml Lab Columns	20

*\*The appropriate filter pore size depends on your matrix. The matrix used should not pass through the filter pores.*

# Lab Columns – Empty Chromatography Columns

## Column Sets

MoBiTec's column sets are a composition of mini columns and syringe-based lab columns. All columns are easy to handle and can be used in a wide variety of chromatography applications.

These sets of empty columns contain four different sizes with three units of each:

- ▶ Mobicol (Mini Column up to 800 µl) (3x)
- ▶ 2.5 ml Lab Columns (3x)
- ▶ 5 ml Lab Columns (3x)
- ▶ 10 ml Lab Columns (3x)

Each column contains a Luer-Lock and a closing cap, a plug, and an inserted lower filter with your choice of pore size (10 µm, 35 µm, or 90 µm). Additionally, two stopcocks are included. An optional upper filter used on top of the matrix can be ordered separately.

### Available Sets

Column Set 1: 10 µm pore size filters inserted

Column Set 2: 35 µm pore size filters inserted

Column Set 3: 90 µm pore size filters inserted



Order No	Description	Amount
 <b>S10011</b>	<b>Column Set 1:</b> Mobicol M1050; 2.5 ml; 5 ml and 10 ml Lab Columns (each with 10 µm pore size filter)	3 each
 <b>S1001</b>	<b>Column Set 2:</b> Mobicol M1050; 2.5 ml; 5 ml and 10 ml Lab Columns (each with 35 µm pore size filter)	3 each
 <b>S10019</b>	<b>Column Set 3:</b> Mobicol M1050; 2.5 ml; 5 ml and 10 ml Lab Columns (each with 90 µm pore size filter)	3 each
 <b>S10031</b>	Set upper filters 10 µm for S10011	1 set
 <b>S1003</b>	Set upper filters 35 µm for S1001	1 set
 <b>S10039</b>	Set upper filters 90 µm for S10019	1 set

**M1050:** Mobicol Gravi Fixed with Luer-Lock and screw cap, fixed outlet plug



## Custom Filling: Mini Columns Packed with the Matrix of Your Choice

### Benefits Contract Filling

- ▶ Save valuable resources – outsource labor- and cost-intensive steps
- ▶ Complete documentation
- ▶ On-time delivery
- ▶ Ready-to-use chromatography columns can be straightforwardly integrated into your kits
- ▶ Renowned for quality:
  - ▶ MoBiTec invented and designed the columns (patented)
  - ▶ Highest quality for more than 30 years, Made in Germany
  - ▶ Our mini columns have already been sold several million times worldwide. They are an integral part of many laboratory diagnostics kits for clinical routine & research

### Target Group

- ▶ Life Science Research
- ▶ Laboratory Diagnostics
- ▶ Food Safety

Pre-Filled Mobicols are an integral part of many kits provided by global leaders in the field of bio-analytical testing solutions. Furthermore, they are used as a complement in laboratory diagnostics kits

or as an asset in assays for Life Science research. Contact our team for working out a solution tailored to your needs.

### Cooperation Example:

#### **Label IT® Nucleic Acid Labeling Kits**



Mirus Bio (Madison, USA) has developed a unique technology to label DNA and RNA. The *Label IT®* Technology is a non-enzymatic labeling method that facilitates the direct covalent attachment of either fluorescent or non-fluorescent labels to any nucleic acid using a simple one-step chemical reaction. After the labeling reaction, the sample should be purified using MoBiTec's G-50 microspin purification columns (Mobicol Spin G-50,  **SCO500** or  **SCO510**). Following microspin column purification, the researcher will get a 100% recovery of the labeled nucleic acid.



# Innovative Tools for Molecular and Cell Biology

## MoBiTec GmbH

MoBiTec develops cutting-edge research tools for the global life science, diagnostic, and pharmaceutical community. These include versatile purification columns and unique vector systems for protein expression in different species. Moreover, MoBiTec is specialized in labeling and detection technologies. Our products are distributed worldwide.

Next to our own portfolio, MoBiTec is a leading distributor for innovative life science reagents and services. We represent a number of well-established international companies:



For the most up-to-date listing of our supplier companies please visit our website

[www.mobitec.com](http://www.mobitec.com)



MoBiTec GmbH

Lotzestrasse 22a | 37083 Goettingen | Phone: +49 (0)551-70722-0 | Fax: +49 (0)551-70722-22 | E-mail: [info@mobitec.com](mailto:info@mobitec.com)