

VirusGEN[®]

Transfection Reagent and Enhancers

Viral Vector Production Platform
for Gene and Cell Therapy Manufacturing



Research Grade | SELECT[™] | GMP



OUR MISSION

Deliver Innovation.

Empower Discovery.

Improve Life.



THE TRANSFECTION EXPERTS


In 1995, Mirus began its journey with the goal of providing the best tools for gene delivery.

With over 25 years of experience and multiple breakthroughs in the field of transfection, Mirus is a global leader and trusted brand in nucleic acid delivery.

Driven by our passion for science, we take pride in the fact that our team of chemists and biologists push the boundaries of transfection technology to enable groundbreaking discoveries and the development of life-changing biotherapeutics throughout the world.

We care about **transfection** as much as you care about your research.

GMPTRANSFECTION.COM



Contents

Product Line Overview..... 1

Performance 5

Quality..... 15

Cost Savings 23

Product Information..... 25





VirusGEN® VIRAL VECTOR PRODUCTION PLATFORM

OVERVIEW

Background: Anatoly P, PhD, Chemist
Foreground: Nick R, PhD, Chemist
Keith P, PhD, Chemist

The Mirus Difference

VirusGEN® Biomimetic Transfection Technology

Many traditional non-viral methods for nucleic acid delivery were characterized by single-component, synthetic formulations, and the majority of first-generation transfection reagents relied on simple, linear polymers such as polyethyleneimine (PEI) or liposomes that transported DNA to cells via non-specific, electrostatic interactions.

In contrast, the technology underlying *TransIT*-VirusGEN® and other Mirus reagents utilizes fully synthetic components to more closely mimic the natural mechanisms that viruses developed over billions of years to deliver nucleic acid cargo to cells and subsequently release them to their nuclear or cytoplasmic destinations. These biomimetic transfection formulations rely on lipid-polymer nanocomplexes (LNPCs) and other components that enable researchers to achieve higher titers than what can be accomplished with PEI-derived polymers or liposomes alone.

**FROM R&D
TO GMP**

The Transfection Experts

With a history spanning more than 25 years of innovation in gene delivery technology, Mirus Bio is a long-standing leader in the field of transfection. Our VirusGEN® product line is a breakthrough in the generation of high-titer viral

vectors for gene and cell therapy and enables biopharmaceutical manufacturers to produce viral vectors at scales that were not previously possible in order to bring more doses to patients in need.

VirusGEN® Transfection Reagents and Enhancers for Production of AAV and Lentivirus

TransIT-VirusGEN® GMP Transfection Reagent is designed to enhance delivery of packaging and transfer vector DNA to suspension and adherent HEK 293 cell types in order to increase production of recombinant lentivirus and adeno-associated virus (AAV). Key benefits of *TransIT-VirusGEN®* GMP Transfection Reagent include:

- ✓ **Performance**—Efficient DNA delivery for large-scale production of high-titer viral vectors
- ✓ **Quality**—Available from research-grade to GMP
- ✓ **Flexibility**—Compatible with different virus production platforms and repeat filtration
- ✓ **Animal Origin Free**—Fully synthetic transfection reagent formulation

With industry-leading performance in high-titer virus manufacturing, *TransIT-VirusGEN®* GMP offers a simplified, cost-effective workflow, making it the superior choice for large-scale therapeutic virus production. Titers can be further increased with optimized enhancers included in the VirusGEN® GMP AAV Transfection Kit and VirusGEN® GMP LV Transfection Kit.



VirusGEN® GMP AAV Transfection Kit



VirusGEN® GMP LV Transfection Kit

TransIT-VirusGEN® Transfection Reagents



TransIT-VirusGEN®

- Research-grade reagent
- For R&D and Discovery



TransIT-VirusGEN® SELECT

- Research-grade reagent
- Additional quality testing
- For preclinical applications



TransIT-VirusGEN® GMP

- GMP-grade reagent
- For clinical trials and commercial manufacturing

VirusGEN® Complex Formation Solution and Enhancers



TransIT-VirusGEN® AAV Complex Formation Solution & Enhancer

- Available as a research-grade or GMP solution
- Optimized for complex formation with TransIT-VirusGEN® transfection reagents
- Enhances AAV titers when used with TransIT-VirusGEN® transfection reagents



TransIT-VirusGEN® LV Enhancer

- Available as a research-grade or GMP solution
- Enhances lentiviral titers when used with TransIT-VirusGEN® transfection reagents



TransIT-VirusGEN® LV Complex Formation Solution

- Available as a research-grade or GMP solution
- Optimized for complex formation with TransIT-VirusGEN® transfection reagents

From Research to Commercialization

TransIT-VirusGEN® Transfection Reagent provides industry-leading virus titers for AAV and lentivirus production and is available as research-grade, SELECT and GMP formats to support seamless transitions from discovery through clinical trials and commercialization.

VirusGEN® AAV and LV complex formation solutions and enhancers are available as research-grade or GMP-grade products.



VirusGEN® VIRAL VECTOR PRODUCTION PLATFORM

PERFORMANCE

Background: Tony L, Biologist
Foreground: Myla G, Production Scientist

VirusGEN® GMP for AAV Production

Adeno-associated virus (AAV) vectors are the most rapidly growing gene therapy delivery platform for the treatment of a vast array of human diseases. The unique properties of AAV (non-pathogenic, non-integrating genome, low immunogenicity and broad tropism) combined with recent advances in capsid optimization have catalyzed rapid growth in the gene therapy field.

The VirusGEN® GMP AAV Transfection Kit was specifically developed for high-titer, industrial scale AAV manufacturing in suspension HEK 293 cells and includes *TransIT-VirusGEN®* GMP Transfection Reagent (150 ml) and VirusGEN® GMP AAV Complex Formation Solution and Enhancer (5 × 1 L).

VirusGEN® GMP AAV Transfection Kit



Kit includes: *TransIT-VirusGEN®* GMP Transfection Reagent (150 ml), VirusGEN® GMP AAV Complex Formation Solution and Enhancer (5 × 1 L)

AAV Performance

TransIT-VirusGEN® Transfection Reagent outperforms commonly used competitor reagents when measured by functional AAV titers or AAV genome copies in suspension 293 cells (**FIGURE 1**). This is also reflected in the higher GC/capsid ratio obtained with *TransIT-VirusGEN®*. The amount of AAV obtained is further increased in suspension HEK 293 cells with the use of VirusGEN® AAV Complex Formation Solution and Enhancer, included in the VirusGEN® AAV Transfection Kit.

The VirusGEN® AAV Transfection Kit demonstrates superior performance across multiple suspension HEK 293 cell lines and media formulations (**FIGURE 2**), and AAV titers obtained with *TransIT-VirusGEN®* Transfection Reagent are significantly higher than what can be obtained with competitor reagents in adherent HEK 293 cells (**FIGURE 3**). Note that VirusGEN® AAV Complex Formation Solution and Enhancer is only recommended for suspension HEK 293 cells.




FIGURE 1. *TransIT-VirusGEN[®]* and the *VirusGEN[®]* AAV Transfection Kit Produce High AAV Titers. (A) Expi293F[™] cells (Thermo Fisher Scientific) grown in Expi293[™] Expression Medium (Thermo Fisher Scientific) were used to generate recombinant AAV2 via transient transfection using *TransIT-VirusGEN[®]* Transfection Reagent (1.5:1 reagent-to-DNA ratio (vol:wt), Mirus Bio LLC), *VirusGEN[®]* AAV Transfection Kit (Mirus Bio LLC), Competitor A1 (1:1) or Competitor A3 (1:1) using the manufacturer's recommended protocol. AAV2 was produced by transfecting pAAV-hrGFP, pAAV-RC and pAAV-Helper plasmids (1:1:1 DNA ratio, 2 µg/ml-2 mg per flask, Agilent Technologies). Cells were transfected at a cell density of 2×10^6 cells/ml, 1 L in a 2.8 L Thomson shake flask. AAV2 was harvested at 72 hours post-transfection using chemical lysis. (A) Functional titers were determined via transduction of HT-1080 cells and GFP expression was measured 48 hours post-transduction using Guava[®] easyCyte[™] 5HT Flow Cytometer. AAV functional titers were measured from virus dilutions with less than 20% GFP positive cells. Genome copies were determined by ddPCR using primers and a probe targeting the CMV promoter. (B) Total assembled capsids were determined using the AAV2 Titration ELISA Kit (Progen) and the GC/capsid ratio was determined by dividing the number of genome copies by total assembled capsids for each condition. The error bars represent the range of duplicate flasks.




FIGURE 2. *TransIT-VirusGEN®* and the *VirusGEN®* AAV Transfection Kit are Compatible with Different Cell Types and Media Formulations. (A) Expi293F™ cells (Thermo Fisher Scientific) or (B) Viral Production cells (Thermo Fisher Scientific) were grown or adapted to Expi293™ Expression Medium (Thermo Fisher Scientific), LV-MAX™ Production (Thermo Fisher Scientific), FreeStyle™ F17 (Thermo Fisher Scientific) or BalanCD HEK293 (Irvine Scientific) media prior to transfection. AAV production using *TransIT-VirusGEN®* Transfection reagent (2:1 reagent-to-DNA ratio (vol:wt), Mirus Bio LLC), *VirusGEN®* AAV Transfection Kit (Mirus Bio LLC) was compared to Competitor A3 (1:1 reagent-to-DNA ratio (vol:wt)) and Competitor C1 (2:1 reagent-to-DNA ratio (vol:wt)) using the manufacturer's recommended protocol. AAV2 was produced by transfecting pAAV-hrGFP, pAAV-RC and pAAV-Helper plasmids (1:1:1 DNA ratio, 1.5 µg/ml = 3 µg/well, Agilent Technologies). Cells were transfected at a cell density of 2 million cells/ml. Harvested virus was used to transduce HT-1080 cells and GFP expression was measured 48 hours post-transduction using Guava® easyCyte™ 5HT Flow Cytometer. AAV functional titers were measured from virus dilutions with less than 20% GFP positive cells. The error bars represent the range of duplicate wells.




FIGURE 3. *TransIT-VirusGEN*® Outperforms PEI-based Reagents for AAV Production in Adherent 293 Cells. AAV2 was produced using adherent 293-AAV cells (Agilent Technologies) grown in DMEM + 10% FBS and transfected using pAAV-hrGFP, pAAV-RC and pAAV-Helper plasmids (1:1:1 DNA ratio, 1.5 µg/ml plasmid DNA, 3 µg/well, Agilent Technologies) using *TransIT-VirusGEN*® Transfection Reagent (2:1, vol:wt, Mirus Bio LLC), Competitor A3 (2:1) or Competitor C1 (2:1) or in 6-well plates. Harvested virus was used to transduce HT-1080 cells and GFP expression was measured 48 hours post-transduction using a Guava® easyCyte™ 5HT Flow Cytometer. Functional AAV titers were measured from virus dilutions with less than 20% GFP positive cells. The error bars represent the standard deviation of triplicate wells.

VirusGEN® GMP for Lentivirus Production

Lentiviral (LV) vectors provide a highly efficient means for modifying mammalian cells and have been extensively utilized for clinical gene and cell therapy applications. In particular, the success of chimeric antigen receptor (CAR) T-cell therapy has been made possible through lentiviral gene transfer.

The VirusGEN® GMP LV Transfection Kit was specifically developed for high-titer, industrial-scale LV manufacturing in suspension or adherent HEK 293 cells and includes *TransIT*-VirusGEN® GMP Transfection Reagent (150 ml), VirusGEN® GMP LV Complex Formation Solution (5 x 1 L) and VirusGEN® GMP LV Enhancer (5 x 1 L).

The VirusGEN® GMP LV Transfection Kit



Kit includes: *TransIT*-VirusGEN® SELECT Transfection Reagent (150 ml), VirusGEN® GMP LV Complex Formation Solution (5 x 1 L), VirusGEN® GMP LV Enhancer (5 x 1 L)

LV Performance

TransIT-VirusGEN® Transfection Reagent outperforms commonly used competitor reagents when measured by functional LV titers obtained from suspension HEK 293 cells (**FIGURE 4**). Lentivirus titers are further increased with the use of VirusGEN® LV Complex Formation Solution and VirusGEN® LV Enhancer, included in the VirusGEN® LV Transfection Kit (red bar).

The VirusGEN® LV Transfection Kit demonstrates high and consistent performance across multiple suspension HEK 293 cell lines and media formulations (**FIGURE 5**) and is compatible with a variety of commercially available packaging and transfer plasmid systems (**FIGURE 6**). The VirusGEN® LV Transfection Kit also achieves high functional LV titers in adherent HEK 293 cells, and the amount of total DNA can be titrated for optimal results, depending on experimental conditions (**FIGURE 7**).




FIGURE 4. TransIT-VirusGEN[®] Transfection Reagent and VirusGEN[®] LV Enhancer Outperform PEI-based Transfection Reagents. Viral Production cells (Thermo Fisher Scientific) were grown in LV-MAX[™] Production Medium (Thermo Fisher Scientific) and seeded at 4×10^6 cells/ml prior to transfection (25 ml/125 ml shake flask). Lentivirus was produced using the TransIT-VirusGEN[®] Transfection Reagent +/- VirusGEN[®] LV Enhancer (Mirus Bio, 3:1 reagent-to-DNA ratio (wt:vol)), Competitor A3 (2:1 reagent-to-DNA ratio (wt:vol)) or Competitor C1 (4:1 reagent-to-DNA ratio (wt:vol)) using the manufacturer's recommended protocol. Lentivirus was produced by transfecting 3rd generation vectors pALD-LentiEGFP-A transfer vector and pALD-VSV-G-A, pALD-Rev-A, pALD-GagPol-A packaging vectors at a 3:0.5:0.5:2 DNA ratio (Aldevron, 25 µg/flask (1 µg/ml final concentration)). The VirusGEN[®] LV Enhancer was added 21 hours post-transfection. (A) Functional titers from virus containing supernatant was measured following transduction of 293T/17 cells and GFP expression was measured at 72 hours post-transduction using Guava[®] easyCyte[™] 5HT Flow Cytometer. Lentivirus functional titers were measured from virus dilutions with less than 20% GFP positive cells. (B) p24 cCapsid protein was measured from virus containing supernatants using a HIV-1 p24 Antigen ELISA kit according to the manufacturer's protocol (Zeptometrix). The error bars represent the standard deviation of triplicate wells.




FIGURE 5. TransIT-VirusGEN[®] Transfection Reagent and VirusGEN[®] LV Enhancer Allow for Cell Line and Media Formulation Flexibility. Viral Production cells (293-VP) grown in LV-MAX[™] Production Medium (Thermo Fisher Scientific), Expi293[™] cells grown in Expi293[™] Expression Medium (Thermo Fisher Scientific) or FreeStyle[™] 293-F cells grown in FreeStyle[™] F17 (Thermo Fisher Scientific) were seeded at 4×10^6 cells/ml prior to transfection (2 ml/non-treated 6-well plate). Lentivirus was produced using the TransIT-VirusGEN[®] Transfection Reagent ± VirusGEN[®] LV Enhancer (Mirus Bio, 3:1 reagent-to-DNA ratio (wt:vol)) according to the recommended protocol. Lentivirus was produced by transfecting 3rd generation vectors pALD-LentiEGFP-A transfer vector and pALD-VSV-G-A, pALD-Rev-A, pALD-GagPol-A packaging vectors at a 3:0.5:0.5:2 DNA ratio (Aldevron, 25 µg/flask (1 µg/ml final concentration)). The VirusGEN[®] LV Enhancer was added 18 hours post-transfection. Functional titers from virus containing supernatant was measured following transduction of 293T/17 cells and GFP expression was measured at 72 hours post-transduction using Guava[®] easyCyte[™] 5HT Flow Cytometer. Lentivirus functional titers were measured from virus dilutions with less than 20% GFP positive cells. The error bars represent the standard deviation of triplicate wells.




FIGURE 6. Lentivirus Commercially Available Vector Competitor Comparison. Lentivirus production using *TransIT-VirusGEN®* Transfection Reagent ± the *VirusGEN®* LV Enhancer was compared to Competitor A3 (2:1 reagent-to-DNA ratio) and Competitor C1 (4:1 reagent-to-DNA ratio) using the manufacturer's recommended protocol. Lentivirus was produced by transfecting Expi293™ cells grown in Expi293™ Expression Medium (Thermo Fisher Scientific) at 4×10^6 cells/ml with two different LV vector mixtures (1 µg/ml, 2 µg/well) including: (A) 3rd generation vectors pALD-LentiEGFP-A transfer vector and pALD-VSV-G-A, pALD-Rev-A, pALD-GagPol-A packaging vectors (3:0.5:0.5:2 DNA ratio, Aldevron) or (B) 3rd generation vectors pLKO.1-puro-CMV-TurboGFP™ transfer vector (Sigma) and ViraSafe™ Pantropic Packaging mix (pRSV-Rev, pCMV-VSV-G, pCgpV, CellBio Labs) at a 3:0.5:0.5:2 DNA ratio. The *VirusGEN®* LV Enhancer was added at 18 hours post-transfection. Virus containing supernatant was used to transduce 293T/17 cells and GFP expression was measured at 72 hours post-transduction using Guava® easyCyte™ 5HT Flow Cytometer. Lentivirus functional titers were measured from virus dilutions with less than 20% GFP positive cells. The error bars represent the range of duplicate wells.




FIGURE 7. Adherent Lentivirus Production With Commercially Available Vector at Different Reagent-to-DNA Ratios. Lentivirus production using *TransIT-VirusGEN*® Transfection Reagent ± the *VirusGEN*® LV Enhancer using the recommended protocol. Lentivirus was produced by transfecting adherent 293T/17 cells grown in serum containing medium with two different LV vector mixtures (1 µg/ml, 0.5 µg/well of 24-well plate) including: (A) 3rd generation vectors pALD-LentiEGFP-A transfer vector and pALD-VSV-G-A, pALD-Rev-A, pALD-GagPol-A packaging vectors (3:0.5:0.5:2 DNA ratio, Aldevron) or (B) 3rd generation vectors pLKO.1-puro-CMV-TurboGFP™ transfer vector (Sigma) and ViraSafe™ Pantropic Packaging mix (pRSV-Rev, pCMV-VSV-G, pCgpV, CellBio Labs) at a 3:0.5:0.5:2 DNA ratio. The *VirusGEN*® LV Enhancer was added at 20 hours post-transfection. Virus containing supernatant was used to transduce 293T/17 cells and GFP expression was measured at 72 hours post-transduction using Guava® easyCyte™ 5HT Flow Cytometer. Lentivirus functional titers were measured from virus dilutions with less than 20% GFP positive cells. The error bars represent the standard deviation of triplicate wells.



VirusGEN® VIRAL VECTOR PRODUCTION PLATFORM

QUALITY

Susan E, Quality Supervisor


Peace of Mind: Mirus GMP Transfection

With the increasing regulatory obligations facing AAV and LV manufacturers, there is a need to utilize more stringent quality processes at earlier stages in the development of gene and cell therapies.

Our commitment to quality is integral to our brand and your success. With our VirusGEN® GMP family of products for recombinant AAV and lentivirus production, Mirus is dedicated to ensuring that customers can rely on Mirus for transfection reagents of exceptional quality from early stages of research and development through clinical trials and commercial manufacturing.



TransIT-VirusGEN® GMP:
Highest Quality Grade TransIT-VirusGEN® Available



VirusGEN® GMP products are manufactured in compliance with 21 CFR 820 and are held to the highest standards of traceability. All products are produced according to a GMP manufacturing process and tested in accordance with documented procedures.

Available Analytical Assays

TransIT-VirusGEN® Transfection Reagents

- ✓ Residual Reagent Assay
- ✓ Identity Assay

VirusGEN® AAV Complex Formation Solution and Enhancer

- ✓ Residual Reagent Assay
- ✓ Identity Assay

VirusGEN® LV Enhancer

- ✓ Identity Assay

QUALITY

Quality Control Documentation

Below is an example of the Certificate of Analysis for *TransIT-VirusGEN*® GMP Transfection Reagent that includes testing for: appearance, endotoxin, identity, mycoplasma and sterility.

CERTIFICATE OF ANALYSIS



Product Name: *TransIT-VirusGEN*® GMP Transfection Reagent, 150 mL
Product Number: MIR 6845-GMP
Lot Number: XXXX
Date of Manufacture: DD/MMM/YYYY
Retest Date: DD/MMM/YYYY
Storage Condition: -10 to -30 °C
Certificate Issue Date: DD/MMM/YYYY

Quality Control Testing and Results:

Test Method	Specification	Result	Pass/Fail
Appearance	Clear, Colorless Solution, Free of Particulates		
Endotoxin; USP <85>	≤ 1 EU/mL		
Identity (NMR) – Proprietary Active Ingredient 1	Conforms		
Identity (NMR) – Proprietary Active Ingredient 2	Conforms		
Identity – Water by GC	Detected		
Identity – Ethanol by GC	Detected		
Mycoplasma; USP <63>	None Detected		
Sterility; USP <71>	No Growth Observed		

The product was GMP manufactured in compliance with 21 CFR 820. This certificate confirms that the product meets the specifications above. The product was tested in accordance with documented procedures and approved as a result of meeting the required specification.

Amanda Rice, Quality Manager

This product is sold to the Buyer with a limited license to use this product for research use and further manufacturing, not for administration into humans. This product, or parts from this product, may not be re-packaged or re-sold without written permission from Mirus Bio LLC.

TransIT-VirusGEN® is a registered trademark of Mirus Bio LLC.

Mirus Bio LLC | 5602 Research Park Blvd, Ste 210 | Madison, WI 53719 USA | 608.441.2852 | FAX 608.441.2849 | www.mirusbio.com

Below is an example of the Certificate of Origin and BSE/TSE Statement for *TransIT-VirusGEN*® GMP Transfection Reagent documenting that it is chemically synthesized and is not manufactured with any animal-derived components.

Mirus

Product Name:	TransIT-VirusGEN® GMP Transfection Reagent, 150 mL
Product Number:	MIR 6845-GMP
Lot Number:	XXXX
Animal Derived Components:	None
Type of Manufacture:	Chemical Synthesis
Country of Manufacture:	USA
Certificate Issue Date:	DD/MMM/YYYY

This product is manufactured via chemical synthesis in the USA entirely from material of non-animal origin. The manufacture, packaging, storage, and transportation of these materials do not involve the use of material of animal origin. This information is to be used for the purpose of determining animal origin only and not to be confused with "country of origin" for import/export purposes.

This product has minimal risk of contamination with Bovine Spongiform Encephalopathy (BSE) or Transmissible Bovine Encephalopathy (TSE). Mirus Bio LLC does not have plans to change the production of this material in any way that would increase the risk of BSE/TSE contamination.




This product is sold to the Buyer with a limited license to use this product for research use and further manufacturing; not for administration into humans. This product, or parts from this product, may not be re-packaged or re-sold without written permission from Mirus Bio LLC.

TransIT-VirusGEN® is a registered trademark of Minus Bio LLC.



Mirus Bio LLC | 5602 Research Park Blvd. Ste 210 | Madison, WI 53719 USA | 608.441.2852 | FAX 608.441.2849 | www.mirusbio.com

QUALITY

Compare *TransIT-VirusGEN*[®] Reagents

Therapeutic Development Pipeline	Research & Development	Preclinical & Early Phase Clinical Trial	Late Phase Clinical Trial & Commercial Manufacturing
<i>TransIT-VirusGEN</i> [®] Reagents	 <i>TransIT-VirusGEN</i> [®]	 <i>TransIT-VirusGEN</i> [®] SELECT	 <i>TransIT-VirusGEN</i> [®] GMP
Composition	Identical <i>TransIT-VirusGEN</i> [®] Reagent Formulation (ready to use, chemically defined, animal origin free in 77.5% ethanol)		
AAV and LV Enhancer Kits Available	Available	Coming Soon	Available
Quality	R&D	Preclinical	GMP
Quality Control	Functional AAV titer assay	<ul style="list-style-type: none"> • Functional AAV titer assay • Formulation identity • Appearance • Sterility: per USP <71> • Bacterial endotoxin: per USP <85> • Mycoplasma: per USP <63> 	<ul style="list-style-type: none"> • Validated formulation identity assay • Appearance • Sterility: per USP <71> • Bacterial endotoxin: per USP <85> • Mycoplasma: per USP <63>
Configuration	0.3 ml; 0.75 ml; 1.5 ml; 5 and 10 pack, 1.5 ml vials	<ul style="list-style-type: none"> • 30 ml bottle • 150 ml bottle 	150 ml bottle
Manufacturing	Research grade manufacturing with aseptic filtration	Research grade manufacturing with aseptic filtration	GMP validated manufacturing process with sterile filtration
Raw Materials	Research grade		GMP grade critical raw materials
Available Documentation	<ul style="list-style-type: none"> • Certificate of Analysis • Certificate of Origin 	<ul style="list-style-type: none"> • Certificate of Analysis • Certificate of Origin Includes: TSE/BSE Statement 	<ul style="list-style-type: none"> • Certificate of Analysis • Certificate of Origin Includes: TSE/BSE Statement • DMF Available in 2022 with Quality Agreement
Available Analytical Assays	—		<ul style="list-style-type: none"> • Validated identity assay • Residual reagent assay

Compare VirusGEN® R&D Complex Formation Solution and Enhancers

Therapeutic Development Pipeline	Research & Development	Preclinical through Clinical Trials and Commercialization
VirusGEN® Enhancers	 <p>VirusGEN® Research Grade Complex Formation Solution and Enhancers</p> <ul style="list-style-type: none"> • VirusGEN® AAV Complex Formation Solution and Enhancer • VirusGEN® LV Complex Formation Solution • VirusGEN® LV Enhancer 	 <p>VirusGEN® GMP Complex Formation Solution and Enhancers</p> <ul style="list-style-type: none"> • VirusGEN® GMP AAV Complex Formation Solution and Enhancer • VirusGEN® GMP LV Complex Formation Solution • VirusGEN® GMP LV Enhancer
Composition	Proprietary aqueous solutions	
Quality Control	<ul style="list-style-type: none"> • Identity • Appearance • Sterility: per USP <71> • Bacterial endotoxin: per USP <85> • Osmolality • pH 	<ul style="list-style-type: none"> • Identity • Appearance • Sterility: per USP <71> • Bacterial endotoxin: per USP <85> • Osmolality • pH • Mycoplasma: per USP <63>
Configuration	100 ml	1 L
Manufacturing	Research grade with sterile filtration	GMP manufacturing process with sterile filtration
Raw Materials	USP-Grade Critical Raw Materials	
Available Documentation	Certificate of Analysis with Material Origin and TSE/BSE Statement	<ul style="list-style-type: none"> • Certificate of Analysis • Certificate of Origin Includes: TSE/BSE Statement • DMF Available in 2022 with Quality Agreement
Available Analytical Assays	—	<ul style="list-style-type: none"> • Residual Reagent Assay Available for VirusGEN® AAV Enhancer • Identity Assay

QUALITY



VirusGEN® VIRAL VECTOR PRODUCTION PLATFORM


WORKFLOW SAVINGS

Background: Austin S, Associate Scientist III
Foreground: Laura J, PhD, Vice President, Scientific Operations

Produce More Doses Per Batch


Gene and cell therapy manufacturers are facing increasing demands in two critical areas:

- 1. Identifying a cost-effective means to produce lentivirus and AAV vectors in sufficient quantities to address the rapidly expanding gene and cell therapy market.
- 2. Minimizing safety and regulatory risks by using GMP raw materials in upstream manufacturing processes.




TransIT-VirusGEN® GMP offers a solution to both of these challenges by providing industry-leading virus titers and compliance with GMP viral vector manufacturing. Higher titers directly translate into cost-savings opportunities by either scaling down biomanufacturing or producing more therapeutic doses per batch.


FROM R&D TO GMP **Product Information**


Research-Grade Reagents & Kits			
Product	Product No.	Description	
 TransIT-VirusGEN® Transfection Reagent	MIR 6703	0.3 ml	
	MIR 6704	0.75 ml	
	MIR 6700	1.5 ml	
	MIR 6705	5 x 1.5 ml	
	MIR 6706	10 x 1.5 ml	


 VirusGEN® AAV Transfection Kit	MIR 6750	Kit includes: <i>TransIT-VirusGEN®</i> Transfection Reagent (2 x 1.5 ml), <i>VirusGEN®</i> AAV Complex Formation Solution and Enhancer (100 ml)	
--	----------	---	--

 VirusGEN® LV Transfection Kit	MIR 6760	Kit includes: <i>TransIT-VirusGEN®</i> Transfection Reagent (2 x 1.5 ml), <i>VirusGEN®</i> LV Complex Formation Solution (100 ml), <i>VirusGEN®</i> LV Enhancer (100 ml)	
---	----------	--	--

SELECT-Grade Reagents & Kits			
Product	Product No.	Description	
 TransIT-VirusGEN® SELECT Transfection Reagent	MIR 6730	30 ml	
	MIR 6735	150 ml	

 TransIT-VirusGEN® SELECT AAV Transfection Kit	MIR 6770	Kit includes: <i>TransIT-VirusGEN®</i> SELECT Transfection Reagent (30 ml), <i>VirusGEN®</i> SELECT AAV Complex Formation Solution and Enhancer (1 L)	
	MIR 6775	Kit includes: <i>TransIT-VirusGEN®</i> SELECT Transfection Reagent (150 ml), <i>VirusGEN®</i> SELECT AAV Complex Formation Solution and Enhancer (1 L)	


 TransIT-VirusGEN® SELECT LV Transfection Kit	MIR 6780	Kit includes: <i>TransIT-VirusGEN®</i> SELECT Transfection Reagent (30 ml), <i>VirusGEN®</i> SELECT LV Complex Formation Solution (1 L), <i>VirusGEN®</i> SELECT LV Enhancer (1 L)	
	MIR 6785	Kit includes: <i>TransIT-VirusGEN®</i> SELECT Transfection Reagent (150 ml), <i>VirusGEN®</i> SELECT LV Complex Formation Solution (1 L), <i>VirusGEN®</i> SELECT LV Enhancer (1 L)	

GMP-Grade Reagents & Kits			
Product	Product No.	Description	
 TransIT-VirusGEN® GMP Transfection Reagent	MIR 6845-GMP	150 ml	

 VirusGEN® GMP AAV Transfection Kit	MIR 6815-GMP	Kit includes: <i>TransIT-VirusGEN®</i> GMP Transfection Reagent (150 ml), <i>VirusGEN®</i> GMP AAV Complex Formation Solution, <i>VirusGEN®</i> GMP AAV Enhancer (5 x 1 L)	
--	--------------	--	--

 VirusGEN® GMP LV Transfection Kit	MIR 6825-GMP	Kit includes: <i>TransIT-VirusGEN®</i> GMP Transfection Reagent (150 ml), <i>VirusGEN®</i> GMP LV Complex Formation Solution (5 x 1 L), <i>VirusGEN®</i> GMP LV Enhancer (5 x 1 L)	
---	--------------	--	--

SAVINGS



Will H, Production Scientist

Terms and Conditions

1. Governing Provisions: All Mirus Bio LLC products are sold and shipped subject to these Terms and Conditions, unless modified by a quotation issued by Mirus Bio LLC or an agreement signed by both Mirus Bio LLC and the customer. These Terms and Conditions supersede any other written document or oral discussion including terms or conditions appearing on a purchase order issued by the customer.

2. Product Safety: Mirus Bio LLC manufactures products for research use by qualified individuals. These products may contain chemicals that may be harmful if misused. Customers are responsible for reviewing and following the instructions included in the product protocols and SDS that accompany the product. Please contact techsupport@mirusbio.com with any questions regarding the safe use of Mirus Bio LLC products. We recommend our products be used in accordance with the NIH guidelines developed for genetic research.

3. Use and Regulatory Restrictions: Certain products produced by Mirus Bio LLC may be labeled for "Further manufacturing, not for administration into humans" and are not intended for diagnostic use or direct administration into humans or animals. The product(s) are not submitted for regulatory review by any governmental body or other organization. Mirus Bio does not validate the product(s) for clinical or diagnostic use, for safety and effectiveness, or for any other specific use or application. The end user is solely responsible for ensuring that use of the product(s) complies with applicable laws, regulations and governmental policies, and that the product(s) are used by technically-qualified individuals. The end user must obtain all necessary approvals and permissions needed. It is solely the end user's responsibility to make sure the product(s) are suitable for their particular use. Products may not be re-packaged or re-sold without written permission from Mirus Bio LLC. The buyer agrees not to infringe upon Mirus patents or to attempt to reverse engineer, reconstruct, synthesize or otherwise modify Mirus Bio products.

4. Warranties: Mirus Bio LLC products are warranted to meet or exceed the stated specifications for the usable life of the product defined on the Certificate of Analysis. This warranty is void if the customer has altered or misused the product or failed to store the product as recommended. The warranty does not include defects or failures associated with materials supplied to Mirus Bio LLC by the customer. Mirus Bio LLC disclaims any and all responsibility for any injury or damage which may be caused by the failure of the buyer or any other person to use these products in accordance with the conditions outlined herein. Mirus Bio's sole obligation and the customer's sole remedy are limited to free replacement of products that fail to meet our minimum performance specifications. THIS WARRANTY IS EXCLUSIVE. MIRUS BIO LLC

SPECIFICALLY DISCLAIMS ANY OTHER WARRANTIES OF ANY KIND OR NATURE, DIRECT OR INDIRECT, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, THE SUITABILITY, FITNESS OR MERCHANTABILITY FOR ANY PURPOSE. MIRUS BIO WILL NOT BE LIABLE FOR ANY CLAIMS OR DAMAGES UNDER ANY LEGAL THEORY INCLUDING, BUT NOT LIMITED TO CONTRACT, NEGLIGENCE, STRICT LIABILITY OR TORT IN CONNECTION WITH THE FAILURE OF MIRUS BIO LLC PRODUCTS TO PERFORM IN ACCORDANCE WITH THE STATED SPECIFICATIONS. Mirus Bio LLC is not liable for any loss, damages or penalties resulting from delays in manufacturing or delivery of products or services. This warranty applies to all products and services offered by Mirus Bio LLC including catalog products, bulk products, custom products and custom services. Unless otherwise agreed, all technical assistance and information is provided free of charge and we make no warranty regarding the accuracy or utility of such information or assistance.

5. Patents: Mirus Bio's products and processes are covered by one or more patents and are subject to other trade secret and proprietary rights. No transfer or grant of rights under any patent or patents is made or is to be implied by any provision under these Terms and Conditions other than a limited license to use the products as designated. You agree not to infringe upon such rights or to attempt to reverse engineer, reconstruct, synthesize or otherwise modify Mirus Bio products. Certain applications of Mirus Bio products may require licenses from other parties. Determining the existence and scope of such third party intellectual property is the responsibility of the customer. Purchase of the product provides the customer with a limited non-transferable license under any Mirus Bio patents or patent applications to use the product for internal research unless there is a written limitation to this license in the product literature. The customer is responsible for carefully reviewing the product literature and respecting any limitations to this license.

6. Customer agrees to indemnify and hold Mirus Bio LLC harmless from and against any and all claims, losses, damages, costs and expenses (including reasonable attorneys' fees and amounts paid in settlement in good faith) which may be suffered or incurred by Mirus Bio LLC (including, without limitation, as a result of any claim or action by a third party) arising in any manner from Customer's use of any Products sold from this Web site or otherwise arising out of any act or omission of Customer, its employees, agents or representatives, or by reason of breach of or failure to perform any of Customers obligations hereunder.

7. Orders, Shipment and Delivery: Orders typically ship from our facility in 1-2 business days. Orders placed on Friday or prior to a holiday will ship on

the following week. Orders for regularly stocked items may only be cancelled prior to shipment of the order. Shipping charges are prepaid and added to the customer's invoice. Orders are shipped by Federal Express. Unless agreed upon in writing, all shipments are Ex Works (EXW), the Mirus facility in Madison, WI, U.S.A., Incoterms 2020 or such other place as may be designated by Mirus. Costs for transportation and risk of loss transfer to the buyer upon delivery to the carrier. The customer is responsible for any tariffs, duties or taxes associated with international shipments.


8. Copyright and Trademarks: All content contained in printed materials or on the Mirus Bio LLC website is the copyright of Mirus Bio LLC unless specifically identified otherwise. Trademarks referenced herein are either registered trademarks of Mirus Bio LLC or third party trademarks, which are the property of their respective owners. All rights are reserved. Please contact license@mirusbio.com to request permission to use copyrighted material or trademarks of Mirus Bio LLC.

9. Payment Terms: All sales are net 30 days of invoice in U.S. dollars unless other terms are agreed to in writing by Mirus Bio LLC. Prices are subject to change without notice. Mirus Bio LLC accepts the following credit cards: VISA, MasterCard, American Express and Discover. Mail payments to:

Mirus Bio LLC
5602 Research Park Blvd, Ste 210
Madison, WI 53719 USA
Attention: Accounts Payable
Fax: 608.441.2849

10. Returns: Mirus Bio LLC must authorize all product returns. Customers wishing to return a product should contact us at techsupport@mirusbio.com. Product that we confirm damaged, incomplete or do not meet our minimum written specifications will be replaced per our warranty (see above). The customer must bear the cost of shipping the product back to Mirus Bio LLC under the conditions specified by Mirus Bio LLC. All such returns may be subject to a 25% restocking fee.

11. These Terms and Conditions constitute the entire agreement between Mirus Bio LLC and Customer with respect to Customer use of the Web site and Customer purchase of products hereunder, except as the foregoing (i) may be amended from time to time by Mirus Bio LLC, or (ii) as related to Customer purchase of products, may be superseded by any express conflicting terms or supplemented by any express additional terms in a separate written contract signed by authorized representatives of Mirus Bio LLC and Customer.



Mirus Bio LLC
5602 Research Park Blvd, Ste 210
Madison, WI 53719 USA