

# Lactococcus lactis strain NZ3910, food grade

Product Information Sheet  
# VS-ELS03910-01



## SUMMARY

shipped on dry ice; store at -80 °C

For research use only!

## Product Description and Application

- *Lactococcus lactis* strain NZ3910, food grade. Host strain for NICE® *Lactococcus* gene expression systems.

For handling carefully read the handbook pertaining to the gene expression system. Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

## Genotype

- *lacF*<sup>-</sup>, *nisRnisk* (in neutral locus)

## Contents

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

## Order Information, Shipping and Storage

| Order#                              | Product  | Quantity |
|-------------------------------------|--|----------|
| VS-ELS03910-01                      | NICE® <i>Lactococcus lactis</i> expression Strain NZ3900, food grade | 1 ml     |
| shipped on dry ice; store at -80 °C |  |          |

## Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: order@mobitec.com

### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: info@mobitec.com

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)