


# Escherichia coli strain MC1061

Product Information Sheet  
# VS-ELS10610-01



## SUMMARY

shipped on dry ice; store at -80 °C  
For research use only!

## Product Description and Application

- *Escherichia coli* strain MC1061. Cloning strain for diverse Gram-positive bacteria, e.g., *Lactococcus lactis* or *Bacillus subtilis*.

For handling carefully read the handbook pertaining to the gene expression system.  
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

## Genotype

- *araD139, Δ(ara, leu)7697, ΔlacX74, galU-, galK-, hsr-, hsm+, strA* (Casadaban MJ and Cohen SN, 1980; Analysis of gene control signals by DNA fusion and cloning in *Escherichia coli*.)

## Contents

- 1 ml cell suspension in growth medium supplemented with 10-30% glycerol.

## Order Information, Shipping and Storage

Order#	Product	Amount
VS-ELS10610-01	<i>Escherichia coli</i> strain MC1061	1 ml
shipped on dry ice; store at -80 °C		

## Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: order@mobitec.com

### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: info@mobitec.com

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)