

SUMMARY

shipped at RT; store at 4 °C

For research use only

MobiFairy stabilizes your proteins!

Product Description

MobiFairy protects and stabilizes biomolecules such as enzymes and antibodies (for example from physical stress like freeze-shock and freeze-thawing). MobiFairy is an effective cryoprotectant which is inert, nontoxic and biologically stable. It can be dialyzed easily. MobiFairy effectively preserves protein function by preventing degradation and irreversible aggregation of many proteins and peptides.

For protein stabilization, MobiFairy has to be diluted 1:5. Enzyme solutions can be stored in a freezer for a longer period of time without significantly losing activity. Moreover, MobiFairy might also be useful for lyophilizing certain proteins.

Order Information, Shipping and Storage

Order#	Product	Quantity
MOBI05	MoBiFairy	2 mL
MOBI06	MoBiFairy	5 x 2 mL
shipped at RT; store at 4 °C		

Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0
fax: +49 (0)551 707 22 22
e-mail: order@mobitec.com

Technical Service – Product information

phone: +49 (0)551 707 22 70
fax: +49 (0)551 707 22 77
e-mail: info@mobitec.com

MoBiTec in your area: Find your local distributor at www.mobitec.com