


# MobiSpinColumn "F" with fixed outlet plug, inserted large 10 um filter and screw cap

Product Information Sheet  
# M105210S

## SUMMARY

shipped at RT; store at RT


For research use only



The MobiSpin Column "F" is designed for a wide variety of applications for nucleic acid and protein extraction and purification. When choosing the MobiSpin Column "F" for a particular application, a suitable resin must be selected and filled into the mini-column. For an appropriate resin selection, impurity versus sample size and the anticipated yield must be considered.

Spin column chromatography offers many advantages over traditional liquid chromatography:

- Comes empty - fill in your own material (matrix, gel slice etc.)
- No sample dilution
- Reproducible results with fast protocols
- One spin application in less than 4 minutes
- Numerous samples can be processed simultaneously
- Large number of applications
- Compatible with laboratory standard



# MobiSpinColumn "F" with fixed outlet plug, inserted large 10 µm filter and screw cap

## Product Information Sheet # M105210S



The MobiSpin Column "F" comes with inserted filter of 10 µm pore size, screw cap and snap-off plug. It can be used without matrix (e.g. to extract trypsinated proteins out of a gel slice for Mass Spec analysis) or filled with a matrix of your choice.

### Order Information, Shipping and Storage

Order#	Product	Quantity
M105210S	MobiSpin Column "F" with fixed outlet plug, inserted large 10 µm filter and screw cap	50
M105010S	MobiSpinColumn "F" with fixed outlet plug, inserted small 10 µm filter and screw cap	50
shipped at RT; store at RT		

### Column Accessories

Order#	Product	Quantity
M2210	Filter (large) 10 µm pore size, for Mobicol M1002, M1003, M1050 & M1053	50
M523515	Filter (large) 35 µm pore size, for Mobicol M1002, M1003, M1050 & M1053	50
M2290	Filter (large) 90 µm pore size, for Mobicol M1002, M1003, M1050 & M1053	50
shipped at RT; store at RT		

### Related Products: Pre-packed ready-to-use columns

Order#	Product	Quantity
SCO200	MobiSpin S-200	20 columns
SCO210	MobiSpin S-200	100 columns
SCO300	MobiSpin S-300	20 columns
SCO310	MobiSpin S-300	100 columns
SCO400	MobiSpin S-400	20 columns
SCO410	MobiSpin S-400	100 columns
SCO234	MobiSpin S-200, S-300, S-400	3 x 10 columns
SCO500	MobiSpin G-50	20 columns
SCO510	MobiSpin G-50	100 columns
shipped at RT; store at 4 °C		

### Contact and Support

MoBiTec GmbH ● Lotzestr. 22a ● D-37083 Goettingen ● Germany

#### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: [order@mobitec.com](mailto:order@mobitec.com)

#### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: [info@mobitec.com](mailto:info@mobitec.com)

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)