

### SUMMARY

shipped at RT

For research use only

#### Absorption / emission max.:

691 nm / 708 nm (in Ethanol)

#### Molar absorbance:

140.000 M<sup>-1</sup>cm<sup>-1</sup>

#### Comments:



- Soluble in water, methanol, DMF, DMSO
- Bright solid state emission
- Suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis
- Negatively charged
- improved water solubility



#### Available Modifications

Available Modifications	Molecular Weight (g • mol <sup>-1</sup> )	Molecular formula
Carboxylic acid	736.88	C <sub>36</sub> H <sub>45</sub> N <sub>2</sub> O <sub>9</sub> S <sub>2</sub> Na
NHS-ester	833.95	C <sub>40</sub> H <sub>48</sub> N <sub>3</sub> O <sub>11</sub> S <sub>2</sub> Na
Amino-derivative	756.99	C <sub>38</sub> H <sub>52</sub> N <sub>4</sub> O <sub>8</sub> S <sub>2</sub>
Maleimide	859.01	C <sub>42</sub> H <sub>51</sub> N <sub>4</sub> O <sub>10</sub> S <sub>2</sub> Na
dUTP	1362.94	C <sub>54</sub> H <sub>70</sub> N <sub>6</sub> O <sub>23</sub> P <sub>3</sub> S <sub>2</sub> Li <sub>5</sub>

#### Diagram, Formel



## Literature

Vollrath A., Schubert S., Windhab N., Biskup C., Schubert U. S., Labeled Nanoparticles Based on Pharmaceutical EUDRAGIT R S 100 Polymers, *Macromol. Rapid Commun.* 31(2010)2053–2058.

Van Noorden C. J. F., *Imaging Enzymes at Work: Metabolic Mapping by Enzyme Histochemistry*, *J Histochem Cytochem* 58(2010)481-497.

Napp J., Dullin C., Müller F., Uhland K., Petri J. B., van de Locht A., Steinmetzer T., Alves F., Time-domain in vivo near infrared fluorescence imaging for evaluation of matriptase as a potential target for the development of novel, inhibitor-based tumor therapies, *International Journal of Cancer*, 127,p1958–1974.

Schuster C. F., Park J.-H., Prax M., Herbig A., Nieselt K., Rosenstein R., Inouye M., Bertram R., Characterization of a mazEF Toxin-Antitoxin Homologue from *Staphylococcus equorum*, *J. Bacteriol.* 195(2013)115-125.

Kaltenpoth M., Strohm E., Gadau J., Polymorphic microsatellite markers for a solitary digger wasp, the European beewolf (*Philanthus triangulum*; Hymenoptera, Sphecidae), *Molecular Ecology Notes*,4(2004)589–592

Ritz C. M., Wissemann V., Microsatellite Analyses of Artificial and Spontaneous Dogrose Hybrids Reveal the Hybridogenic Origin of *Rosa micrantha* by the Contribution of Unreduced Gametes, *J Hered.* 102(2011)217-227.

Harloff H.-J., Lemcke S., Mittasch J., Frolov A., Wu J. G., Dreyer F., Leckband G., Jung C., A mutation screening platform for rapeseed (*Brassica napus* L.) and the detection of sinapine biosynthesis mutants, *Theor Appl Genet* 124(2012)957–969.

Arndt V., Daniel C., Nastainczyk W., Alberti S., Höfeld J., BAG-2 Acts as an Inhibitor of the Chaperone-associated Ubiquitin Ligase CHIP, *Molecular Biology of the Cell* 16(2005)5891–5900.

Frerichmann S. L. M., Kirchoff M., Müller A. E., Scheidig A. J., Jung C., Kopisch-Obuch F. J., EcoTILLING in *Beta vulgaris* reveals polymorphisms in the FLC-like gene BvFL1 that are associated with annuality and winter hardiness, *BMC Plant Biology* 13(2013)52.

Owczarek-Lipska M., Dolf G., Guziewicz K. E., Leeb T., Schelling C., Posthaus H., Braunschweig M. H., Bovine cardiac troponin I gene (TNNI3) as a candidate gene for bovine dilated cardiomyopathy, *Archiv Tierzucht* 52(2009)113-123.

de Almeida R. A., Heuser T., Blaschke R., Bartram C. R., Janssen J. W. G., Control of MYEOV Protein Synthesis by Upstream Open Reading Frames, 281(2006)695-704.

Padirac A., Fujii T., Rondelez Y., Quencher-free multiplexed monitoring of DNA reaction circuits, *Nucl. Acids Res.* (2012) 40 (15): e118.

Millera O. J., Harrakb A. E., Mangeatb T., Bareta J.-C., Frenza L., Debsa B. E., Mayota E., Samuels M. L., Rooneyd E. K., Dieue P., Galvand M., Link D. R., Griffiths A. D., High-resolution dose–response screening using droplet-based microfluidics; *PNAS* 109(2012)378-383.

Franz M., Hilger I., Grün K., Kossatz S., Richter P., Petersen I., Jung C., Gummert J., Figulla H. R., Kosmehl H., Neri D., Berndt A., Renner A., Selective imaging of chronic cardiac rejection using a human antibody specific to the alternatively spliced EDA domain of fibronectin, *The Journal of Heart and Lung Transplantation*, 32(2013)641–650.

Ferreira D. N., Katayama I. A., Oliveira I. B., Rosa K. T., Furukawa L. N. S., Coelho M. S., Casarini D. E., Heimann J. C., Salt-Induced Cardiac Hypertrophy and Interstitial Fibrosis Are Due to a Blood Pressure-Independent Mechanism in Wistar Rats, *J. Nutr.* 140(2010)1742-1751.

Pogacica T., Mancinia A., Santarellia M., Bottaria B., Lazzia C., Neviana E., Gattia M., Diversity and dynamic of lactic acid bacteria strains during aging of a long ripened hard cheese produced from raw milk and undefined natural starter, *Food Microbiology*, 36(2013)207–215.

Kandaswamy R., McQuillin A., Sharp S. I., Fiorentino A., Anjorin A., Blizard R. A., Curtis D., Gurling H. M. D., Genetic Association, Mutation Screening, and Functional Analysis of a Kozak Sequence Variant in the Metabotropic Glutamate Receptor 3 Gene in Bipolar Disorder, *JAMA Psychiatry*. 70(2013)591-598.

Mancini A., Lazzi C., Bernini V., Neviani E., Gatti M., Identification of dairy lactic acid bacteria by tRNAAla-23S rDNA-RFLP, *Journal of Microbiological Methods* 91(2012)380–390.

Glatzer S., Merten N. J., Dierks C., Wöhlke A., Philipp U., Distl O., A Single Nucleotide Polymorphism within the Interferon Gamma Receptor 2 Gene Perfectly Coincides with Polledness in Holstein Cattle, *PLoS ONE* 8(6): e67992.

van de Kooij B., Rooswinkel R. W., Kok F., Herrebut M, de Vries E., Paauwe M., Janssen G. M. C., van Veelen P. A., Borst J., Polyubiquitination and proteasomal turnover controls the anti-apoptotic activity of Bcl-B, *Oncogene* (2013)1–10.

Koplin A., Preissler S., Ilina Y., Koch M., Scior A., Erhardt M., Deuerling E., A dual function for chaperones SSB-RAC and the NAC nascent polypeptide-associated complex on ribosomes, *JCB*, 189(2010)57-68.

### Order Information, Shipping and Storage

Order#	Product	Quantity
MFP-D681-00-1	MFP™-DY-681-Carboxylic Acid	1 mg
shipped at RT; store at RT		
MFP-D681-01-1	MFP™-DY-681-NHS-Ester	1 mg
shipped at RT; store at 4 °C		
MFP-D681-02	MFP™-DY-681-Amino function	1 mg
MFP-D681-03-2	MFP™-DY-681-Maleimide	1 mg
shipped at RT; store at RT		
MFP-D681-04	MFP™-DY-681-Hydrazide	1 mg
MFP-D681-05	MFP™-DY-681-STP-Ester	1 mg
shipped at RT; store at 4 °C		
MFP-D681-30-2	MFP™-DY-681-Biotin	1 mg
shipped at RT; store at RT		
MFP-D681-34	MFP™-DY-681-dUTP	100 nmol
shipped on blue ice; store at -20 °C		
MFP-D681-10	MFP™-DY-681-Azide	1 mg

### Contact and Support

MoBiTec GmbH ● Lotzestr. 22a ● D-37083 Goettingen ● Germany

#### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: [order@mobitec.com](mailto:order@mobitec.com)

#### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: [info@mobitec.com](mailto:info@mobitec.com)

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)