

### SUMMARY

shipped on dry ice; store at -80 °C

For research use only!

## Product Description and Application

- *Bacillus subtilis*, strain NZ8900. Host strain for *Bacillus subtilis* SURE gene expression system.
- *Bacillus subtilis* 168 strain with *spaR* and *spaK* integrated in the chromosome in the *amyE* locus.

For handling carefully read the handbook supplied with the gene expression system.  
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

## Genotype

- *amyE::spaRK*, Km<sup>R</sup>

## Contents

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

## Order Information, Shipping and Storage

Order#	Product	Quantity
PBS024	<i>Bacillus subtilis</i> strain NZ8900	1 ml
shipped on dry ice; store at -80 °C		

## Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: [order@mobitec.com](mailto:order@mobitec.com)

### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: [info@mobitec.com](mailto:info@mobitec.com)

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)