

# Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns ( 90µm pore size)

## Product Information Sheet # S10019

### SUMMARY



shipped at RT; store at RT

For research use only

Easy-To-Handle Lab Columns with Luer-Lock Connections for Syringes and Tubing.

### Product Description and Application

Set of empty columns for e.g. affinity and gel chromatography, ion exchange columns, gel filtration, etc.



Mobicol: 1 ml (M1050)

### Advantages:

- easy handling
- compatible with laboratory standard
- smallest and largest volumes can be treated easily
- Luer-lock connectors for tubings and syringes
- can be closed tightly
- come with lower filters of 10, 35 or 90 µm pore size inserted
- low price
- complete system

### Set contents:


**3 columns each** of 1 ml column (Mobicol), 2.5 ml, 5 ml and 10 ml columns & **two stopcocks**. All columns come with Luer-lock cap, plug, closing cap and inserted lower filters. Upper filter (used on top of the matrix) available separately in sets as well.

**Set # 1:** 10 µm pore size filters inserted

**Set # 2:** 35 µm pore size filters inserted

**Set # 3:** 90 µm pore size filters inserted

# Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns ( 90µm pore size)



## Product Information Sheet # S10019

### Order Information, Shipping and Storage

#### Column Sets

Order#	Product	Quantity
S10011	Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns (10 µm pore size)	3 each
S1001	Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns (35 µm pore size)	3 each
S10019	Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns (90 µm pore size)	3 each
S10031	Set upper filters 10 µm for S10011	1 set
S1003	Set upper filters 35 µm for S1001	1 set
S10039	Set upper filters 90 µm for S10019	1 set

shipped at RT; store at RT

#### Column Accessories

Order#	Product	Quantity
M3001	Tubing-to-Luer-lock adapter (material outside: metal; inside: glass and PTFE; with 1 m PTFE tubing) for Mobicols or laboratory columns	1
M3002	Luer adapter, (M+F) 10 each	20
M3003	Luer adapter for the top-cap (M)	20
M3004	Luer adapter for outlet (F)	20
M3005	Bottom plugs for M1002	50
M3006	Stopcocks, one-way, Luer-lock	4
M3007	Luer-lock adapter (M/M)	1

shipped at RT; store at RT

#### Laboratory Columns

Order#	Product	Quantity
S10121	2.5 ml columns (10 µm filter pore size)	20
S1012	2.5 ml columns (35 µm filter pore size)	20
S10129	2.5 ml columns (90 µm filter pore size)	20
S10131	5.0 ml columns (10 µm filter pore size)	20
S1013	5.0 ml columns (35 µm filter pore size)	20
S10139	5.0 ml columns (90 µm filter pore size)	20
S10141	10 ml columns (10 µm filter pore size)	20
S1014	10 ml columns (35 µm filter pore size)	20
S10149	10 ml columns (90 µm filter pore size)	20
S1210	Upper filter 10 µm for 2.5 ml columns	20
S1235	Upper filter 35 µm for 2.5 ml columns - (thickness 1.5 mm)	20
S123532	Upper filter 35 µm for 2.5 ml columns - (thickness 3.2 mm)	20
S1290	Upper filter 90 µm for 2.5 ml columns	20
S1310	Upper filter 10 µm for 5.0 ml columns	20
S1335	Upper filter 35 µm for 5.0 ml columns - (thickness 1.5 mm)	20
S133532	Upper filter 35 µm for 5.0 ml columns - (thickness 3.2 mm)	20
S1390	Upper filter 90 µm for 5.0 ml columns	20
S1410	Upper filter 10 µm for 10 ml columns	20
S143532	Upper filter 35 µm for 10 ml columns - (thickness 3.2 mm)	20
S1490	Upper filter 90 µm for 10 ml columns	20

shipped at RT; store at RT

### Contact and Support

MoBiTec GmbH ● Lotzestr. 22a ● D-37083 Goettingen ● Germany

#### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: order@mobitec.com

#### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: info@mobitec.com

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)