

# Lactococcus lactis expression strain MG1363

Product Information Sheet  
# VS-ELS01363



## SUMMARY

shipped on dry ice; store at -80 °C

For research use only!

## Product Description and Application

- *Lactococcus lactis* subsp. *cremoris* expression strain MG1363. Host strain for *Lactococcus* gene expression systems.

For handling carefully read the handbook pertaining to the gene expression system.  
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

## Genotype

- Plasmid-free derivative of SH4109 (Gasson M. J.: Plasmid complements of *Streptococcus lactis* NCDO 712 and other lactic streptococci after protoplast-induced curing, Journal of Bacteriology (1983), Vol. 154, No. 1, 1 - 9)

## Contents

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

## Order Information, Shipping and Storage

Order#	Product	Quantity
VS-ELS01363	<i>Lactococcus lactis</i> expression strain MG1363	1 ml
shipped on dry ice; store at -80 °C		

## Contact and Support

MoBiTec GmbH ● Lotzestr. 22a ● D-37083 Goettingen ● Germany

### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: [order@mobitec.com](mailto:order@mobitec.com)

### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: [info@mobitec.com](mailto:info@mobitec.com)

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)