


Bacillus megaterium protoplasts, strain WH320

Product Information Sheet
BMEG02



SUMMARY

shipped on dry ice; store at -80 °C
For research use only!

Product Description and Application

- *Bacillus megaterium* protoplasts, strain WH320. Host strain for *Bacillus megaterium* gene expression systems

For handling carefully read the handbook supplied with the gene expression system.
Care must be taken not to interrupt the cooling chain. Always store frozen protoplasts at -80 °C or lower.

Genotype

- n/a – WH320 is a chemical mutant of DSM319, lacZ⁻ (Rygus T, Hillen W.: Inducible high-level expression of heterologous genes in *Bacillus megaterium* using the regulatory elements of the xylose-utilization operon, Appl. Microbiol. Biotechnol. (1991), 35, 594-599)

Contents

- 5x500 µl ready-to-use protoplasts

Order Information, Shipping and Storage

Order#	Product	Quantity
BMEG02	<i>Bacillus megaterium</i> protoplasts, strain WH320	5x500 µl
shipped on dry ice; store at -80 °C		

Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0
fax: +49 (0)551 707 22 22
e-mail: order@mobitec.com

Technical Service – Product information

phone: +49 (0)551 707 22 70
fax: +49 (0)551 707 22 77
e-mail: info@mobitec.com

MoBiTec in your area: Find your local distributor at www.mobitec.com