

# Bacillus megaterium protoplasts, strain YYBm1

Product Information Sheet  
# BMEG04



## SUMMARY

shipped on dry ice; store at -80 °C

For research use only!

## Product Description and Application

- *Bacillus megaterium* protoplasts, strain YYBm1. Host strain for *Bacillus megaterium* gene expression systems. Ready for transformation.
- YYBm1 is unable to utilize xylose as carbon source.

For handling carefully read the handbook supplied with the gene expression system.  
Care must be taken not to interrupt the cold chain. Always store frozen protoplasts at -80 °C or lower.

## Genotype

- Defined deletion mutant of wild type strain DSM319:  $\Delta nprM$ ,  $\Delta xyIA$

## Contents

- 5 x 500  $\mu$ l protoplasts of *Bacillus megaterium* YYBm1 supplemented with 11% glycerol.

## Order Information, Shipping and Storage

Order#	Product	Quantity
BMEG04	<i>Bacillus megaterium</i> protoplasts, strain YYBm1	5 x 500 $\mu$ l
shipped on dry ice; store at -80 °C		

## Contact and Support

MoBiTec GmbH ● Lotzestr. 22a ● D-37083 Goettingen ● Germany

### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: [order@mobitec.com](mailto:order@mobitec.com)

### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: [info@mobitec.com](mailto:info@mobitec.com)

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)