

# Bacillus megaterium protoplasts, strain MS941

Product Information Sheet  
# BMEG50



## SUMMARY

shipped on dry ice; store at -80 °C  
For research use only!

## Product Description and Application

- *Bacillus megaterium* protoplasts, strain MS941. Host strain for *Bacillus megaterium* gene expression systems

For handling carefully read the handbook supplied with the gene expression system.  
Care must be taken not to interrupt the cooling chain. Always store frozen protoplasts at -80 °C or lower.

## Genotype

- Defined deletion mutant of wildtype strain DSM319,  $\Delta nprM$

## Contents

- 5x500 µl ready-to-use protoplasts

## Order Information, Shipping and Storage

Order#	Product	Quantity
BMEG50	<i>Bacillus megaterium</i> protoplasts, strain MS941	5x500 µl
shipped on dry ice; store at -80 °C		

## Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

### Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0  
fax: +49 (0)551 707 22 22  
e-mail: order@mobitec.com

### Technical Service – Product information

phone: +49 (0)551 707 22 70  
fax: +49 (0)551 707 22 77  
e-mail: info@mobitec.com

**MoBiTec in your area:** Find your local distributor at [www.mobitec.com](http://www.mobitec.com)