


Bacillus subtilis strain NZ8963

Product Information Sheet
PBS023



SUMMARY

shipped on dry ice; store at -80 °C
For research use only!

Product Description and Application

- *Bacillus subtilis*, strain NZ8963. Subtilin producing strain; suitable for *Bacillus subtilis* SURE gene expression system.

For handling carefully read the handbook supplied with the gene expression system.
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

Genotype

- Wild type

Contents

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

Order Information, Shipping and Storage

Order#	Product	Quantity
PBS023	<i>Bacillus subtilis</i> strain NZ8963	1 ml
shipped on dry ice; store at -80 °C		

Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0
fax: +49 (0)551 707 22 22
e-mail: order@mobitec.com

Technical Service – Product information

phone: +49 (0)551 707 22 70
fax: +49 (0)551 707 22 77
e-mail: info@mobitec.com

MoBiTec in your area: Find your local distributor at www.mobitec.com