

Bacillus subtilis strain NZ8901

Product Information Sheet
PBS025



SUMMARY

shipped on dry ice; store at -80 °C
For research use only!

Product Description and Application

- *Bacillus subtilis*, strain NZ8901. Negative control strain; suitable for *Bacillus subtilis* SURE gene expression system.
- *Bacillus subtilis* 168 strain with kanamycin gene integrated in the chromosome in the *amyE* locus.

For handling carefully read the handbook supplied with the gene expression system.
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

Genotype

- *amyE::Km^R*

Contents

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

Order Information, Shipping and Storage

Order#	Product	Quantity
PBS025	<i>Bacillus subtilis</i> strain NZ8901	1 ml
shipped on dry ice; store at -80 °C		

Contact and Support

MoBiTec GmbH ● Lotzestr. 22a ● D-37083 Goettingen ● Germany

Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0
fax: +49 (0)551 707 22 22
e-mail: order@mobitec.com

Technical Service – Product information

phone: +49 (0)551 707 22 70
fax: +49 (0)551 707 22 77
e-mail: info@mobitec.com

***Bacillus subtilis* strain NZ8901**

Product Information Sheet
PBS025



MoBiTec in your area: Find your local distributor at www.mobitec.com