

Lactococcus lactis strain NZ9100

Product Information Sheet
VS-ELS09100



SUMMARY

shipped on dry ice; store at -80 °C

For research use only!

Product Description and Application

- *Lactococcus lactis* strain NZ9100. Host strain for NICE® *Lactococcus* gene expression systems.

For handling carefully read the handbook pertaining to the gene expression system.
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

Genotype

- *nisRnisK* (in neutral locus)

Contents

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

Order Information, Shipping and Storage

Order#	Product	Quantity
VS-ELS09100	NICE® <i>Lactococcus lactis</i> expression Strain NZ9100	1 ml
shipped on dry ice; store at -80 °C		

Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

Customer Service – General inquiries & orders

phone: +49 (0)551 707 22 0
fax: +49 (0)551 707 22 22
e-mail: order@mobitec.com

Technical Service – Product information

phone: +49 (0)551 707 22 70
fax: +49 (0)551 707 22 77
e-mail: info@mobitec.com

MoBiTec in your area: Find your local distributor at www.mobitec.com